Virtual[®] 380

EN Instructions for Use

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For dental use only.

Description

Virtual® impression materials are addition-reaction silicones (vinylpolysiloxanes) used to create fine detail impressions of dentition. Virtual impression materials are available in a variety of viscosities allowing dental professionals to choose the material and technique best suited for each individual case.

The Virtual line of addition-silicone (vinylpolysiloxane) impression materials is recommended for use to create highly detailed impressions of the hard and soft tissues of the oral cavity.

- Final impressions used for the fabrication of indirect restorations (crowns, bridges, inlays, onlays and veneers)
- Dental implant impressions
- Matrix from “wax-ups” or for the treatment planning, study models
- Edentulous impressions
- Matrix used to create temporary restorations

Each viscosity provides specific attributes to meet the needs of the dental professional.

Colours and working times

Technical data

	Monophase	Heavy Body
Farbe	Blue	Blue
Classification ISO 4823 / ADA Spec. No. 19	Type 2 Medium-bodied consistency	Type 1 Heavy-bodied consistency
Mixing Ratio [Base:Catalyst]	5:1	5:1

Material	Setting Time	Total Work Time [min:sec]	Mouth Set Time* [min:sec]
Monophase [tray/wash]	Fast Set	1:00	3:30
Heavy Body [tray]	Regular Set Fast Set	1:30 1:00	4:30 2:30

* Minimum time the impression material should remain in the mouth before removal.
The processing times are reduced or prolonged depending on the prevalent room temperature.

Indication

Impressions of dentulous and edentulous jaw sections within the framework of restorative, prosthetic, orthodontic or functional treatment measures.

Contraindication

The use of Virtual Monophase 380, Virtual Heavy Body 380 Regular/Fast is contraindicated if the patient is known to be allergic to any of the products' ingredients.

Limitations of use

The product is not intended to be reprocessed or reused.

Side effects

There are no known side effects to date.

Interactions

The setting of vinylpolysiloxanes is inhibited by latex gloves.

Do not touch teeth, preparations, retractions cords, etc. with latex gloves. It is recommended that operators wash their hands thoroughly or use vinyl gloves in order to eliminate all traces of impurities. Other products which may inhibit the set of impression materials include rubber dams, retraction cords and certain substances.

If the operator suspects that the preparation has been contaminated, the preparation must be rinsed and dried to eliminate all traces of impurities.

Composition

– Virtual Monophase 380 Fast, Virtual Heavy Body 380 Regular/Fast

Vinylpolysiloxane, silicon dioxide, inorganic aluminium filler (CAS 1318-02-1), methyl hydrogen siloxane

Application

– **Heavy Body (blue):** High viscosity, hydrophilic vinylpolysiloxane suitable for final impressions. This material is recommended for use as a tray material when using the Tray/Wash impression technique.

– **Monophase (blue):** Medium viscosity, hydrophilic vinylpolysiloxane suitable for final impressions.

This material is recommended for use as a tray and wash material when using the Double-Arch, Single-Phase impression technique.

– Instructions for Virtual 380 cartridge materials

The Virtual 380 cartridge is suitable for use in conjunction with e.g. Pentamix[®], Pentamix[®] 2, Pentamix[®] 3, MixStar-eMotion, Sympress, Duomix and Duomix II, PowerMix, Dynamix Speed and Modulmix mixing units (except Pentamix Lite).

Carefully remove the Virtual 380 cartridges from the packaging.

IMPORTANT: Avoid dropping the heavy Virtual 380 cartridges, as this may severely damage the cartridges, so that they can no longer be used.

Pentamix, Pentamix 2, Pentamix 3, Pentamix Lite, MixStar-eMotion, Sympress, Duomix and Duomix II, PowerMix, Dynamix Speed and Modulmix are not registered trademarks of Ivoclar Vivadent.

- **Inserting the Virtual 380 cartridge:**

(Please also refer to the corresponding directions for use of the manufacturer of the mixing unit.)

Fig. 1

Fig. 2

Fig. 3

Fig. 4

Fig. 5

1. Hold the cartridge with the sealed outlet openings facing upwards.
2. Grasp the flap on the sealing cap and bend it upward 90 degrees. Do not apply excessive force in the process.
3. While holding the cartridge securely in one hand, take the flap between thumb and index finger and completely remove the sealing cap. (Fig. 1)
4. Insert the cartridge into the mixing unit according to the manufacturer's instructions. (Fig. 2)
5. To control the flow of the pastes, start the mixing unit according to the manufacturer's instructions and wait until a small amount of material is extruded (This should always be done prior to starting the mixing process). Carefully wipe off excess material in a vertical direction to avoid that the two components become intermixed inside the openings. (Fig. 3)
6. Attach one of the dynamic mixers included in the delivery to the cartridge following the instructions of the mixing unit manufacturer. (Fig. 4)
(NOTE: Only use the dynamic mixers from Ivoclar Vivadent – As the dynamic mixers from other manufacturers do not fit the openings precisely, leakage may occur.)
7. Once the dynamic mixer is correctly seated, place the yellow locking ring over the mixer and push it down completely. Turn the ring 1/4 turn clockwise to lock the mixer securely in place. (Fig. 5)
8. Load a separate intra-oral elastomer syringe (only Virtual 380 Monophase) or a suitable impression tray with the mixed Virtual impression material.
9. To ensure an optimum mixing speed and the extrusion of an optimum amount of material, select the quickest program.
10. Immediately after extrusion of the material, fully retract the plungers of the mixing unit to avoid continual pressure on the cartridge and prevent the escape of excess material through the dynamic mixer.
11. Leave the used dynamic mixer on the cartridge. It serves as a seal until next use. In this way, cross-contamination of base and catalyst and clogging of the outlet openings with set material is avoided.

NOTE: It is recommended that the cartridge be removed from the mixing unit at the end of each working day or when the mixing unit is not used for a longer period of time. It should be stored in an upright position, with the outlet openings (sealed either by a dynamic mixer or the sealing cap) facing upwards.

- **Pre-treatment of impression trays (tray adhesives)**

It is recommended that a tray adhesive for silicone-based impression materials be used to reduce the chance of distortion when removing impressions from the mouth.

- **Disinfection**

Impressions made with Virtual impression materials can be immersed in a disinfectant solution (e.g. glutaraldehyde 0.5%, benzalkonium chloride 0.5%). Please observe the instructions of the manufacturer. The disinfection does not affect the surface or the dimension.

- **Model fabrication**

The impression may be poured immediately after disinfection, or up to two weeks later, provided that the impression is stored at room temperature. Dimensional stability is guaranteed for 14 days, if stored appropriately. Virtual impression materials are compatible with all popular dental plasters on the market, e.g. Type 3: Elite® Model (Zhermack), Type 4: Fujirock® (G.C. International).

- **Galvanization**

Virtual impression materials can be silver or copper plated in a galvanic bath.

Special notes

Virtual impression materials should be at room temperature (23 °C / 73 °F) when processed. Lower temperatures, e.g. when stored in the refrigerator, will prolong the working time and the time the impression material should remain in the mouth, while higher temperatures will reduce it. Vinylpolysiloxanes are chemically resistant. Unpolymerized material may stain clothing.

Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclarvivadent.com, and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclarvivadent.com).

Warnings

If uncured material comes into contact with the eyes, rinse the eye affected immediately with plenty of water. If irritation persists, consult a physician. If the material comes into contact with the skin, rinse with plenty of water.

Disposal information

- Remaining stocks must be disposed of according to the corresponding national legal requirements.

Shelf life and storage

- Storage temperature: 2-28 °C (36-82 °F)
- Do not use the product after the indicated date of expiration.
- Expiration date: See note on cartridge

Additional information

Keep material out of the reach of children!

The product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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