Virtual®

EN Instructions for Use

- Vinylpolysiloxane impression material (addition-reaction)

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For dental use only.

Description

Virtual impression materials are addition-reaction silicones (vinylpolysiloxanes) used to create fine detail impressions of dentition. Virtual impression materials are available in a variety of viscosities allowing dental professionals to choose the material and technique best suited for each individual case.

The Virtual line of addition-silicone (vinylpolysiloxane) impression materials is recommended for use to create highly detailed impressions of the hard and soft tissues of the oral cavity.

- Final impressions used for the fabrication of indirect restorations (crowns, bridges, inlays, onlays and veneers)
- Dental implant impressions
- Matrix from "wax-ups" or for the treatment planning, study models
- Edentulous impressions
- Matrix used to create temporary restorations

Each viscosity provides specific attributes to meet the needs of the dental professional.

Colours and working times

Technical data

	Extra Light Body	Light Body	Monophase	Heavy Body	Putty
Colour	Beige	Beige	Blue	Blue	Blue
Classification ISO 4823 / ADA Spec. No. 19	Type 3 light-bodied consistency	Type 3 light-bodied consistency	Type 2 medium- bodied consistency	Type 1 heavy-bodied consistency	Type 0 putty consistency
Mixing Ratio [Base:Catalyst]	1:1	1:1	1:1	1:1	1:1

Material	Setting Time	Total Work Time [min:sec]	Mouth Set Time* [min:sec]
Extra Light Body [wash]	Regular Set	1:30	4:30
	Fast Set	1:00	2:30
Light Body [wash]	Regular Set	1:30	4:30
	Fast Set	1:00	2:30
Monophase [tray/wash]	Regular Set	1:30	4:30
	Fast Set	1:00	2:30
Heavy Body [tray]	Regular Set	1:30	4:30
	Fast Set	1:00	2:30
Putty [tray]	Regular Set	1:30	4:30
	Fast Set	1:00	2:30

* Minimum time the impression material should remain in the mouth before removal.

The processing times are reduced or prolonged depending on the prevalent room temperature.

Indication

Impressions of dentulous and edentulous jaw sections within the framework of restorative, prosthetic, orthodontic or functional treatment measures.

Contraindication

The use of Virtual Extra Light Body Regular/Fast, Virtual Light Body Regular/Fast, Virtual Monophase Regular/Fast, Virtual Heavy Body Regular/Fast, Virtual Putty Regular/Fast is contraindicated if the patient is known to be allergic to any of the products' ingredients.

Limitations of use

The product is not intended to be reprocessed or reused.

Side effects

There are no known side effects to date.

Interactions

The setting of vinylpolysiloxanes is inhibited by latex gloves.

Do not touch teeth, preparations or retractions cords with latex gloves. It is recommended that operators wash their hands thoroughly or use vinyl gloves in order to eliminate all traces of impurities, specifically when hand mixing putty. Other products which may inhibit the set of impression materials include rubber dams, retraction cords and certain substances. If the operator suspects that the preparation has been contaminated, the preparation must be rinsed and dried to eliminate all traces of impurities.

Composition

- **Virtual Extra Light Body Regular/Fast, Virtual Light Body Regular/Fast, Virtual Monophase Regular/Fast, Virtual Heavy Body Regular/Fast, Virtual Putty Regular/Fast**
Vinylpolysiloxane, silicon dioxide, inorganic aluminium filler (CAS 1318-02-1), methyl hydrogen siloxane

Application

– Putty (Regular and Fast Set)

Very high viscosity vinylpolysiloxane that provides easy mixing and the hydraulic forces necessary to extend wash material into the sulcus providing a detailed impression of prepared margins. The putty viscosity maintains outstanding detail reproduction when used alone. Recommended impression techniques include Putty/Wash One-Step and Putty/Wash Two-Step. Important: Avoid wearing latex gloves when handling putty, putty jars, and dosing spoons (see note under interactions).

1. Take equal amounts of base (blue) and catalyst (white) using the colour-coded dosing spoons.
2. Hand mix the identical portions of Virtual Putty base and catalyst until you obtain an evenly coloured mixture (approximately 30 seconds of mixing time). A proper mixture should have no streaks. Note that equal proportions of base and catalyst must be used in order to obtain proper setting times. A larger quantity of catalyst will not accelerate the setting time. However, the setting time and curing depend on the temperature of the fingers/hand and the mixing intensity.
3. Place mixture in impression tray. It is recommended that a tray adhesive for silicone-based impression materials be used to reduce the chance of distortion when removing impressions from the mouth.

Make sure to re-seal the Virtual Putty jars immediately after use using the lid with the correct colour. Moreover, use the spoon with the correct colour to measure out the material. Failure to do so can cause discolouration on the surface of the white Putty Catalyst paste. However, this discolouration does not impair the functions of the product.

– Extra-Light Body, Light Body, Heavy Body, Monophase (Regular and Fast Set)

Extra-Light Body (beige): Extra low viscosity, hydrophilic vinylpolysiloxane suitable for final impressions. This material is recommended for use as a wash material when using the Tray/Wash and Putty/Wash impression techniques.

Light Body (beige): Low viscosity, hydrophilic vinylpolysiloxane suitable for final impressions. This material is recommended for use as a wash material when using the Tray/Wash and Putty/Wash impression techniques.

Heavy Body (blue): High viscosity, hydrophilic vinylpolysiloxane suitable for final impressions. This material is recommended for use as a tray material when using the Tray/Wash impression technique.

Monophase (blue): Medium viscosity, hydrophilic vinylpolysiloxane suitable for final impressions. This material is recommended for use as a tray and wash material when using the Double-Arch, Single-Phase impression technique.

– Instructions for cartridge materials

Loading the cartridge dispenser

Fig. 1

Fig. 2

Fig. 3

Fig. 4

Fig. 5

1. Press the black release lever located below the plunger on the back of the dispenser and pull plunger as far back as possible. (Fig. 1)
2. Lift the cartridge lock and insert the cartridge with the "V" shape on the cartridge base turned downwards. Lower the cartridge lock. (Fig. 2)
3. Remove cartridge cap by turning it 1/4 turn counter clockwise. Discard the cap. (Fig. 3)
4. Important: It is very important to clear or bleed the cartridge prior to applying the mixing tip! This will assure a proper mix and set time. Gently press on the dispenser lever until both components (base and catalyst) begin to flow out of the cartridge evenly onto a mixing pad.
5. Select the appropriate mixing tip (Fig. 4):
 - a. The large mixing tips (blue-coloured base) are recommended for Heavy Body and Monophase impression materials.
 - b. The small mixing tips (yellow-coloured base) are recommended for Extra-Light Body and Light Body impression materials.(The colour of the mixing tip base and the colour of the cartridge cap are identical). Insert the mixing tip. Push the tip down until the notch on the mixing tip touches the notch on the cartridge. While gripping the coloured base, and not the mixing tip, turn the mixing tip base 1/4 turn clockwise.
6. If necessary apply an intra-oral tip to the mixing tip (Fig. 5):
 - a. The large intra-oral tips (clear) are to be used with the large mixing tips (blue base).
 - b. The small intra-oral tips (yellow) are to be used with the small mixing tips (yellow base)
7. Begin mixing by pressing gently on the lever of the dispenser.

Removing the cartridge

1. Lift the release lever to pull the slide back. Lift the cartridge lock and remove cartridge.
2. Leave the mixing tip on the cartridge. The material will polymerize within the mixing tip, functioning as a natural cap.

- Pre-treatment of impression trays (tray adhesives)

It is recommended that tray adhesive for silicone-based impression materials be used to reduce the chance of distortion when removing impressions from the mouth.

- Disinfection

Impressions made with Virtual impression materials can be immersed in a disinfectant solution (glutaraldehyde 0.5%, benzalkonium chloride 0.5%). (Please observe the instructions of the manufacturer). The disinfection does not affect the surface or the dimension.

- Model fabrication

The impression may be poured immediately after disinfection, or up to two weeks later, provided that the impression is stored at room temperature. Dimensional stability is guaranteed for 14 days, if stored appropriately. Virtual impression materials are compatible with all popular dental plasters on the market, e.g. Type 3: Elite® Model (Zhermack), Type 4: Fujirock® (G.C. International).

- Galvanization

Virtual impression materials can be silver or copper plated in a galvanic bath.

Special notes

Virtual impression materials should be at room temperature (23 °C / 73.4 °F) when processed. Lower temperatures, e.g. when stored in the refrigerator, will prolong the working time and the time the impression material should remain in the mouth, while higher temperatures will reduce it. Vinylpolysiloxanes are chemically resistant. Unpolymerized material may stain clothing.

Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclarvivadent.com, and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclarvivadent.com).

Warnings

If uncured material comes into contact with the eyes, rinse the eye affected immediately with plenty of water. If irritation persists, consult a physician. If the material comes into contact with the skin, rinse with plenty of water.

Disposal information

- Remaining stocks must be disposed of according to the corresponding national legal requirements.

Shelf life and storage

- Storage temperature: 2 – 28 °C (36 – 82 °F).
- Do not use the product after the indicated date of expiration.
- Expiration date: See note on cartridge or jar.

Additional information

Keep material out of the reach of children!

The product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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