# Compoglass<sup>®</sup> Flow

English

## EN Instructions for Use

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## Description

Compoglass<sup>®</sup> Flow is a flowable, light-curing, radiopaque compomer-based restorative material. Compoglass Flow cures with light in the wavelength range of 400–500 nm (blue light). Compoglass Flow combines the advantages of glass ionomer cements with those of light-curing composite materials.

## Shades

Compoglass Flow is available in the following shades (Chromascop/A-D shades):  
140/A2, 210/A3, 230/A3.5, 340/A4

## Composition

The monomer matrix consists of urethane dimethacrylate, tetraethylene glycol dimethacrylate and cycloaliphatic dicarboxylic acid dimethacrylate (32.9 wt.%). The inorganic fillers include ytterbium trifluoride, Ba-Al-fluorosilicate glass and spheroid mixed oxide (66.8 wt.%). Additional contents are catalysts, stabilizers, and pigments (0.3 wt.%). The total content of inorganic fillers is 41.8 Vol.% or 66.8 wt.% respectively. The particle size is between 0.2 and 3.0 µm.

## Indication

- Restoration of deciduous teeth
- Class V restorations (cervical caries, root erosion, wedge-shaped defects)
- Anterior Class III restorations
- Temporary restorations

## Contraindication

- The placement of Compoglass Flow restorations is contraindicated if a patient is known to be allergic to any of the ingredients of Compoglass Flow;
- if the stipulated working technique cannot be applied;
- for direct and indirect pulp capping;
- for permanent, stress-bearing occlusal restorations in teeth of the secondary dentition.

## Side effects

In rare cases, ingredients of Compoglass Flow may cause a sensitizing reaction. The product should not be used in such cases. To avoid possible irritation of the pulp, areas close to the pulp should be protected with a suitable pulp/dentin protector (apply a calcium hydroxide-based preparation to areas in close proximity of the pulp) and covered with an adequate lining.

## Interactions

Substances containing eugenol/clove oil may inhibit the polymerization of compomers. These materials must not be used in conjunction with Compoglass Flow. Discolouration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.

CE 0123

Rx ONLY

Date information prepared:  
2021-06-05/Rev. 1

573128/WE3

 **Manufacturer**  
Ivoclar Vivadent AG  
FL-9494 Schaan/Liechtenstein  
www.ivoclarvivadent.com

  
ivoclar  
vivadent<sup>®</sup>  
clinical

## Usage

### 1. Shade Determination

Clean the teeth prior to shade determination. The shade is selected with the tooth still moist.

### 2. Isolation

Appropriate isolation is required, preferably with a rubber dam (e.g. OptraDam® Plus).

### 3. Cavity preparation

The cavity is prepared according to the principles of the adhesive technique, i.e., by preserving as much of the tooth structure as possible. Do not prepare sharp internal edges or angles or additional undercuts in caries-free areas. The dimensions of the cavity are generally dictated by the extent of the caries or the size of the previous restorations. Bevel enamel edges of anterior teeth. In the posterior region, only the sharp enamel edges should be rounded. Do not prepare caries-free cervical defects. Clean them with pumice or a suitable cleaning paste (e.g. Proxyl RDA 36, medium) and with the help of rubber polishers or rotary brushes. Ablate sclerotised dentin around the tooth neck in a thin layer with a rotary instrument (diamond finisher, coarse disc). Subsequently, rinse the cavity with water to remove all the residue and dry it with water- and oil-free air.

### 4. Pulp protection / Base

Do not apply a base if an enamel-dentin bonding agent is used. Only cover very deep areas close to the pulp with a calcium hydroxide liner and subsequently use a pressure-resistant cement (e.g. Vivaglass® Liner). Do not cover the remaining cavity walls so that they can be used to generate a bond with an enamel-dentin adhesive.

### 5. Apply matrix / interdental wedge

Use a wrap around matrix for cavities affecting the proximal area or a sectional matrix and wedge it.

### 6. Conditioning / Application of the bonding agent

Condition and apply the bonding agent according to the instructions for use of the product in use. Ivoclar Vivadent recommends using the self-etching adhesive AdheSE® One F, for example.

### 7. Application of Compoglass Flow

- In order to achieve optimum results, Compoglass Flow should be applied in layers of max. 3 mm (light shades) or 2 mm (dark shades = 340/A4).
- For complete polymerization, sufficient exposure time is required. For the recommendations regarding exposure time and light intensity see Table 1.
- Hold the light emission window of the curing light as close to the restorative as possible.
- When using a metal matrix, additionally polymerize from the buccal or the lingual/palatal aspect after removing the matrix.

Program	Unit	Bluephase C8	Bluephase Style	Bluephase 20i
Turbo Program		–	–	2x5 s
High Power Program		40 s	20 s	20 s
Soft Start Program		40 s	–	20 s

Light intensity	Exposure Time
≥ 500 mW/cm <sup>2</sup>	40 s
≥ 1000 mW/cm <sup>2</sup>	20 s

Table 1

### 8. Finishing / Checking the occlusion / Polishing

After polymerization, remove excess material with suitable finishers (e.g. Astropol F) or fine diamonds. Remove proximal excess material with flexible discs, finishing strips or diamond-coated files. Check the occlusion and articulation and grind in the necessary adjustments to prevent premature contacts or undesired occlusal paths on the surface of the restoration. Use silicone polishers (e.g. OptraPol® Next Generation) as well as polishing discs and polishing strips to polish the restoration to a high gloss.

### Additional notes

1. Use Compoglass Flow at room temperature. Cold material may be difficult to dispense.
2.  For single use only. If Compoglass Flow is directly applied from the Cavifill in the mouth of the patient, the Cavifill should only be used for one patient for hygienic reasons (prevention of cross-contamination between patients).
3. Do not disinfect Cavifills with oxidizing disinfectants.

### Warning

Prevent any contact of uncured Compoglass Flow with skin / mucous membrane and eyes. Unpolymerized Compoglass Flow may cause slight irritation and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

### Storage and hygiene

- Do not use Compoglass Flow after the indicated date of expiration.
- Expiry date: see note on packaging.
- Storage temperature 2–28 °C (36–82 °F).
- Close Cavifills immediately after usage. Exposure to light causes premature polymerization.

### Keep out of the reach of children.

#### For dental use only.

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damage resulting from failure to observe the instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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## Ivoclar Vivadent Limited
## Ivoclar Vivadent, Inc.