ProBase[®] Hot

Instructions for Use

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Mode d'emploi

Istruzioni d'uso

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Bruksanvisning

Brugsanvisning

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Bruksanvisning

Productinformatie

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ivoclar[®]
vivadent[®]
technical

Description

ProBase® Hot sets quality standards for heat-curing denture base materials with regard to working properties, shape and shade stability and comfort of fit. The material is available in a variety of shades. Various methods of polymerization render the material suitable for a number of application possibilities.

Composition

Powder

Polymethyl methacrylate, softening agent, benzoyl peroxide, pigments

Liquid

Methyl methacrylate, dimethacrylate (linking agent), catalyst

Indication

- Complete dentures
- Partial dentures
- Combination dentures
- Relining

Contraindication

- Direct intraoral contact of unpolymerized material.
- If the patient is known to be allergic to any of the ingredients in ProBase Hot

Side effects

In individual cases, local allergic reactions to polymethyl methacrylate materials have been reported.

Application

Preparation

Isolate boiled-out, well-wetted plaster surfaces with two layers of Ivoclar Separating Fluid and allow to dry. To ensure an adequate bond with the denture base, well roughen the teeth and wet with monomer.

- **Isolate plaster surfaces twice.**
- **Invest wax model with plaster in the flask.**

Dosage

- **Ideal mixing ratio for one denture**
22.5 g polymer (powder) : 10 ml monomer (liquid)
- **With dosage system**
1 graduation mark polymer : 10 ml monomer

Integrated dosage system

The integrated dosage system ensures an ideal mixing ratio and, therefore, minimum polymerization shrinkage of ProBase Hot. The measuring cylinder for the polymer indicates the quantity of material used for one or two medium-sized dentures. The graduation on the measuring cylinder for the monomer is in millilitres. Use the appropriate graduation mark.

Mixing

Thoroughly mix polymer and monomer in the given ratio by means of a spatula. Mix thoroughly. Subsequently, leave the material to mature in the closed mixing cup at room temperature (23 °C / 73 °F) for approx. 8 to 10 minutes.

Working time

When the material has matured sufficiently and is no longer sticky, it can be worked for approx. 20 minutes at 23 °C/73 °F.

- **Thoroughly mix powder and liquid.**
- **Dough time and working time depend on the temperature.**

Pressing

Place a sufficient quantity of the resin dough in the hand-warm (approx. 40 °C / 104 °F), isolated flask halves. Carefully close flask, load with 80 bar pressure and fix with a clamp. **Maintain pressure.**

Polymerization

Heat-polymerization can be carried out in different ways:

Standard procedure (recommended method)

Place closed flask in cold water. Heat up to 100 °C / 212 °F and let boil for 45 minutes.

Alternative methods

- Place flask in cold water, heat up to 70 °C / 158 °F and leave it for 30 minutes. Then heat up to 100 °C / 212 °F and let boil for 30 minutes.
- Place flask in water of 70 °C / 158 °F and leave it for 60 minutes. Subsequently, heat up to 100 °C / 212 °F and let boil for 30 minutes.
- Place flask in boiling water. Bring the water to the boil again and then let boil for 40 minutes. This procedure is only suitable for medium-sized dentures.
- Place flask in cold water, heat up to 80 °C / 176 °F and polymerize for 10 hours. Switch off heat source and leave the flask to cool in the same water bath over night.
- Polymerize the contents of the flask for 10 hours at 80 °C / 176 °F in the drying cabinet.

Residual monomer content can be reduced by increasing the polymerization temperature and by prolonging the polymerization length. We recommend using the standard procedure to keep the residual monomer content at minimum levels. The residual monomer content after polymerizing the material according to the standard procedure is <2.2%.

Cooling

Let the flask cool at room temperature for 30 minutes. Subsequently, completely cool the flask with cold water.

Deflasking and finishing

Open the completely cooled flask and remove plaster. Check occlusion of the denture. Subsequently, finish and polish in the usual manner.

Repair and correction possibilities of ProBase Hot

Corrections and repairs can be carried out with the cold-curing ProBase Cold material, using the pouring technique. Thoroughly roughen the corresponding surfaces and wet them with monomer.

Warnings

- The monomer contains methyl methacrylate (MMA).
- Methyl methacrylate is easily flammable and irritating (flash point +10 °C / 50 °F).
- Irritating to eyes, skin, and respiratory system.
- May cause sensitization by skin contact.
- Avoid contact of the skin with monomer or uncured material. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.
- Do not breathe vapour.
- Keep away from sources of ignition – no smoking.
- Do not empty into drains.
- Take precautionary measures against static discharges.

Storage

- Store material in a cool, dark, well-ventilated place. Storage temperature: 2–28 °C (36–82 °F).
- Do not use the materials after the indicated date of expiration.
- Keep out of the reach of children.

Date information prepared: 08/2012

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Pulver:

Polymethylmethacrylat, Weichmacher, Benzoylperoxid, Pigmente

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.