Stratos® 300

Operating Instructions

Page 3

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Seite 23

Mode d'emploi

Page 43

Istruzioni d'uso

Pagina 63

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Página 83

ivoclar
vivadent®

Stratos® 300

	Page
Views of the Unit and Parts List	4
– Description of components	4
– Accessories	6
– Accessories for average model transfer	6
– Accessories for personalized model transfer	6
1 Introduction / Signs and Symbols	8
1.1 Preface	8
1.2 Signs and symbols	8
1.3 Information on the Operating Instructions	8
2 Safety First	9
2.1 Field of application	9
2.2 Health and safety instructions	9
3 Product Description	9
3.1 Functional description	9
3.2 Indication, contraindication	9
4 Installation and Initial Start-up	10
4.1 Unpacking and checking the contents	10
4.2 Assembly and initial set-up	10
5 Handling and Operation	10
5.1 Centric position	10
5.2 Protrusion	11
5.3 Lateral and Bennett movement	12
5.4 Immediate side-shift movement	13
5.5 Retrusion	13
6 Model Orientation in the Stratos 300	14
6.1 Average orientation of dentulous and edentulous cases with a rubber band	14
6.2 Average orientation of dentulous cases with the set-up table	14
6.3 Average orientation of edentulous cases with the horizontal guide	15
6.4 Average orientation of dentulous cases with the horizontal guide	15
6.5 Personalized model transfer with the registration joint holder	16
6.6 Personalized model transfer with the UTS transferbow	17
6.7 Personalized setup of anterior guidance	17
7 Maintenance and Cleaning	18
7.1 Monitoring and maintenance	18
7.2 Cleaning	18
8 What If ...	19
8.1 Technical malfunctions	19
8.2 Repairs	19
9 Product Specifications	20
9.1 Delivery form	20
9.2 Technical data	20
10 Miscellaneous	21
10.1 Tips on the coordination of articulators	21

Parts List

Stratos® 300

- 1.1 Incisal plate
- 1.2 Incisal guide rod
- 1.3 Incisal indicator
- 1.4 Guide rod holder
- 1.5 Guide rod retaining screw
- 1.6 Guide rod initial set position (red)
- 1.7 FI knurled screw
- 1.8 Hole for incisal indicator

- 2.1 GI knurled screw
- 2.2 Lower frame
- 2.3 Upper frame
- 2.4 Incline supports
- 2.5 Notch for rubber band
- 2.6 Non-slip base
- 2.7 Joint bolt

- 3.1 Bennett flap
- 3.2 B knurled screw
- 3.3 Bennett angle scale

- 4.1 Centric fixation lock (activable centric locking system)
- 4.2 Centric rubber band
- 4.3 Centric rubber band holder

- 5.1 P knurled screw
- 5.2 P angle indicator
- 5.3 Protrusion angle scale

- 6.1 Locating hole for transferbow

- 7.1 Holding device for incline support holder
- 7.2 Incline support holder

- 8.0 Articulator number






- 9.0 PR scale
- 9.1 Screw for PR stop
- 9.2 PR stop

- 10.0 ISS screw
- 10.1 ISS screw thread
- 10.2 ISS scale

- 11.0 Opening stop


- 12.0 Retrusion screw

- 33 Retentive base block


Accessories for the Stratos® 300

- 30 15° incisal plate made of plastic
- 31 30° incisal plate made of plastic
- 32 Incisal plate made of MMA-soluble material for personalized anterior guidance

- 33 Retentive base block
- 34 Plaster protection plate
- 35 Magnetic base block
- 36 Collar
- 37 Retention disk


Accessories for average model transfer

40 Instrument carrier for horizontal guide, set-up table, 2-D setting up template and bite fork support

41 Horizontal guide
41.1 Symphysis fork

42 2-D setting up template

43 Set-up table

Accessories for personalized model transfer

50 Plane indicator

51 Bite fork support

52 Registration joint holder (CP)

53 3-D setting up template


54 FH registration joint holder

60 Adjustable support pins (Type 2) for the UTS transferbow

61 UTS transferbow system

1. Introduction / Signs and Symbols

1.1 1.1 Preface

Dear Customer

Thank you for having purchased the Stratos 300 articulator. The Stratos 300 is an advanced technology product that is characterized by high quality and precision.

The Stratos 300 has been designed according to the latest industry standards. Inappropriate use may be hazardous. Please observe the relevant safety instructions and do read the Operating Instructions carefully.

We wish you much success and satisfaction with the Stratos 300.

1.2 Signs and symbols

The signs and symbols in these Operating Instructions facilitate the finding of important information. They have the following meanings:

Risks and dangers

Important information

Contraindication

1.3 Information on the Operating Instructions

These Instructions apply to the following apparatus: Stratos 300.
Target group: dentists, dental technicians, dental lab professionals

These Operating Instructions facilitate the correct, safe and economic use of the Stratos 300 articulator.

If you have lost the Operating Instructions, you may order an extra copy at a nominal fee from your local Ivoclar Vivadent Service Center.

2. Safety First

This chapter is particularly important for staff who work with the apparatus or who have to carry out maintenance or repair work.

2.1 Field of application


The Stratos 300 must only be used for the indications stipulated in Chapter 3. Further instructions to assure the proper use of the Stratos 300 articulator are as follows:

- The instructions, regulations and notes in these Operating Instructions must be observed.
- The unit must be properly maintained (see Chapter 7).


2.2 Health and safety instructions

If the centric rubber bands are removed from their retainers and the centric fixation lock is open, the upper and lower frames may be accidentally separated.

Risk of crushing between centric locking catch and incline support.

3. Product Description

3.1 Functional description

The Stratos 300 is an individually adjustable articulator built according to the Camper's plane (CP).

The device is suitable to simulate different human jaw movements. The protrusion angle and the Bennett angle can be set individually.

Furthermore, a retrusive movement of 35° as well as an "immediate side shift" are integrated in the Stratos 300.

An activable centric locking system permits the articulator to be set to a precise, reproducible initial set position and facilitates separation and connection of the upper and lower frames. Even when the centric lock is open, the upper and lower frames are still connected. Only after removing the centric rubber bands from their retainers can the two parts be separated.

The locating holes for transferbows, automatic centric return and a non-slip base are all standard equipment. In addition, the Stratos 300 features a dirt-and-wear-resistant finish.

A wide range of accessories permits the Stratos 300 to meet the specific needs and requirements of every user. Therefore, the apparatus can be used in the field of removable as well as fixed denture prosthetics.

3.2 Indication and contraindication

Indication

The articulator is suitable for the spatial fixation of models for lab-fabricated reconstructions. The basic equipment permits simulation of average movement patterns of the human temporomandibular joint.

Contraindication

No contraindications are known to date provided that the articulator is used strictly according to the Operating Instructions.

4. Installation and Initial Start-up

4.1 Unpacking and checking the contents

Remove the components of the Stratos 300 from their packaging and check the delivery for completeness.

Stratos 300 basic model:


- 1 x Stratos 300
- 1x Incisal plate 0°
- 1x Incisal indicator
- 2x Retentive base blocks
- 1x Plaster protection plate
- 1x Incline support holder

If certain parts are missing or damaged, contact your local Ivoclar Vivadent Customer Service. We recommend keeping the original packaging for possible future transportation purposes.

4.2 Assembly and initial set-up

Incline support holder

The incline support holder enables ergonomic handling of the articulator. If mounted, it maintains the articulator at approx. a 45° angle. This position permits a superb overview of the model mounted in the articulator.


Insert the incline support holder (7.2) into the holding device (7.1) for the incline support holder. If necessary, the incline support holder can be removed from the apparatus.

Secure the base block (33) to the articulator by means of the knurled screw (2.1).

5. Handling and Operation


5.1 Centric position


The Stratos 300 features a precise centric locking system that can be activated in the following positions:

5.1.1 Centric fixation (4.1) open

This position permits the simulation of jaw movements.


If the rubber bands (4.2) are removed, the two frames can be separated.


5.1.2 Centric fixation lock (4.1) closed

In this position, the two frames of the Stratos 300 can no longer be separated and the articulator is fixed in the centric position. The centric fixation can be opened by applying slight pressure. This greatly facilitates handling for the user.

5.1.3 Centric fixation (4.1) locked


To lock the centric fixation completely, e.g. if the articulator has to be transported, or if no movements are desired, pull the joint axis into the precisely defined, final position by turning the centric screw to the right.

Loosen the centric fixation before carrying out masticatory movements.


5.2 Protrusion

5.2.1 Protrusive movements can be simulated if the centric fixation (4.1) is open, according to the schematic below. The angle of the protrusion path can be adjusted to be more acute or flat (0–60°) by setting the respective protrusion angle.


Protrusion angle:

By loosening the P knurled screw (5.1), the protrusion angle may be adjusted and fixed with the P angle indicator (5.2).



The set protrusion angle can be read off at the top edge of the P angle indicator (5.2).

5.2.2 PR stop (9.2) for centric shift

The PR stop permits each condyle to be incrementally advanced from 0 to 4 mm into a protrusive position.

Loosen the screw for the PR stop (9.1) and move the PR stop (9.2) into the desired position on the PR scale (9.0). By fastening the screw for the PR stop (9.1), the position can be fixed.


Important:

In this case, the centric fixation cannot be used.

5.3 Lateral and Bennett movement


5.3.1 Lateral and Bennett movements

The lateral movements can be simulated when the centric fixation (4.1) is open, according to the schematic below. To perform this movement, press with your thumb unilaterally on the joint of the corresponding side.


Bennett angle

By loosening the B knurled screw (3.2), the Bennett angle can be read off, adjusted and fixed with the Bennett angle scale (3.3).


5.4 Immediate side-shift movement

The ISS screw permits the fixation of the immediate side-shift.


By turning the ISS screw (10.0), the desired distance from the Bennett flap (3.1) to the joint bolt can be fixed. The distance can be continuously fixed between 0 and 1.5 mm and read off the ISS scale (10.2).


5.5 Retrusion

The retrusion screw (9.1) permits the fixation of the retrusive movements.

To simulate retrusive movements, the centric fixation lock (4.1) must be open and the screw to fix the protrusive movement (9.1) loosened. Moreover, the Bennett flap (3.1) must be set to 0°. This allows retrusive movements to be carried out properly.



The distance may be continuously set from 0 to 2 mm with the retrusion screw (12.0) and read off the PR scale (9.0).

6. Model Orientation in the Stratos 300

6.1 Average model orientation of dentulous and edentulous cases with a rubber band

6.1.1 Return articulator to the initial set position

- Check if Bennett inserts (3.1) are secured in place.
- Mount centric rubber bands (4.2).
- Bring centric fixation lock (4.1) into the upper position and secure in place.
- Secure incisal guide rod (1.2) in the initial position (red mark 1.6) in the guide rod holder.


6.1.2. Incisal indicator (1.3)

Slide the incisal indicator completely into the incisal guide rod.

Insert the longer part of the incisal indicator (1.3) into the upper hole (1.8) of the incisal guide rod (note notch in the rod).

6.1.3 Next, attach the thin rubber band to the incline supports (2.4), notches (2.5), and under the longer part of the incisal indicator (1.3).

6.1.4 Adjust occlusal plane according to the rubber band. Adjust mesial contact area either between the mandibular central incisors or on the wax bite block, according to the tip of the incisal indicator (1.3).

Example: edentulous case

6.1.5 For best results, mount mandibular model on modelling material and adjust. After that, cast maxillary model.

Example: dentulous case

6.2 Average orientation of dentulous cases with the set-up table

6.2.1 Return the articulator to the initial set position

- See 6.1.1 for description

Remove the incisal indicator (1.3).

6.2.2 Mount instrument carrier to the upper frame (2.3) and fasten it with the GI knurled screw (2.1).

6.2.3 Insert set-up table completely into the instrument carrier and secure it.

6.2.4 Place articulator with the upper frame (2.3) facing down on the table. Orientate the dentulous mandibular model according to the mesial contact area and symmetry in the molar region. If necessary, use modelling material to hold it in place.

6.2.5 Apply plaster on model and base block and close articulator slowly.

6.2.6 Then, orientate maxillary model in the usual manner.

6.3 Average orientation of edentulous cases with the horizontal guide

6.3.1 Return articulator to the initial set position


- See 6.1.1 for description

Remove incisal indicator (1.3).

6.3.2 Mark and bisect the Trigonum retromolare on the mandibular model.

Bisect the distance between the upper and lower mucolabial folds and set the symphysis fork (41.1) of the horizontal plane on the resulting value.

Using the resulting values, orientate the mandibular model to the horizontal guide and hold it in place, e.g. with a rubber band.

6.3.3 Mount instrument carrier to the upper frame (2.3)

Fasten horizontal guide in place by means of the knurled screw of the instrument carrier.

Insert horizontal guide completely into the instrument carrier. After that, apply some plaster on the model and base block and close articulator slowly.

6.3.4 Next, orientate maxillary model as usual.

6.4 Average orientation of dentulous cases with the horizontal guide

6.4.1 Return articulator to the initial set position

- See 6.1.1 for description

Remove incisal indicator (1.3).

6.4.2 Transfer maxillary mesial contact area vertically to the mandibular model. Mark the tips of both distobuccal cusps of the lower second molars. If these molars are missing, the first molars may be used instead.

- Put mandibular mesial contact area behind the incisal tip of the symphysis fork.
- Adjust horizontal guide in such a way that the rear edges barely cover the marked cusps of the second molars and align symmetry.
- The retromolar pads are used as reference points for free-end surfaces, similarly to the Trigonum retromolare in edentulous cases.

Subsequently, attach horizontal guide to the model with e.g. wax or a rubber band.

6.4.3 Mount instrument carrier to the upper frame (2.3).

Fasten horizontal guide in place by means of the knurled screw of the instrument carrier.

Insert horizontal guide completely into the instrument carrier.

After that, apply some plaster on the model and base block and close articulator slowly.

6.4.4 Next, orientate maxillary model as usual.

6.5 Personalized model transfer with the registration joint holder

6.5.1 Return articulator to the initial set position.


– See 6.1.1 for description

Remove incisal indicator (1.2).

6.5.2 Mount the plane indicator (50) instead of the incisal guide rod (1.2).

6.5.3 Mount bite fork support (51) to the lower frame (2.2) by means of the instrument carrier (40).

6.5.4 Insert registration joint holder (52) instead of the incisal plate (1.1) and secure it.

6.5.5 Mount UTS registration to the registration joint holder and secure bite fork with the bite fork support.

Example: Centric Tray

Example: UTS bite fork (dentulous case)


6.5.6 Mount maxillary model to the bite fork and cast. Close articulator until the plane indicator lies on the registration joint holder.

Example: UTS bite fork (edentulous case)

FH registration joint holder (54)

The procedure is the same as with the CP registration joint holder. The plane indicator (50) lies on the relevant support.

6.5.7 Subsequently, orientate mandibular model in the usual manner.

6.6 Personalized model transfer with the UTS transferbow

6.6.1 Return articulator to the initial set position

- See 6.1.1 for description


Remove incisal guide rod (1.2).

6.6.2 Mount the plane indicator (50) instead of the incisal guide rod (1.2).

6.6.3 Mount bite fork support to the lower frame (2.2) by means of the instrument carrier.

6.6.4

- Remove nose rest from the transferbow and insert support pins (Type II) (60) from above.

- Insert the UTS into the locating holes for transferbows (6.1).
- Using the support pins, adjust the transferbow so that it is parallel to the table top.
- Finally, secure bite fork with the bite fork support.

6.6.5 Mount maxillary model to the bite fork and cast. Close the articulator until the plane indicator rests in the notch in the cross-beam of the support pin.

If Type-I support pins are used, close articulator until the plane indicator touches the SME indicator of the transferbow.


6.6.6 Then, orientate mandibular model in the usual manner.

6.7 Personalized setup of anterior guidance

For the fabrication of personalized incisal plates, Ivoclar Vivadent offers a special transparent incisal plate. The plate is made of a MMA-soluble plastic that bonds with the set-up resin (SR Ivolen). Anterior guidance plates thus fabricated can always be exactly repositioned in the Stratos 300.

Procedure

First, mount the transparent 0° incisal plate. Then, position the models of the patient for which the anterior guidance should be formed in the Stratos 300.

Subsequently, mix appropriate resin and spread it on the incisal plate. Finally, simulate the registered masticatory movements (protrusion, retrusion, laterotrusion) with the upper frame of the Stratos 300. The movement patterns will thus be scratched into the resin by the incisal guide rod and recorded after polymerization of the resin.

The forming of the natural facets should be done from excursion to centric. Otherwise, the resin might be forced out of the incisal plate. Make sure that no vertical increase of occlusion occurs.

7. Maintenance and Cleaning

This chapter describes the user maintenance and cleaning procedure for the Stratos 300. Only those tasks that can be carried out by qualified dental lab experts are listed. All other tasks must be performed by qualified service personnel at a certified Ivoclar Vivadent Service Center.

7.1 Monitoring and maintenance

The time for these maintenance procedures depends on the frequency of use and the working habits of the user. For these reasons, the recommended times are only approximates.

What	Part	Frequency
Check centric rubber band for damage and cracks and replace, if necessary.	Centric rubber band	Monthly or as necessary
Check if knurled screws are dirty.	Knurled screw	Weekly
Check if joint bolt is sufficiently lubricated.	Joint bolt	Monthly or as necessary
Check if incline support holder and holding device are dirty.	Holding device for incline support holder	Weekly or as necessary

7.2 Cleaning

What	Frequency	Cleaning agent/Measure
Knurled screws contaminated with wax or plaster	Weekly or as necessary	Rinse with warm water
Joint bolt contaminated with dust	Weekly or as necessary	Rinse with warm water

Protect the Stratos 300 from acids and solvents (e.g. MMA) to prevent the finish from being damaged.

8. What If ...?

This chapter will help you to recognize malfunctions and take appropriate measures or, if possible, to perform some repairs.

8.1 Technical malfunctions

Error	Cause/Description	Action
Upper frame exceeds opening stop and slips backwards.	The opening stop is not in place or defective.	Replace opening stop or mount it properly.
Articulator cannot be opened when the centric fixation lock is closed.	Centric fixation has been fixed too tightly.	Fix centric fixation less tightly.
Joint bolt is damaged or grooved.	The articulator was forcibly opened and closed without the centric fixation being loosened.	Loosen centric fixation lock before opening and closing the articulator.
The surface of the articulator is damaged or etched.	Surface was cleaned with acid, solvent or lye.	Do not use acid, solvents or lye to clean the apparatus.
The surface of the apparatus is scratched or damaged.	Surface was scratched with a sharp instrument.	Do not scratch the surface with a sharp instrument.

8.2 Repairs

Repairs may only be carried out by a certified Ivoclar Vivadent Service Center. Please refer to the addresses mentioned in Chapter 10. If repairs during the warranty period are not carried out by a certified Ivoclar Vivadent Service Center, the warranty will expire immediately.

9. Product Specifications

9.1 Delivery forms

Stratos 300 basic model

- 1 Stratos 300
- 1 Incisal plate 0°
- 1 Incisal indicator
- 2 Retentive base blocks
- 1 Plaster protection plate
- 1 Incline support holder

"Average" Accessories Assortment

- 1 Instrument carrier
- 1 Setup table
- 1 Horizontal guide
- 1 2-D setting up template

"Personalized" Accessories Assortment

- 1 Complete set of joint inserts
- 1 Incisal plate 15° and 30° each, as well as 1 personalized plate
- 1 Instrument carrier
- 1 Plane indicator
- 1 Registration joint holder
- 1 Bite fork support

Separately available:

- Incisal plate 15° or 30°
- Incisal plate for personalized anterior guidance, package of 5 plates
- Instrument carrier
- Plane indicator
- Horizontal guide
- 2-D setting up template
- 3-D setting up template
- Set-up table
- Adjustable support pins (Type II)
- Registration joint holder
- Bite fork support
- Magnetic base blocks, package of 2 or 10 blocks
- Retentive base blocks, package of 10 or 50 blocks
- Plaster protection plate, package of 5 plates
- FH registration joint holder
- Incline support holder
- Retention disk for magnetic base block

9.2 Technical data

- Bonnwill triangle 108 mm
- Balkwill angle 15°
- Working height 118 mm
- Retrusion path angle 35°
- Protrusion screw 0-4 mm
- ISS: 0-1.5 mm
- Protrusion angle: 0°-60°
- Bennett angle: 0°-30°
- Exchangeable, colour-coded 0° incisal plate. Other angles are available as accessories.
- Weight: 950 g
- Colour: white (RAL 9016), incline support: gold

The delivery forms may vary from country to country.

10. Miscellaneous

10.1 Tips on the coordination of articulators

The companies below offer split-cast systems for the Stratos 300.

These systems permit users to coordinate their Stratos 300 articulators. For further information please contact:

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Please note that the standard accessories can no longer be used after the apparatus have been coordinated.

This apparatus has been developed for use in dentistry. Setup and operation should be carried out strictly according to the Operating Instructions. Liability cannot be accepted for damages resulting from misuse or failure to observe the Instructions. The user is solely responsible for testing the apparatus for its suitability for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent S.A.
Ivoclar Vivadent AB
Ivoclar Vivadent, Inc.