Bluephase® C8

Licence to cure

LED for every use

EN Instructions for Use

Page 6

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Pagina 42

Appendix

Page 78

DE Gebrauchsinformation

Seite 18

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Página 54

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Page 30

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Página 66

CE 0123

For dental use only!

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist


Made in Austria

Ivoclar Vivadent AG, FL-9494 Schaan / Liechtenstein

ivoclar
vivadent
clinical

Bluephase® C8

LED for every use

Technical data

Wavelength range:	385–515 nm
Light intensity:	800 mW/cm ² ± 10%
Curing modes:	High Power, Soft Start, Low Power
Exposure time:	20 sec 15 sec (Tetric EvoCeram / IPS Empress Direct)
Light probe:	10 mm
Power supply:	Mains operated
Operating voltage:	100–240 VAC / 50–60 Hz
Warranty	3 years

English	Deutsch	Français	Italiano	Español	Português
List of parts	Teileverzeichnis	Liste des pièces détachées	Elenco componenti	Despiece	Lista de partes
① Handpiece	Handstück	Pièce à main	Manipolo	Pieza de mano	Peça de mão
② Start button	Starttaster	Bouton de mise en marche	Tasto di accensione	Interruptor de puesta en marcha	Botão de ligar
③ Display	Display	Affichage	Display	Pantalla	Display
④ Program selection button	Programmwahl-taster	Bouton de sélection des programmes	Tasto di selezione programma	Tecla de selección de programa	Tecla de seleção do programa
⑤ Time selection button	Zeitwahl-taster	Bouton de sélection du temps	Tasto di selezione tempo	Tecla de selección de tiempo	Tecla de seleção do tempo
⑥ Volume button	Lautstärketaster	Bouton de réglage du volume	Tasto volume	Tecla de volumen	Tecla de volume
⑦ Handpiece holder	Handstückhalter	Socle pour pièce à main	Base del manipolo	Soporte de la Pieza de mano	Suporte da peça de mão
⑧ Power cord	Netzkabel	Cordon d'alimentation	Cavo di alimentazione	Cable eléctrico	Cabo de força
⑨ Power pack	Netzgerät	Power pack	Alimentatore	Unidad de alimentación	Unidade de alimentação
⑩ Light probe 10 mm	Lichtleiter 10 mm	Conducteur de lumière 10 mm	Conduttore ottico 10 mm	Conducto de luz 10 mm	Conductor de luz 10 mm
⑪ Anti-glare cone	Blendschutzkegel	Embout de protection oculaire	Coni antiriflesso	Conos antideslumbrantes	Cone anti-ofuscante

English	Deutsch	Français	Italiano	Español	Português
Indicators on the handpiece	Anzeigen des Handstücks	Indicateurs sur la pièce à main	Indicatori sul manipolo	Indicadores en la pieza de mano	Indicações da peça de mão
HIGH LOW SOFT	Curing program and operating state	Belichtungsprogramm bzw. Betriebszustand	Programme de polymérisation, temps de polymérisation, état de marche	Programma di polimerizzazione e stato operativo	Programa de polimerización y estado de funcionamiento
	Selected indications and curing time	Angewählte Indikation bzw. Belichtungsdauer	Indications spécifiques du programme de polymérisation sélectionné	Indicazione selezionata e tempo di polimerizzazione	Indicaciones seleccionadas y tiempo de polimerización
	Correct power supply	Netzanschluss korrekt	Alimentation électrique correcte	Alimentación eléctrica correcta	Correto suprimento de força
	Display is in Stand-by mode	Display ist im Stand-by Modus	L'écran est en mode Stand-by	Display en modo Stand-by	Display está no modo „Stand-by“.

Signs and symbols	Zeichen-erklärung	Signes et symboles	Segni e simboli	Señales y símbolos	Sinais e símbolos
	Contraindication	Nicht zulässige Anwendung	Ne pas utiliser	Uso no permitido	Uso não autorizado
	Double insulation (Apparatus complies with safety class II)	Doppelt isoliert (Gerät der Schutzklasse II)	Double protection (l'appareil répond aux normes de sécurité de classe II)	Doppio isolamento (apparecchio conforme alla Classe di sicurezza II)	Doble aislamiento (El dispositivo cumple con la Clase de Seguridad II)
	Protection against electrical shock (BF type apparatus)	Schutz gegen elektrischen Schlag (Gerätetyp BF)	Protection contre les décharges électriques (type d'appareil BF)	Protezione da scossa elettrica (apparecchio di tipo BF)	Protección frente al choque eléctrico (aparato tipo BF)
	Observe Instructions for Use	Gebrauchsinformation beachten	Respecter le mode d'emploi	Attenersi alle istruzioni d'uso	Observar las instrucciones de uso
	Observe Instructions for Use	Gebrauchsinformation beachten	Respecter le mode d'emploi	Attenersi alle istruzioni d'uso	Observar las instrucciones de uso
	Caution	Vorsicht	Mise en garde	Attenzione	Precaución
	The apparatus must not be disposed of in the normal domestic waste	Das Gerät darf nicht über den normalen Hausmüll entsorgt werden	Ne pas jeter l'appareil avec les ordures ménagères	L'apparecchio non può essere smaltito come normale rifiuto urbano	El aparato no se puede eliminar con la basura doméstica normal
	Recyclable	Recyclbar	Recyclable	Riciclabile	Reciclable

Bluephase® C8


Light intensity	800 mW/cm ² ± 10%
Exposure time for composites	20 s
Exposure time for Tetric EvoCeram / IPS Empress Direct	15 s

HIGH POWER Program

Consistently high light intensity for the polymerization of restorative and cementation materials for direct and indirect restorations.

Filling materials	Curing time
Composite	
• 2 mm ¹⁾	
IPS Empress Direct / IPS Empress Direct Flow / Tetric EvoCeram / Tetric EvoFlow / Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	15 s
Heliomolar / Heliomolar HB / Heliomolar Flow / Others	20 s
• 4 mm ²⁾	
Tetric EvoCeram Bulk Fill	15 s
Tetric Basic White	20 s
Compomer ³⁾	
Compoglass F / Compoglass Flow	40 s
Indirect restorations / Luting materials	
Variolink II Base ⁴⁾ / Variolink Veneer / Variolink N Base ⁴⁾ / Variolink N Clear Veneer / Variolink Ultra ³⁾ / Dual Cement ⁵⁾ / Variolink II ⁵⁾	per mm ceramic 10 s per segment
Miscellaneous	
Helioseal / Helioseal F / Helioseal Clear	20 s
Monopaque	30 s
MultiCore Flow / Multicore HB	40 s
Heliosit Orthodontic	20 s
Telio Add-On Flow	20 s
Telio Stains	15 s
IPS Empress Direct Color	15 s
IPS Empress Direct Opaque	30 s

LOW POWER Program

Reduced light intensity with reduced heat development for the polymerization of adhesives, liners, and restorative materials in areas near the pulp when restoring Class V cavities.

Adhesive	Curing time
AdheSE / AdheSE One F Excite / Excite DSC Tetric N-Bond / Tetric N-Bond Self-Etch Heliobond (Syntac)	10 s
Temporary materials	
Telio CS Inlay/Onlay	10 s
System.inlay/onlay	10 s
Fermit / Fermit N	10 s
Telio CS Link / System.link	20 s per segment
Miscellaneous	
Heliosit Orthodontic	20 s
Vivaglass Liner	20 s

SOFT START Program

Step-by-step increase of the light intensity with reduced shrinkage stress and reduced heat development for the polymerization of restorative materials.

Filling materials	Curing time
Composite	
• 2 mm ¹⁾	
IPS Empress Direct / IPS Empress Direct Flow / Tetric EvoCeram / Tetric EvoFlow / Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	20 s
Heliomolar / Heliomolar HB / Heliomolar Flow / Others	25 s
• 4 mm ²⁾	
Tetric EvoCeram Bulk Fill	20 s
Tetric Basic White	25 s
Compomer ³⁾	
Compoglass F / Compoglass Flow	40 s

- 1) Applies to a maximum layer thickness of 2 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case, e.g. with dentin shades)
 2) Applies to a maximum layer thickness of 4 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case, e.g. with dentin shades)
 3) Applies to a maximum layer thickness of 3 mm 4) Applies to light-curing (use of the base past only) 5) Applies to dual-curing

Dear Customer,

Optimum polymerization is an important requirement for all light cured materials in order to consistently produce high quality restorations. The polymerization light selected also plays a decisive role in this respect. Therefore, we would like to thank you for having purchased Bluephase® C8. Bluephase C8 is a high-quality medical device which has been designed according to the latest standard of science and technology in compliance with the relevant industry standards.

These Instructions for Use will help you safely to start up the device, make full use of its capabilities and ensure a long service life.

Should you have any further questions, please do not hesitate to contact us (see addresses on the back cover).

Your Ivoclar Vivadent Team

Table of Contents

Safety	
Intended use	8
Indication	8
Safety notes	9
Start-up	
Handpiece holder	10
Handpiece	11
Operation	
Selecting the curing program and the curing time	12
Cure Memory function	14
Start	
Acoustic signals	14
Light intensity	14
Measuring the light intensity	14
Maintenance and cleaning	
Housing	15
Light probes	15
Disposal	15
What if ...	16
Warranty / Procedure in case of repair	16
Product specifications	17

Safety

Intended use


Bluephase C8 is an LED polymerization light that produces energy-rich blue light. It is used for the polymerization of light-curing dental materials immediately at the dental unit. The intended use also includes the observation of the notes and regulations in these Operating Instructions.

Indication

With its "Polywave[®]" broadband spectrum, Bluephase C8 is suitable for the polymerization of all light-curing dental materials curing in the wavelength range of 385–515 nm. These materials include restoratives, bonding agents/adhesives, bases, liners, fissure sealants, temporary materials as well as luting materials for brackets and dental-lab materials such as those used in the fabrication of ceramic inlays.

Safety notes

The Bluephase C8 is an electronic device and a medical product which is subject to IEC 60601-1 (EN 60601-1) and EMC directives, as well as the 93/42/EEC Medical Device Directive. The appliance complies with the relevant EU regulations and is classified as an LED Class 2 product.

CE 0123

The curing light has been shipped from the manufacturer in a safe and technically sound condition. In order to maintain this condition and to ensure risk-free operation, the notes and regulations in these Instructions for Use have to be observed. To prevent damage to equipment and risks for patients, users and third parties, the following safety instructions have to be observed.

Contraindications

Materials, the polymerization of which is activated outside the wavelength range of 385–515 nm (no materials known to date). If you are not sure about certain products, please ask the manufacturer of the corresponding material.

Do not use the appliance near flammable or combustible substances.

Portable and mobile high-frequency communication devices may interfere with medical equipment. The use of mobile phones during operation is not allowed.

Usage and liability

- Bluephase C8 must only be employed for the intended use. Any other uses are contraindicated. Liability cannot be accepted for damage resulting from misuse or failure to observe the Operating Instructions.
- The user is responsible for testing Bluephase C8 for its use and suitability for the intended purposes. This is particularly important if other equipment is used in the immediate vicinity of Bluephase at the same time.
- Use only original spare parts and accessories from Ivoclar Vivadent (see Accessories). The manufacturer does not accept any liability for damage resulting from the use of other spare parts or accessories.
- The light probe is an applied part and may warm up to a maximum of 45 °C (113 °F) at the interface to the hand-piece during operation.

Operating voltage

Before switching on, make sure that

- a) the voltage indicated on the rating plate complies with the local power supply and
- b) the unit has acquired the ambient temperature.

Assumption of impaired safety

If it has to be assumed that safe operation is no longer possible, the power must be disconnected to avoid accidental operation. This may be the case, for example, if the apparatus is visibly damaged or no longer works correctly. The apparatus is only completely disconnected from the power supply if the power cord has been pulled out of the plug-in socket.


Eye protection

Direct or indirect exposure of the eyes must be prevented. Prolonged exposure to the light is unpleasant to the eyes and may result in injury. Therefore, using the supplied anti-glare cones is recommended. Individuals who are generally sensitive to light, who take photosensitizing drugs, have undergone eye surgery, or people who work with the apparatus or in its vicinity for long periods of time should not be exposed to the light of this device. Therefore, they should wear protective goggles (orange) that absorb light below 515 nm.

Heat development (burn hazard)

As it is the case with all high-performance lights, the high light intensity results in a certain heat development. Prolonged exposure of areas near the pulp and soft tissues may result in irreversible or reversible damage. Therefore, the recommended curing times have to be observed. Uninterrupted curing times of more than 40 seconds on the same tooth surface, as well as direct contact with the gingiva, oral mucous membrane or skin have to be prevented. Polymerize indirect restorations at intermittent intervals of 20 seconds each or use external cooling with an air stream.

Start-up

Delivery form

Check the delivery for completeness and any possible transportation damage (see Delivery form). If components are missing or damaged, immediately contact your dealer or your Service Center.

Before you switch on the device, make sure that the voltage mentioned on the rating plate complies with your local power supply. The rating plate is attached to the underside of the power pack.


Connect the power cord with the power supply and the power pack.

Handpiece

Before you insert the light probe, remove the protective foil on the corresponding opening of the handpiece.

Slightly turn the light probe while attaching it to the handpiece.

After that, mount the anti-glare cone on the light probe.


Operation

Disinfect or autoclave contaminated surfaces of the curing light as well as of the light probes and anti-glare cones before each use. Furthermore, make sure that the stipulated light intensity permits adequate polymerization. For this purpose, check the light probe for contamination and damage as well as the light intensity at regular intervals (e.g. with the Ivoclar Vivadent radiometer "Bluephase Meter").

Selecting the curing program and the curing time

The curing program and the curing time can be individually set. Bluephase C8 is equipped with the following 3 curing programs for the different indications. The desired curing program is selected with the program selection buttons. The display changes accordingly (see Indicators on the handpiece). The device comes equipped with the following preset programs:

Factory settings

HIGH
(High Power)

15 seconds

LOW
(Low Power)

10 seconds

SOFT
(Soft start)

20 seconds

The intended curing time is selected using the time selection buttons. Users may choose between 5, 10, 15, 20 and 30 seconds.



Observe the Instructions for Use of the material applied when selecting the curing time.

The curing recommendations for composite materials apply to all shades and, if not mentioned otherwise in the Instructions for Use, to a maximum layer thickness of 2 mm. Generally, these recommendations apply to situations where the emission window of the light probe is placed directly over the material to be polymerized. Increasing the distance between the light source and the material will require the curing time to be extended accordingly. For instance, if the distance to the material is 9 mm, the effective light output is reduced by approx. 50%. In this case, the recommended curing time has to be doubled.

Bluephase® C8


Light intensity	800 mW/cm ² ± 10%
Exposure time for composites	20 s
Exposure time for Tetric EvoCeram / IPS Empress Direct	15 s

HIGH POWER Program

Consistently high light intensity for the polymerization of restorative and cementation materials for direct and indirect restorations.


Filling materials	Curing time
Composite	
• 2 mm ¹⁾	
IPS Empress Direct / IPS Empress Direct Flow / Tetric EvoCeram / Tetric EvoFlow / Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	15 s
• 4 mm ²⁾	
Heliomolar / Heliomolar HB / Heliomolar Flow / Others	20 s
• 4 mm ²⁾	
Tetric EvoCeram Bulk Fill	15 s
Tetric Basic White	20 s
Compomer ³⁾	
Compoglass F / Compoglass Flow	40 s
Indirect restorations / Luting materials	
Variolink II Base ⁴⁾ / Variolink Veneer / Variolink N Base ⁴⁾ / Variolink N Clear Veneer / Variolink Ultra ⁴⁾ / Dual Cement ⁵⁾ / Variolink II ⁶⁾	per mm ceramic: 10 s per segment
Miscellaneous	
Helioseal / Helioseal F / Helioseal Clear	20 s
Monopaque	30 s
MultiCore Flow / Multicore HB	40 s
Heliosit Orthodontic	20 s
Telio Add-On Flow	20 s
Telio Stains	15 s
IPS Empress Direct Color	15 s
IPS Empress Direct Opaque	30 s

LOW POWER Program

Reduced light intensity with reduced heat development for the polymerization of adhesives, liners, and restorative materials in areas near the pulp when restoring Class V cavities.

Adhesive	Curing time
AdheSE / AdheSE One F	
ExcTE / ExcTE D5C	10 s
Tetric N-Bond / Tetric N-Bond Self-Etch	
Heliobond (Syntac)	
Temporary materials	
Telio CS Inlay/Onlay	10 s
Systemp.inlay/Onlay	10 s
Fermit / Fermit N	10 s
Telio CS Link / Systemp.link	20 s per segment
Miscellaneous	
Heliosit Orthodontic	20 s
Vivaglass Liner	20 s

SOFT START Program

Step-by-step increase of the light intensity with reduced shrinkage stress and reduced heat development for the polymerization of restorative materials.

Filling materials	Curing time
Composite	
• 2 mm ¹⁾	
IPS Empress Direct / IPS Empress Direct Flow / Tetric EvoCeram / Tetric EvoFlow / Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	20 s
• 4 mm ²⁾	
Heliomolar / Heliomolar HB / Heliomolar Flow / Others	25 s
• 4 mm ²⁾	
Tetric EvoCeram Bulk Fill	20 s
Tetric Basic White	25 s
Compomer ³⁾	
Compoglass F / Compoglass Flow	40 s

- 1) Applies to a maximum layer thickness of 2 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case, e.g. with dentin shades)
- 2) Applies to a maximum layer thickness of 4 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case, e.g. with dentin shades)
- 3) Applies to a maximum layer thickness of 3 mm
- 4) Applies to light-curing (use of the base past only)
- 5) Applies to dual-curing

Start up

Cure Memory function

The last settings used, together with the combination of curing program and curing time, are automatically saved.

Start

The light is switched on by means of the start button. Once the selected curing time has elapsed, the curing program is automatically terminated. If desired, the light can be switched off before the set curing time has elapsed by pressing the start button again. The fan is activated simultaneously to the light. Once the curing time has elapsed, the fan continues to run for a certain time to cool the apparatus.

Acoustic signals

Acoustic signals can be heard for the following functions:

- Start (Stop)
- Every 10 seconds
- Program change
- Curing time change
- Error message

If desired, the volume of the acoustic signals can be adjusted.

For this purpose, press the blue volume button (see Product overview) while the light is off and reduce the volume or turn the acoustic signals off with the left time selection button. If you want to turn the acoustic signals back on or increase the volume, press the blue volume button and then the right time selection button.

Light intensity

The light intensity is maintained at a consistent level during operation. If the supplied 10 mm light probe is used, the light intensity has been calibrated to $800 \text{ mW/cm}^2 \pm 10\%$.

The use of light probe other than the one provided has a direct influence on the light intensity emitted.

In parallel-walled light probes (10 mm), the diameter is equal at both ends. In focussing light probes (10>8 mm light probe, Pin-Point light probe 6>2 mm), the diameter of the rear end is larger than that of the light emission window. The incident blue light is thus bundled to a smaller surface, which increases the light intensity emitted.

Pin-Point light probes are suitable for the polymerization of confined areas, such as the attachment of veneers prior to excess removal. For thorough curing, it is necessary to change the light probe.

Measuring the light intensity

The light intensity of the Bluephase C8 and the enclosed 10 mm light probe can be checked by means of the Bluephase Meter.

If the measured value does not correspond with the expected light intensity, proceed as follows:

- Check the selected curing program
- Clean light sensor if it is contaminated
- Remove the light probe and clean the light emission window of the handpiece with a cotton swab dipped in alcohol.
- Clean light probe if it is contaminated (see Maintenance and cleaning)
- Replace a damaged light probe with a new one

If these measures do not improve the results, please contact your dealer or your local Service Center.

Maintenance and cleaning

For reasons of hygiene, we recommend using a disposable protective sleeve for each patient. Make sure to fit the protective sleeve closely to the light probe. Disinfect contaminated surfaces of the device and anti-glare cones (e.g. FD 366/Dürr Dental, Incidin Liquid/Ecolab) and sterilize the light probe before each use if disposable protective sleeves are not used. Make sure that no liquids or other foreign substances enter the handpiece, the handpiece holder and particularly the power pack during cleaning (risk of electrical shock). Disconnect the handpiece from the power source when cleaning it.

Housing

Wipe the handpiece and handpiece holder with a customary aldehyde-free disinfecting solution. Do not clean with highly aggressive disinfecting solutions (e.g. solutions based on orange oil or with an ethanol content of more than 40%), solvents (e.g. acetone), or pointed instruments, which may damage or scratch the plastic. Clean soiled plastic parts with a soap solution.

Light probe

Before cleaning and/or disinfecting the light probe, it should be pretreated. This applies to both automated and manual cleaning and disinfection.

Pretreatment

- Remove gross contamination immediately after use or within 2 hours. For this purpose, thoroughly rinse the light probe under running water (for at least 10 seconds). Alternatively, use a suitable aldehyde-free disinfecting solution to prevent the fixation of blood.
- To remove contamination manually, use a soft brush or soft cloth. Partially polymerized composite can be removed with alcohol and a plastic spatula if necessary. Do not use sharp or pointed objects, as they may scratch the surface.

Cleaning and disinfection

For cleaning, immerse the light probe in a cleaning solution and ensure that it is sufficiently covered with liquid (ultrasound or careful brushing with a soft brush may support the effect). A neutral-enzymatic cleaning agent is recommended. When cleaning and disinfecting, please make sure that the agents used are free of

- organic, mineral and oxidizing acids (the minimum admissible pH value is 5.5)
- alkaline solutions (the maximum admissible pH value is 8.5)
- oxidizing agents (e.g. hydrogen peroxide)

Afterwards, remove the light probe from the solution and thoroughly rinse it under running fresh water (for at least 10 seconds). Cleaning in a thermal disinfectant is an effective alternative.

Sterilization

Thorough cleaning and disinfecting is imperative to ensure that the subsequent sterilization is effective. Use only autoclave sterilization for this purpose. The sterilization time (exposure time at sterilization temperature) is 4 minutes at 134 °C (273 °F); pressure should be 2 bar (29 psi). Dry the sterilized light probe using either the special drying program of your steam autoclave or hot air. The light probe has been tested for up to 200 sterilization cycles.


After that, check the light probe for damage. Hold it against light. If individual segments appear black, glass fibres are broken. If this is the case, replace the light probe with a new one.

Disposal

The curing light must not be disposed of as urban waste. Dispose unserviceable polymerization lights according to the corresponding legal requirements in your country.

What if ...

Symbol	Causes	Error rectification
	Apparatus is overheated or under-cooled	Allow the apparatus to cool down (or assume room temperature if it is undercooled) and try again after a certain time. If the error persists, please contact your dealer or your local Service Center.
	Electronic component of the power pack defective	Disconnect power pack from the power supply and then connect it again. If the error reappears, please contact your dealer or your local Service Center.
	Electronic component of the hand-piece defective	Please contact your dealer or your local Service Center.

Warranty / Procedure in case of repair

The warranty for Bluephase C8 is three years from the date of purchase. Malfunctions resulting from faulty material or manufacturing errors are repaired free of charge during the warranty period. The warranty does not provide the right to recover any material or non-material damage other than the ones mentioned. The apparatus must only be used for the intended purposes. Any other uses are contraindicated. The manufacturer does not accept any liability resulting from misuse and warranty claims cannot be accepted in such cases. This is particularly valid for:

- damage resulting from improper handling.
- damage of components resulting from wear under standard operating conditions.
- damage resulting from external influences, e.g. blows, drop to the floor.
- damage resulting from incorrect set-up or installation.
- damage resulting from connecting the unit to a power supply, the voltage and frequency of which do not comply with the ones stated on the rating plate.

- damage resulting from improper repairs or modifications that have not been carried out by certified Service Centres.

In case of a claim under warranty, the complete apparatus (handpiece, handpiece holder, power cord and light probe) must be returned, carriage paid, to the dealer or directly to Ivoclar Vivadent, together with the purchase document. Use the original packaging with the corresponding cardboard inserts for transportation. Repair work may only be carried out by a certified Ivoclar Vivadent Service Center.

In case of a defect that cannot be rectified, please contact your dealer or your local Service Center (see addresses on the back cover). A clear description of the defect or the conditions under which the defect occurred will facilitate locating the problem. Please enclose this description when returning the apparatus.

Product specifications

Delivery form

- 1 Handpiece holder,
- 1 Handpiece including power cord and power pack,
- 1 10 mm light probe,
- 3 Anti-glare cones,
- 1 Set of Instructions for Use

Accessories

The following accessories are available for Bluephase C8:

REF.	Description
608537	Light probe 10 mm, black (G2)
608538	Pin-point light probe 6/2 mm, black (G2)
551756	Anti-glare cones
592496	Anti-glare shield
627436	Handpiece holder C8 (G2)
607922	Bluephase Meter
608554	sleeves

Technical data

Operating voltage	5 VDC
Power Pack	100–240 VAC / 50-60 Hz / max. 0.4 A Output 5 VDC TYP 15.2630 Manufacturer: Friwo

Operating conditions

Temperature	+10 °C to +35 °C (+50 °F to 95 °F)
Relative humidity	30 % to 75 %
Ambient pressure	700 hPa to 1060 hPa
Dimensions of handpiece holder	L=205mm; B=150mm; H=85mm
Weight of the handpiece holder	190 g
Light source	Polywave® LED
Wavelength range	385–515 nm
Light intensity	800 mW/cm ² ± 10 %
Operating	5 min on / 6 min off (intermittently)
Dimensions of Handpiece:	L=260mm; B=42mm; H=120mm
Weight of the handpiece	200 g

Transportation and storage conditions:

Temperature	-20 °C to +60 °C (-4 to 140 °F)
Relative humidity	10 % to 75 %
Ambient pressure	500 hPa to 1060 hPa

Bluephase C8 has to be stored in closed, roofed rooms.
Protect the device from severe jarring.

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Appendix

Guidance and manufacturer's declaration - electromagnetic emission

The Bluephase C8 is intended for use in the electromagnetic environment specified below. The customer or the user of the Bluephase C8 should assure that it is used in such an environment.

Emision test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Bluephase C8 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Bluephase C8 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	P < 75W (pass without test)
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	P < 75W (pass without test)

Guidance and manufacturer's declaration - electromagnetic immunity

The Bluephase C8 is intended for use in the electromagnetic environment specified below. The customer or the user of the Bluephase C8 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of typical commercial or dental environment.
Surge IEC 61000-4-5	± 1 kV line - line ± 2 kV line - earth	± 1 kV line - line no prot. earth	Mains power quality should be that of typical commercial or dental environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U, (>95 % dip in U.) for 0.5 cycle 40 % U, (60 % dip in U.) for 5 cycles 70 % U, (30 % dip in U.) for 25 cycles <5 % U, (>95 % dip in U.) for 5 sec	<5 % U, (>95 % dip in U.) for 0.5 cycle 40 % U, (60 % dip in U.) for 5 cycles 70 % U, (30 % dip in U.) for 25 cycles <5 % U, (>95 % dip in U.) for 5 sec	Mains power quality should be that of typical commercial or dental environment. If the user of the Bluephase C8 requires continued operation during power mains interruptions, it is recommended that the Bluephase C8 be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or dental environment.

NOTE: U_i is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity

The Bluephase C8 is intended for use in the electromagnetic environment specified below. The customer or the user of the Bluephase C8 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should not be used closer to any part of the Bluephase C8, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			Recommended separation distance
Conducted RF IEC 61000-4-6	3 V _{ms} 150 kHz to 80 MHz	10 V	$d = 0.35 \sqrt[3]{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m	$d = 0.35 \sqrt[3]{P}$ 80 MHz to 800 MHz $d = 0.70 \sqrt[3]{P}$ 800 MHz to 2.5 GHz
			<p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic

environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Bluephase C8 is used exceeds the applicable RF compliance level above, the Bluephase C8 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Bluephase C8.

^b Over the frequency range 150 kHz to 80 MHz, field strength should be less than 10 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Bluephase C8

The Bluephase C8 is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled.

The customer or the user of the Bluephase C8 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Bluephase C8 as recommended below, according to the maximum output power of the communication equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 0.4 \sqrt[3]{P}$	80 MHz to 800 MHz $d = 0.4 \sqrt[3]{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt[3]{P}$
0.01	0.04	0.04	0.07
0.1	0.13	0.13	0.22
1	0.40	0.40	0.70
10	1.3	1.3	2.2
100	4.0	4.0	7.0

Ivoclar Vivadent – worldwide

Ivoclar Vivadent AG

Bendererstrasse 2
FL-9494 Schaan
Liechtenstein
Tel. +423 235 35 35
Fax +423 235 33 60
www.ivoclarvivadent.com

Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive
P.O. Box 367
Noble Park, Vic. 3174
Australia
Tel. +61 3 979 595 99
Fax +61 3 979 596 45
www.ivoclarvivadent.com.au

Ivoclar Vivadent Ltda.

Alameda Caiapós, 723
Centro Empresarial Tamboré
CEP 06460-110 Barueri – SP
Brazil
Tel. +55 11 2424 7400
Fax +55 11 3466 0840
www.ivoclarvivadent.com.br

Ivoclar Vivadent Inc.

1-6600 Dixie Road
Mississauga, Ontario
L5T 2Y2
Canada
Tel. +1 905 670 8499
Fax +1 905 670 3102
www.ivoclarvivadent.us

Ivoclar Vivadent (Shanghai)

Trading Co., Ltd.
2/F Building 1, 881 Wuding Road,
Jing An District
200040 Shanghai
China
Tel. +86 21 6032 1657
Fax +86 21 6176 0968
www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

Calle 134 No. 7-B-83, Of. 520
Bogotá
Colombia
Tel. +57 1 627 33 99
Fax +57 1 633 16 63
www.ivoclarvivadent.co

Ivoclar Vivadent SAS

B.P. 118
F-74410 Saint-Jorioz
France
Tel. +33 450 88 64 00
Fax +33 450 68 91 52
www.ivoclarvivadent.fr

Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2
D-73479 Ellwangen, Jagst
Germany
Tel. +49 (0) 79 61 / 8 89-0
Fax +49 (0) 79 61 / 63 26
www.ivoclarvivadent.de

Wieland Dental + Technik GmbH & Co. KG

Schwenninger Strasse 13
D-75179 Pforzheim
Germany
Tel. +49 (0) 72 31 / 37 05-0
Fax +49 (0) 72 31 / 35 79 59
www.wieland-dental.com

Ivoclar Vivadent Marketing (India) Pvt. Ltd.

503/504 Raheja Plaza
15 B Shah Industrial Estate
Veera Desai Road, Andheri (West)
Mumbai, 400 053
India
Tel. +91 (22) 2673 0302
Fax +91 (22) 2673 0301
www.ivoclarvivadent.in

Ivoclar Vivadent s.r.l.

Via Isonzo 67/69
40033 Casalecchio di Reno (BO)
Italy
Tel. +39 051 611 35 55
Fax +39 051 611 35 65
www.ivoclarvivadent.it

Ivoclar Vivadent K.K.

1-28-24-4F Hongo
Bunkyo-ku
Tokyo 113-0033
Japan
Tel. +81 3 6903 3535
Fax +81 3 5844 3657
www.ivoclarvivadent.jp

Ivoclar Vivadent Ltd.

12F W-Tower, 1303-37
Seocho-dong, Seocho-gu,
Seoul 137-855
Republic of Korea
Tel. +82 (2) 536 0714
Fax +82 (2) 596 0155
www.ivoclarvivadent.co.kr

Ivoclar Vivadent S.A. de C.V.

Av. Insurgentes Sur No. 863,
Piso 14, Col. Napoles
03810 México, D.F.
México
Tel. +52 (55) 50 62 10 00
Fax +52 (55) 50 62 10 29
www.ivoclarvivadent.com.mx

Ivoclar Vivadent Ltd.

12 Omega St, Rosedale
PO Box 303011 North Harbour
Auckland 0751
New Zealand
Tel. +64 9 914 99 99
Fax +64 9 914 99 90
www.ivoclarvivadent.co.nz

Ivoclar Vivadent Polska Sp. z o.o.

Al. Jana Pawla II 78
00-175 Warszawa
Poland
Tel. +48 22 635 54 96
Fax +48 22 635 54 69
www.ivoclarvivadent.pl

Ivoclar Vivadent Marketing Ltd.

Prospekt Andropova 18 korp. 6/
office 10-06
115432 Moscow
Russia
Tel. +7 499 418-03-00
Fax +7 499 418-03-10
www.ivoclarvivadent.ru

Ivoclar Vivadent Marketing Ltd.

Qlaya Main St.
Siricon Building No.14, 2nd Floor
Office No. 204
P.O. Box 300146
Riyadh 11372
Saudi Arabia
Tel. +966 1 293 83 45
Fax +966 1 293 83 44
www.ivoclarvivadent.com

Ivoclar Vivadent Pte. Ltd.

171 Chin Swee Road
#02-01 San Centre
Singapore 169877
Tel. +65 6535 6775
Fax +65 6535 4991
www.ivoclarvivadent.com

Ivoclar Vivadent S.L.U.

c/ Emilio Muñoz Nº 15
Entrada c/ Albaracin
E-28037 Madrid
Spain
Tel. + 34 91 375 78 20
Fax + 34 91 375 78 38
www.ivoclarvivadent.es

Ivoclar Vivadent AB

Dalvägen 14
S-169 56 Solna
Sweden
Tel. +46 (0) 8 514 93 930
Fax +46 (0) 8 514 93 940
www.ivoclarvivadent.se

Ivoclar Vivadent Liaison Office

: Tesvikiye Mahallesi
Sakayik Sokak
Nisantas' Plaza No:38/2
Kat:5 Daire:24
34021 Sisli – Istanbul
Turkey
Tel. +90 212 343 08 02
Fax +90 212 343 08 42
www.ivoclarvivadent.com

Ivoclar Vivadent Limited

Ground Floor Compass Building
Feldspar Close
Warrens Business Park
Enderby
Leicester LE19 4SE
United Kingdom
Tel. +44 116 284 78 80
Fax +44 116 284 78 81
www.ivoclarvivadent.co.uk

Ivoclar Vivadent, Inc.

175 Pineview Drive
Amherst, N.Y. 14228
USA
Tel. +1 800 533 6825
Fax +1 716 691 2285
www.ivoclarvivadent.us