Bluephase® Meter

Operating Instructions
Bedienungsanleitung
Mode d'emploi
Istruzioni d'uso
Instrucciones de uso
Instruções de Uso
Bruksanvisning
Brugsanvisning
Käyttöohjeet
Bruksanvisning
Productinformatie
Οδηγίες Χρήσεως
Kullanım Kılavuzu

CE

For dental use only!
 Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist.

Made in Austria
 Ivoclar Vivadent AG, FL-9494 Schaan / Liechtenstein

In USA distributed by
 Ivoclar Vivadent Inc.,
 175 Pineview Drive, Amherst,
 NY 14228, USA

 ivoclar
 vivadent®
 clinical

611504/WE3/2012-10-08/Rev4

English

Dear Customer,

Thank you for having purchased the Bluephase Meter from Ivoclar Vivadent. The light intensity of the polymerization light used is the decisive factor when it comes to the sufficient curing of restorations. In order to achieve adequate polymerization and thus durable high-quality composite restorations, it is recommended to check the light intensity of the polymerization lights regularly with a radiometer.

Intended use

The Bluephase Meter is a radiometer for measuring energy-rich blue light in the wavelength range of 385-515 nm emitted by LED polymerization units for dental use. The intended use also includes the observation of the notes and regulations in these Operating Instructions.

Indications

The Bluephase Meter is suitable for measuring the light intensity of LED polymerization units for use in dentistry. In order to achieve accurate measuring results, only light probes with a circular emission window must be used.

Signs and symbols

Contraindication

Symbols on the curing light

Double insulation (Apparatus complies with safety class II)

Observe Operating Instructions

The apparatus must not be disposed in the normal domestic waste. Please visit your local Ivoclar Vivadent website for detailed information about the disposal of the unit.

Safety notes

The Bluephase Meter is a medical electronic device which is subject to IEC 601-1 (EN 60601-1) and EMC standards. It complies with the relevant EU regulations (see appendix).

The apparatus has been shipped from the manufacturer in a safe and technically sound condition. In order to maintain this condition and to ensure risk-free operation, the notes and regulations in these Operating Instructions have to be observed.

NRTL (Nationally Recognized Testing Laboratory)

Warnings

- Protect the measuring device from heavy shock or falling down.
- Cleaning agents and other liquids must not enter the apparatus, since they may damage it.
- Do not expose the radiometer to sunlight.
- Both the storage and use of other electronic devices in the immediate vicinity of the apparatus may impair proper operation.
-  Portable and mobile high-frequency communication devices may affect medical apparatus. The use of mobile phones during operation is not allowed.
- Repair work may only be carried out by a certified Ivoclar Vivadent Service Center.
- Use only original spare parts.
- For dental use only! Suitable for measuring energy-rich blue light in the wavelength range of 385-515 nm and not for white light.

Operation

- 1) Place the light probe directly to the centering gauge and flush on the sensor.
- 2) Switch on the polymerization light.
Note: In order to achieve accurate measuring results, only light probes with a circular emission window must be used.
- 3) The Bluephase Meter automatically switches on at a light intensity of at least 300 mW/cm². A lower light intensity is indicated with "----" on the display.
- 4) Hold the polymerization light above the sensor for about 5 seconds.
- 5) The measuring result shows the power density in milliwatt per square centimeter.
- 6) Switch off the polymerization light.
- 7) The Bluephase Meter automatically switches off after 20 seconds of light exposure or 3 seconds after removing the light source.

What if...?

- ... **the Bluephase Meter shows a value below 400 mW/cm²?** 400 mW/cm² is the absolute minimum of light intensity for polymerization lights. If the light intensity is lower, composites can no longer be completely polymerized. Therefore, it is recommended to replace a curing light with a light intensity of less than 400 mW/cm².

- ... **the measuring results correspond with the data provided by the manufacturer?**

The polymerization light achieves the stipulated light intensity and can therefore be used according to the curing recommendations of the respective manufacturer or stipulations of the corresponding composite manufacturer.

- ... **the measuring results deviate strongly?**

If the measuring results strongly differ from the manufacturer's data, the following procedure is recommended:

- Check the selected curing program.
- Clean the sensor of the Bluephase Meter (see Maintenance and Cleaning).
- Clean the light probe with a disinfectant and a soft cloth. Carefully remove any possible material residue, such as composite material, from the light probe with e.g. fingernail or a plastic spatula. Do not use sharp or pointed instruments, since they may scratch the surface of the light probe and thus reduce light transmission.
- Check light probe for damage. If the light probe is held against a light source and certain segments appear black, then the glass fibres are broken. In this case, the light probe has to be replaced with a new one.

If these measures do not improve the results, the curing times have to be prolonged respectively in order to achieve sufficient polymerization. The following standard applies: With half the light intensity, the curing time has to be doubled (e.g. the measured light intensity is 500 mW/cm² while the data stipulated by the manufacturer is 1000 mW/cm² with 20 seconds of curing time; in this case, the curing time has to be doubled to 40 seconds).

- ... **"LO" is shown on the display?**

A low residual capacity is repeatedly shown by "LO" on the display. The batteries have to be replaced as soon as possible.

- ... **the display does not show anything?**

- After 20 seconds of exposure to light (e.g. also sunlight), the Bluephase Meter automatically changes into the stand-by mode. Remove the light source before you conduct a new measuring. If still no measuring is possible, cover the sensor with the hand for a short time.
- If these measures do not improve the results, check the batteries and replace them, if required.
- If there is still no display upon exposure to the light of a curing light despite these measures, send the Bluephase Meter to a certified Ivoclar Vivadent Service Center for repair.

Change of batteries

Batteries are consumables and must be exchanged regularly. For this purpose, remove the cover of the battery case at the rear of the unit and exchange the batteries. Subsequently, remount the cover of the battery case. Dispose of unserviceable batteries according to the national rules and regulations.

Maintenance and cleaning

Clean the radiometer with a soft cloth and customary aldehyde-free disinfecting solution. Do not clean with highly aggressive disinfecting solutions (e.g. solutions based on orange oil or with an ethanol content of more than 40%), solvents (e.g. acetone), or pointed instruments, which may damage or scratch the plastic.

- Do not spray disinfecting solution directly on the apparatus and do not hold it under running water in order to prevent liquids from entering the device.
- Do not sterilize the radiometer in the autoclave.

Warranty

The warranty for the Bluephase Meter is three years from the date of purchase. Malfunctions resulting from faulty material or manufacturing errors are repaired free of charge during the warranty period. The warranty does not provide the right to recover any material or non-material damage other than the ones mentioned. The apparatus must only be used for the intended purposes. Any other uses are contraindicated. The manufacturer does not accept any liability resulting from misuse and warranty claims cannot be accepted in such cases. This is particularly valid for:

- damage resulting from improper handling.
- damage to components resulting from wear under standard operating conditions (e.g. sensor).
- damage resulting from external influences, e.g. blows, drop to the floor;
- damage resulting from improper repairs or modifications that have not been carried out by certified Service Centers.

In case of a claim under warranty, the apparatus must be returned, carriage paid, to the dealer or directly to Ivoclar Vivadent in its original packaging together with the purchase document.

Procedure in case of repair

Repair work may only be carried out by a certified Ivoclar Vivadent Service Center. In case of a defect that cannot be rectified, please contact your dealer or your local Service Center (see addresses on the reverse side). A clear description of the defect or the conditions under which the defect occurred will facilitate locating the problem. Please enclose this description when returning the apparatus.

Product specifications

Delivery form

- 607922 Bluephase Meter
- 611504 Operating Instruction
- 3x Batteries (Type LR6/AA/1.5V)

Technical data

Operating voltage 4.5 VDC
 Battery operation 3x Type LR6/AA/1.5V

Measuring range:

Wavelength 385–515 nm
 Power density 300–2500 mW/cm²
 Only for light probes with a circular emission window Ø 7-13 mm
 Measuring accuracy +/- 20%

Operating conditions:

Temperature +10 °C to +40 °C
 Relative humidity 30% to 75%
 Ambient pressure 700 hPa to 1060 hPa

Dimensions L = 126 mm W = 80 mm H = 41 mm
 Weight 170 g

Transportation and storage conditions

- Temperature -20 °C to +70 °C
- Relative humidity 10% to 75%
- Ambient pressure 500 hPa to 1060 hPa
- bluephase radiometer has to be stored in closed, roofed rooms.
 Protect the device from severe jarring.

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Appendix

Guidance and manufacturer's declaration - electromagnetic emission

The "Bluephase Meter" is intended for use in the electromagnetic environment specified below. The customer or the user of the "Bluephase Meter" should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The "Bluephase Meter" uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The "Bluephase Meter" is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	No power network input	
Voltage fluctuations / flicker emissions IEC 61000-3-3	No power network input	

Guidance and manufacturer's declaration - electromagnetic immunity

The "Bluephase Meter" is intended for use in the electromagnetic environment specified below. The customer or the user of the "Bluephase Meter" should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	n.a. (no cables connected)	
Surge IEC 61000-4-5	± 1 kV line - line ± 2 kV line - earth	n.a. (no cables connected)	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	n.a. (no mains power)	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or dental environment.

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Appendix

Guidance and manufacturer's declaration - electromagnetic immunity

The "Bluephase Meter" is intended for use in the electromagnetic environment specified below. The customer or the user of the "Bluephase Meter" should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should not be used closer to any part of the "Bluephase Meter", including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			Recommended separation distance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	n.a.	No cables connected to the device.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted

theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the "Bluephase Meter" is used exceeds the applicable RF compliance level above, the "Bluephase Meter" should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the "Bluephase Meter".

^b Over the frequency range 150 kHz to 80 MHz, field strength should be less than 0 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the “Bluephase Meter”

The “Bluephase Meter” is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the “Bluephase Meter” can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the “Bluephase Meter” as recommended below, according to the maximum output power of the communication equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	n.a.	0.12	0.23
0.1	n.a.	0.38	0.73
1	n.a.	1.2	2.3
10	n.a.	3.8	7.3
100	n.a.	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 4: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Ivoclar Vivadent – worldwide

Ivoclar Vivadent AG

Bendererstrasse 2
FL-9494 Schaan
Liechtenstein
Tel. +423 235 35 35
Fax +423 235 33 60
www.ivoclarvivadent.com

Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive
P.O. Box 367
Noble Park, Vic. 3174
Australia
Tel. +61 3 979 595 99
Fax +61 3 979 596 45
www.ivoclarvivadent.com.au

Ivoclar Vivadent Ltda.

Alameda Caiapós, 723
Centro Empresarial Tamoré
CEP 06460-110 Barueri – SP
Brazil
Tel. +55 11 2424 7400
Fax +55 11 3466 0840
www.ivoclarvivadent.com.br

Ivoclar Vivadent Inc.

1-6600 Dixie Road
Mississauga, Ontario
L5T 2Y2
Canada
Tel. +1 905 670 8499
Fax +1 905 670 3102
www.ivoclarvivadent.us

Ivoclar Vivadent (Shanghai) Trading Co., Ltd.

2/F Building 1, 881 Wuding Road,
Jing An District
200040 Shanghai
China
Tel. +86 21 6032 1657
Fax +86 21 6176 0968
www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

Calle 134 No. 7-B-83, Of. 520
Bogotá
Colombia
Tel. +57 1 627 33 99
Fax +57 1 633 16 63
www.ivoclarvivadent.com

Ivoclar Vivadent SAS

B.P. 118
F-74410 Saint-Jorioz
France
Tel. +33 450 88 64 00
Fax +33 450 68 91 52
www.ivoclarvivadent.fr

Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2
D-73479 Ellwangen, Jagst
Germany
Tel. +49 (0) 79 61 / 8 89-0
Fax +49 (0) 79 61 / 63 26
www.ivoclarvivadent.de

Ivoclar Vivadent Marketing (India) Pvt. Ltd.

503/504 Raheja Plaza
15 B Shah Industrial Estate
Veera Desai Road, Andheri (West)
Mumbai, 400 053
India
Tel. +91 (22) 2673 0302
Fax +91 (22) 2673 0301
www.ivoclarvivadent.in

Ivoclar Vivadent s.r.l.

Via Isonzo 67/69
40033 Casalecchio di Reno (BO)
Italy
Tel. +39 051 611 35 55
Fax +39 051 611 35 65
www.ivoclarvivadent.it

Ivoclar Vivadent K.K.

1-28-24-4F Hongo
Bunkyo-ku
Tokyo 113-0033
Japan
Tel. +81 3 6903 3535
Fax +81 3 5844 3657
www.ivoclarvivadent.jp

Ivoclar Vivadent Ltd.

12F W-Tower, 1303-37
Seocho-dong, Seocho-gu,
Seoul 137-855
Republic of Korea
Tel. +82 (2) 536 0714
Fax +82 (2) 596 0155
www.ivoclarvivadent.co.kr

Ivoclar Vivadent S.A. de C.V.

Av. Insurgentes Sur No. 863,
Piso 14, Col. Napoles
03810 México, D.F.
México
Tel. +52 (55) 50 62 10 00
Fax +52 (55) 50 62 10 29
www.ivoclarvivadent.com.mx

Ivoclar Vivadent Ltd.

12 Omega St, Rosedale
PO Box 303011 North Harbour
Auckland 0751
New Zealand
Tel. +64 9 914 99 99
Fax +64 9 914 99 90
www.ivoclarvivadent.co.nz

Ivoclar Vivadent Polska Sp. z o.o.

Al. Jana Pawła II 78
00-175 Warszawa
Poland
Tel. +48 22 635 54 96
Fax +48 22 635 54 69
www.ivoclarvivadent.pl

Ivoclar Vivadent Marketing Ltd.

Prospekt Andropova 18 korp. 6/
office 10-06
115432 Moscow
Russia
Tel. +7 499 418-03-00
Fax +7 499 418-03-10
www.ivoclarvivadent.ru

Ivoclar Vivadent Marketing Ltd.

Qlaya Main St.
Siricon Building No.14, 2nd Floor
Office No. 204
P.O. Box 300146
Riyadh 11372
Saudi Arabia
Tel. +966 1 293 83 45
Fax +966 1 293 83 44
www.ivoclarvivadent.com

Ivoclar Vivadent Pte. Ltd.

171 Chin Swee Road
#02-01 San Centre
Singapore 169877
Tel. +65 6535 6775
Fax +65 6535 4991
www.ivoclarvivadent.com

Ivoclar Vivadent S.L.U.

c/ Emilio Muñoz Nº 15
Entrada c/ Albarracín
E-28037 Madrid
Spain
Tel. + 34 91 375 78 20
Fax + 34 91 375 78 38
www.ivoclarvivadent.se

Ivoclar Vivadent AB

Dalvägen 14
S-169 56 Solna
Sweden
Tel. +46 (0) 8 514 93 930
Fax +46 (0) 8 514 93 940
www.ivoclarvivadent.se

Ivoclar Vivadent Liaison Office

: Tesvikkiye Mahallesi
Sakayik Sokak
Nisantas' Plaza No:38/2
Kat:5 Daire:24
34021 Sisli – Istanbul
Turkey
Tel. +90 212 343 08 02
Fax +90 212 343 08 42
www.ivoclarvivadent.com

Ivoclar Vivadent Limited

Ground Floor Compass Building
Feldspar Close
Warrens Business Park
Enderby
Leicester LE19 4SE
United Kingdom
Tel. +44 116 284 78 80
Fax +44 116 284 78 81
www.ivoclarvivadent.co.uk

Ivoclar Vivadent, Inc.

175 Pineview Drive
Amherst, N.Y. 14228
USA
Tel. +1 800 533 6825
Fax +1 716 691 2285
www.ivoclarvivadent.us

ivoclar
vivadent®
clinical