**ivoclar vivadent**

ISO 13485  
  
 CE  

Recommended for inlays, onlays, 3/4 crowns, telescopes, crown crowns, posts, short and long span bridges, PFM crowns, implant superstructures, partial dentures.

**Indications**

Composition

40.3 Au  
 39.3 Pt  
 9.7 Ag  
 9.3 Cu  
 7.2 Zn  
 6.3 Ni  
 5.8 In  
 4.8 Sn  
 3.8 Pb  
 2.8 Bi  
 1.8 Fe  
 0.8 Co

**ivoclar vivadent**

Recommended for inlays, onlays, 3/4 crowns, telescopes, crown crowns, posts, short and long span bridges, PFM crowns, implant superstructures, partial dentures.

**Indications**

Composition

40.3 Au  
 39.3 Pt  
 9.7 Ag  
 9.3 Cu  
 7.2 Zn  
 6.3 Ni  
 5.8 In  
 4.8 Sn  
 3.8 Pb  
 2.8 Bi  
 1.8 Fe  
 0.8 Co

**ivoclar vivadent**

Recommended for inlays, onlays, 3/4 crowns, telescopes, crown crowns, posts, short and long span bridges, PFM crowns, implant superstructures, partial dentures.

**Indications**

Composition

40.3 Au  
 39.3 Pt  
 9.7 Ag  
 9.3 Cu  
 7.2 Zn  
 6.3 Ni  
 5.8 In  
 4.8 Sn  
 3.8 Pb  
 2.8 Bi  
 1.8 Fe  
 0.8 Co

**ivoclar vivadent**

Recommended for inlays, onlays, 3/4 crowns, telescopes, crown crowns, posts, short and long span bridges, PFM crowns, implant superstructures, partial dentures.

**Indications**

Composition

40.3 Au  
 39.3 Pt  
 9.7 Ag  
 9.3 Cu  
 7.2 Zn  
 6.3 Ni  
 5.8 In  
 4.8 Sn  
 3.8 Pb  
 2.8 Bi  
 1.8 Fe  
 0.8 Co

**ivoclar vivadent Worldwide**

**Australia**  
 Ivoclar Vivadent Pty. Ltd.  
 1 - 5 Overseas Drive  
 P.O. Box 367  
 Noble Park, Vic. 3174  
 Tel. +61 3 979 595 99  
 Fax +61 3 979 596 45  
 www.ivoclarvivadent.com.au

**Brazil**  
 Ivoclar Vivadent Ltda.  
 Alameda Caiapós, 723  
 Centro Empresarial Tamboré  
 CEP 06460-110  
 Barueri - SP  
 Tel. +55 11 2424 7400  
 Fax +55 11 2424 7440  
 www.ivoclarvivadent.com.br

**Canada**  
 Ivoclar Vivadent Inc.  
 1-6600 Dixie Road  
 Mississauga, Ontario L5T 2Y2  
 Tel. +1 905 670 8459  
 Fax +1 905 670 3102  
 www.ivoclarvivadent.us

**China**  
 Ivoclar Vivadent (Shanghai)  
 Trading Co., Ltd.  
 21F Building 1,  
 881 Wuding Road  
 Jing An District  
 Shanghai 200040  
 Tel. +86 21 6032 1657  
 Fax +86 21 6176 0968  
 www.ivoclarvivadent.cn

**Colombia**  
 Ivoclar Vivadent Marketing Ltd.  
 Calle 134 No. 7-B-83, Of. 520  
 Bogotá  
 Tel. +57 1 623 37 99  
 Fax +57 1 633 16 63  
 www.ivoclarvivadent.com

**France**  
 Ivoclar Vivadent SAS  
 B.P. 118  
 74410 Saint-Jorioz  
 Tel. +33 450 88 64 00  
 Fax +33 450 88 91 52  
 www.ivoclarvivadent.fr

**Germany**  
 Ivoclar Vivadent GmbH  
 Dr. Adolf-Schneider-Str. 2  
 73479 Ellwangen, Jagst  
 Tel. +49 79 61 8 89 9  
 Fax +49 79 61 63 26  
 www.ivoclarvivadent.de

**India**  
 Ivoclar Vivadent Marketing  
 (India) Pvt. Ltd.  
 881 Wuding Road  
 Jing An District  
 Shanghai 200040  
 Tel. +86 21 6032 1657  
 Fax +86 21 6176 0968  
 www.ivoclarvivadent.cn

**Japan**  
 Ivoclar Vivadent K.K.  
 B.P. 118  
 74410 Saint-Jorioz  
 Tel. +33 450 88 64 00  
 Fax +33 450 88 91 52  
 www.ivoclarvivadent.fr

**Mexico**  
 Ivoclar Vivadent S.A. de C.V.  
 Av. Insurgentes Sur No. 863  
 03810 Mexico, D.F.  
 Tel. +52 55 5062 1000  
 Fax +52 55 5062 1029  
 www.ivoclarvivadent.com.mx

**New Zealand**  
 Ivoclar Vivadent Ltd.  
 12 Omega Street, Rosedale  
 PO Box 303011 North Harbour  
 Auckland 0751  
 Tel. +64 9 914 9999  
 Fax +64 9 914 9990  
 www.ivoclarvivadent.co.nz

**Poland**  
 Ivoclar Vivadent S.r.l.  
 ul. Jana Pawla II 78  
 00175 Warszawa  
 Tel. +48 22 635 54 96  
 Fax +48 22 635 54 69  
 www.ivoclarvivadent.pl

**Russian Federation**  
 Ivoclar Vivadent Marketing Ltd.  
 Prospekt Andropova 18, Korus 6  
 Office 10-06  
 115432 Moscow  
 Tel. +7 499 418 0300  
 Fax +7 499 418 0310  
 www.ivoclarvivadent.ru

**Spain**  
 Ivoclar Vivadent S.L.U.  
 C/ta. Fuencarral, 24 Portal 1 Bajó  
 28108 Alcobendas (Madrid)  
 Tel. +34 913 75 78 20  
 Fax +34 913 75 78 38  
 www.ivoclarvivadent.es

**Sweden**  
 Ivoclar Vivadent AB  
 Djalnågen 14  
 169 56 Solna  
 Tel. +46 8 514 93 930  
 Fax +46 8 514 93 940  
 www.ivoclarvivadent.se

**Turkey**  
 Ivoclar Vivadent Liaison Office  
 Tesvikliye Mahallesi  
 Sakayik Sokak  
 Nisantasi Plaza No: 40  
 Kat: 6 Daire: 31-32  
 34021 Sisi Istanbul  
 Tel. +90 212 343 0802  
 Fax +90 212 343 0842  
 www.ivoclarvivadent.com.tr

**UK**  
 Ivoclar Vivadent Limited  
 Compass Building  
 Fiddigap Road  
 Warrens Business Park  
 Enderby  
 Leicestershire LE19 4SE  
 Tel. +44 116 284 78 80  
 Fax +44 116 284 78 81  
 www.ivoclarvivadent.co.uk

**EN INSTRUCTIONS FOR USE**

**MODELLATION**  
 Design the framework in a reduced anatomical shape taking the planned veneer into consideration. Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm. Make sure the framework demonstrates adequate stability of shape. Avoid sharp angles. Design the connector areas to be adequate for the position and alloy being used.

**SPRUNG**  
 Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center. The connection sprues between the reservoir and the coping should be 2.5-3.0 mm in length and width.

**INVESTMENT**  
 Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy). Use investment following the manufacturer's instruction.

**BURN-OUT**  
 The suggested burnout temperature: 750-820C/1380-1510F

**MELTING AND CASTING**  
 Use a separate ceramic crucible for each alloy and preheat the crucible in the burnout furnace. Used and new alloys must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturer's instructions for use. If you use Ivoclar Vivadent Magic Wand, set propane to 0.35 bar/5 psi and oxygen to 0.7 bar/10 psi pressure. Keep the alloy in the reducing atmosphere of the flame between the inner and outer cones. Do not use flux. After casting bench cool to room temperature.

**CASTING TEMPERATURE:** 1315-1375C/2400-2510F

**METAL PREPARATION**  
 Carefully divest and clean the object with Al<sub>2</sub>O<sub>3</sub>. Do not use a hammer for divesting the object to prevent deformation. Finish the framework with carbide burs or with ceramic-bonded grinding instruments. Avoid inhalation of dust during grinding!

**OXIDATION**  
 Blast surface with 50-100 micron Al<sub>2</sub>O<sub>3</sub> at max. 5.2 bar/75 psi pressure before oxidation. Subsequently, clean with ultrasonic or steam. Place the object on the firing tray and provide adequate support. Place the tray in a porcelain furnace set at a low temperature of 650C/1200F and increase the furnace temperature to 950C/1740F without vacuum with 5 min. hold time at top temperature.

**SOLDERS AND FLUXES**  
 Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/1112F. The soldering gap should be the same thickness as the soldering strip. Allow the object to cool slowly after soldering.

**Pre-Solder:** HGPKF 1015 Y, HGPKF 1030 Y  
**Post-Solder:** Low Fusing White Gold Solder, .585 Fine Solder  
**Laser weld material:** Laser Ceramic White

**POLISHING**  
 After glazing remove oxide and flux residue and finish/polish the framework with rubber finishers/polishers.

**INDICATIONS**  
 Recommended for inlays, onlays, 3/4 crowns, telescopes, crown crowns, posts, short and long span bridges, PFM crowns, implant superstructures, partial dentures.

**CONTRAINDICATIONS**  
 For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

**SIDE EFFECTS**  
 In individual cases, sensitivity or allergies to components of this alloy may occur.

**INTERACTIONS**  
 Galvanic effects may occur between different alloys in the same oral environment.  
 For additional information look into the alloy property chart.

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