Short and long span bridges, PFM crowns, posts.

Indications

Recommended for inlays, onlays, 3/4 crowns, telescope crowns, cone crowns, posts.

Table with 2 columns: Element, Content (Au, Pt, Pd, Ag, Cu, Sn, Zn, Ni, Pb, Cr, Ni, Ti, etc.)

Composition

Use a separate carbon/ceramic crucible for each alloy and preheat the crucible (ceramic) in the burnout furnace.

INSTRUCTIONS FOR USE

MODELLATION Design the framework in a reduced anatomical shape taking the planned veneer into consideration. Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm.

SPRUNG Provide the modelled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center. The connection sprues between the reservoir and the coping should be 2.5-3.0 mm in length and width.

INVESTMENT Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy). Use investment following the manufacturer's instruction.

BURN-OUT The suggested burnout temperature: 750-820C/1380-1510F MELTING AND CASTING Use a separate carbon/ceramic crucible for each alloy and preheat the crucible (ceramic) in the burnout furnace.

OXIDATION Blast surface with 50-100 micron Al2O3 at max. 4.5 bar/65 psi pressure before oxidation. Subsequently, clean with ultrasonic or steam.

HEAT TREATMENT Hardening: 540C/1000F for 15 minutes: air cool. SOLDERS AND FLUXES Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/1112F.

POLISHING After glazing remove oxide and flux residue and finish/polish the framework with rubber finishers/polishers. INDICATIONS Recommended for inlays, onlays, 3/4 crowns, crowns, telescope crowns, cone crowns, posts, short and long span bridges, PFM crowns.

CONTRAINDICATIONS For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended. SIDE EFFECTS In individual cases, sensitivity or allergies to components of this alloy may occur.

INTERACTIONS Galvanic effects may occur between different alloys in the same oral environment. For additional information look into the alloy property chart.

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Aquarius Hard Rich yellow, gold-based dental metal ceramic casting alloy, Type 4. Ivoclar Vivadent Worldwide. Australia, Brazil, Canada, China, Colombia, Germany, India, Italy, Mexico, New Zealand, Poland, Russia, Spain, Sweden, Switzerland, Taiwan, Thailand, Turkey, USA, Vietnam, UK.

Polishing and Oxidation instructions. After glazing remove oxide and flux residue and finish/polish the framework with rubber finishers/polishers. INDICATIONS Recommended for inlays, onlays, 3/4 crowns, crowns, telescope crowns, cone crowns, posts, short and long span bridges, PFM crowns.

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Ivoclar Vivadent logo and contact information. ISO 13485 Quality Management System Certified. Made in U.S.A. www.ivoclarvivadent.com

SPRUNG METHOD RECOMMENDATIONS DIRECT: single crowns, inlays and onlays. INDIRECT: multiple units and multiple single crowns.

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