Harmony® 2

Ag-based dental casting alloy, Type 3

EN Instructions for use
Ag-based dental casting alloy, Type 3

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Ag-based dental casting alloy, Type 3

■ INDICATIONS *

Telescope Crowns, Conus Crowns, Bridges

■ WAXING / MODELLATION

Design the framework in a volume total or in a reduced anatomical shape taking the planned veneer into consideration. For Lab Composite veneering material use mechanical retainers. Small crowns require a minimum thickness of 0.3 mm. Abutment crowns require a minimum thickness of 0.5 mm. Ensure the framework design provides adequate support for the veneering material. Avoid sharp angles. Connectors must have the required dimensions to provide resistance to deformation. Create large surface areas for planned soldering, with a gap of 0.05–0.2 mm.

■ SPRUING

Provide the modelled single-toothrestoration or bridge framework with sprues of a suitable size. Reservoir, spacer leads, connectors and other parts shall be placed according to the specific technique used. Waxing must be done directly on the framework. Be sure to position the connector sprues between the reservoir and the casting so that a maximum of 2.5–3.0 mm in length and width. The wax pattern including the sprues must be weighed in grams in order to determine the needed amount of alloy. Wax conversion formula: wax weight (gram) x alloy density = grams of alloy required.

■ INVESTING

Use a phosphate-/gypsum-bonded investment material. Follow the manufacturer's instructions.

■ PREHEATING / BURN-OUT

Recommended burn-out temperature: 650–705 °C / 1200–1300 °F

■ MELTING AND CASTING

Torch: Propane 0.15 bar / 2 psi; Oxygen 0.35 bar / 5 psi
Other specifics may be required by the type of casting machine. It is recommended to use a separate and clean carbon/ceramic crucible for each alloy. The recommended ratio of used material to new material is 1:1. Use casting flux if necessary.

Casting Temperature: 1170–1230 °C / 2140–2245 °F

■ SPARK PLUG FINISHING

After bench cooling, carefully dress and clean the casting with aluminum oxide (Al₂O₃). Do not use a hammer for dressing. Finish the casting with carbide burs and/or with ceramic-bonded grinding instruments. For veneering with Lab Composite material, the framework must be finished and polished. Subsequently, steam clean or ultrasonic clean with distilled water or ethanol and dry the framework. Use the appropriate Lab Composite material, following the manufacturer's instructions.

■ HEAT TREATMENT

Soften: 705 °C / 1300 °F for 15 min; bench cool

Hardening: 350 °C / 660 °F for 30 min; bench cool

■ SOLDERING AND LASER WELDING

The soldering gap should not be wider than the thickness of the soldering material. Allow the soldered casting to cool before use. Use flux sparingly.

Post Solder: .585 Fine Gold Solder, Porta OP Solder 710 Flux: Bondal Flux

■ POLISHING

Carefully remove any oxide and flux residue. Smooth the metal surfaces with rubber polishers. Polish to a high gloss finish using polishing paste. Subsequently, clean using ultrasonic cleaning equipment or careful steam cleaning.

ADDITIONAL SAFETY CONCERN AND INSTRUCTIONS

■ CONTRAINDICATIONS

For patients with known allergy/sensitivity to any major or minor elements of this alloy, consultation with a physician is recommended. Alloy is not to be used for any application not included within the indications.

■ SIDE EFFECTS

In individual cases, sensitivity or allergies to elements of this alloy may occur. Ivoclar Vivadent makes no claims regarding the MRI-compatibility of dental alloys. It is recommended that the patient be made aware of the possibility for dental alloys to affect MRI results and to disclose the presence of dental alloys to the MRI technician prior to conducting a test.

■ INTERACTIONS

Galvanic effects may occur between different or dissimilar alloys in the same oral environment.

■ CAUTION

Metal vapors and metal dust are harmful if inhaled. Therefore, the use of extraction equipment and/or suitable protective masks is advised!

■ STORAGE CONDITION

Store in a dry environment at room temperature.

■ DISCLAIMER

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the instructions for use. Use only where indicated for dentures resulting from failure to observe instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and for any purpose not explicitly stated in the instructions. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

■ PROCESSING DATA

Investment Material: Phosphate-bonded / Gypsum-bonded

Preheating / Burn-out Temperature: 650–705 °C / 1200–1300 °F

Crucible: Carbon crucible / Ceramic crucible

Casting Temperature: 1170–1230 °C / 2140–2245 °F

Recommended Lab Composite: SR Nexo®

Softening: 705 °C / 1300 °F for 15 min; bench cool

Hardening: 350 °C / 660 °F for 30 min; bench cool

Post Solder / Flux: .585, Porta OP Solder 710 Bondal Flux

Laser Welding Wire: Laser C&B Yellow

■ TECHNICAL DATA (ISO 22674:2016)

Type / Color: 3 Yellow

Density (g/cm³): 9.8

Melting Range (Solidus/Liquidus): 875–1120 °C / 1605–2050 °F

Elastic Modulus (GPa/psi): 80 As Cast 11,603,000

Vickers Hardness: 175 160 180

Tensile Strength (MPa/psi): 540/78,320 490/71,070 520/75,420

0.2% Proof Stress (MPa/psi): 300/43,510 270/39,160 270/39,160

Elongation (%): 5 5 5

* See TYPE CLASSIFICATION DUE TO PHYSICAL PROPERTIES

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