
# Cervitec<sup>®</sup> Plus

## Dossier Scientifique

# Sommaire

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
*Dubielecka M, Slotwinska SM (2005): Suppression of caries in mothers and caries risk in offspring [15]*

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*Duskova J, Broukal Z, Kratky M (2000): Inhibition of the oral Streptococcus mutans transfer in the mother and child care: results in mothers and infants in the 3rd year of study [16]*

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*Plotzitza B, Kneist S, Berger J, Hetzer G (2005): Efficacy of chlorhexidine varnish applications in prevention of early childhood caries [17]*

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*Baca B, Munoz M, M Bravo, P Junco, AP Baca (2004): Effectiveness of chlorhexidin-thymol varnish in preventing caries lesions in primary molars [18]*

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
*Sköld-Larsson K, Fornell AC, Lussi A, Twetman S (2004): Effect of topical applications of a chlorhexidine / thymol-containing varnish on fissure caries assessed by laser fluorescence [21]*

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
*Baca P, Muñoz MJ, Junco P, Baca AP (2002): Effectiveness of chlorhexidine-thymol varnish for caries reduction in permanent first molars of 6-7-year-old children: 24-month clinical trial [7]*

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*Araujo A, Naspitz G, Chelotti A, Cai S (2002): Effect of Cervitec on mutans streptococci in plaque and on caries formation on occlusal fissures of erupting permanent molars [22]*

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*Joharji RM, Adenubi JO (2001): Prevention of pit and fissure caries using antimicrobial varnish: 9 month clinical evaluation [6]*

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*Heintze SD, Twetman S: Interdental mutans streptococci suppression in vivo (2002): a comparison of different chlorhexidine regimens in relation to restorative material [31]*

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*Twetman S, Petersson L (1999): Interdental caries incidence and progression in relation to mutans streptococci suppression after chlorhexidine / thymol varnish treatments in schoolchildren [25]*

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*Twetman S, Petersson L (1997): Effect of different chlorhexidine varnish regimes on mutans streptococci levels in interdental plaque and saliva [29]*

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*Øgaard B, Larsson E, Henriksson T, Birkhed D, Bishara S (2001): Effects of combined application of antimicrobial and fluoride varnishes in orthodontic patients [32]*

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*Eronat C, Alpöz AR (1997): Effect of Cervitec varnish on the salivary Streptococcus mutans levels on the patients with fixed orthodontic applications [34]*

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
*M. Madléna, S. Marton, G. Nagy (2000): Effect of chlorhexidine varnish on bacterial levels in plaque and saliva during orthodontic treatment [5]*

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
*Brailsford SR, Fiske J, Gilbert S, Clark D, Beighton D (2002): The effects of the combination of chlorhexidine/thymol- and fluoride-containing varnishes on the severity of root caries lesions in frail institutionalised elderly people [42]*

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*Ekenbäck BS, Linder LE, Lönnies H (2000): Effect of four dental varnishes on the colonization of cariogenic bacteria on exposed sound root surfaces [43]*

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
Wallman C, Birkhed D (2002): *Effect of chlorhexidine varnish and gel on mutans streptococci in margins of restorations in adults [48]*

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#### 4.2.8.3 Implants

*Johansson LA, Ekfeldt E, Petersson LG, Edwardsson S (1994): Antimicrobial effect of a chlorhexidine varnish, Cervitec in healing caps in osseointegrated prosthetic treatment [51]*

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*Besimo CE, Guindy JS, Lewetag D, Meyer J (1999): Prevention of bacterial leakage into and from prefabricated screw-retained crowns on implants in vitro [52; 53]*

*Duarte AR, Rossetti PH, Rossetti LM, Torres SA, Bonachela WC (2006): In vitro sealing ability of two materials at five different implant-abutment surfaces [54]*

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## 6. Bibliographie

1. Addy M. Dentine hypersensitivity: new perspectives on an old problem. *Int Dent J* 2002;367- 375.
2. Rudhart A, Rompola E, Hopfenmüller W, Bernimoulin JP. Effectiveness of Cervitec and Fluor Protector in patients with dentin hypersensitivity. *J Dent Res* 1998;77:746.
3. Panduric V, Knezevic A, Tarle Z, Sutalo J. The efficiency of dentine adhesives in treating noncaries cervical lesions. *J Oral Rehabil* 2001;28:1168-1174.
4. Zhang Q, van Palenstein Helderma WH, van't Hof MA, Truin GJ. Chlorhexidine varnish for preventing dental caries in children, adolescents and young adults: a systematic review. *Eur J Oral Sci* 2006;114:449-455.
5. Madléna M, Vitalyos G, Marton S, Nagy G. Effect of chlorhexidine varnish on bacteril levels in plaque and saliva during orthodontic treatment. *J Clin Dent* 2000;11:42-46.
6. Joharji RM, Adenubi JO. Prevention of pit and fissure caries using an antimicrobial varnish: 9 month clinical evaluation. *J Dent* 2001;29:247-254.
7. Baca P, Munoz MJ, Bravo M, Junco P, Baca AP. Effectiveness of chlorhexidine-thymol varnish for caries reduction in permanent first molars of 6-7-year-old children: 24-month clinical trial. *Community Dent Oral Epidemiol* 2002;30:363-368.
8. Haukali G, Poulsen S. Effect of a varnish containing chlorhexidine and thymol (cervitec) on approximal caries in 13- to 16-year-old schoolchildren in a low caries area. *Caries Res* 2002;37:185-189.
9. Fennis-le YL, Verdonschot EH, Burgersdijk RC, Konig KG, van 't Hof MA. Effect of 6-monthly applications of chlorhexidine varnish on incidence of occlusal caries in permanent molars: a 3- year study. *J Dent* 1998;26:233-238.
10. de Soet JJ, Gruythuysen RJ, Bosch JA, van Amerongen WE. The effect of 6-monthly application of 40% chlorhexidine varnish on the microflora and dental caries incidence in a population of children in Surinam. *Caries Res* 2002;36:449-455.
11. Zhang Q, van 't Hof MA, Truin GJ, Bronkhorst EM, van Palenstein Helderma WH. Cariesinhibiting effect of chlorhexidine varnish in pits and fissures. *J Dent Res* 2006;85:469-472.
12. Forgie AH, Paterson M, Pine CM, Pitts NB, Nugent ZJ. A randomised controlled trial of the caries-preventive efficacy of a chlorhexidine-containing varnish in high-caries-risk adolescents. *Caries Res* 2000;34:432-439.
13. Bratthall D, Serinirach R, Rapisuwon S, Kuratana M, Luangjarmekorn V, Luksila K, Chaipanich P. A study into the prevention of fissure caries using an antimicrobial varnish. *Int Dent J* 1995;45:245-254.
14. Baca P, Junco P, Bravo M, Baca AP, Munoz MJ. Caries incidence in permanent first molars after discontinuation of a school-based chlorhexidine-thymol varnish program. *Community Dent Oral Epidemiol* 2003;31:179-183.
15. Dubielecka M, Slotwinska SM. Suppression of caries in mothers and caries risk in offspring. *J Dent Res* 2005;84 (Spec Iss B):Abstract # 0068.
16. Duskova J, Broukal Z, Kratky M. Inhibition of the oral steptococcu mutans - transfer in the mother and child care: results in mothers and infants in the 3rd year of study. *Caries Res* 2000;34:347.
17. Plotzitza B, Kneist S, Berger J, Hetzer G. Efficacy of chlorhexidine varnish applications in the prevention of early childhood caries. *Eur J Paediatr Dent* 2005;6:149-154.
18. Baca P, Munoz MJ, Bravo M, Junco P, Baca AP. Effectiveness of chlorhexidine-thymol varnish in preventing caries lesions in primary molars. *J Dent Child (Chic)* 2004;71:61-65.
19. Sanderink R, Bernhard H, Knoke M, Meyer J, Weber C, Weiger R. Curriculum – Orale Mikrobiologie und Immunologie. Quintessenz; Berlin 2004.

20. Gokalp S, Baseren M. Use of laser fluorescence in monitoring the durability and cariostatic effects of fluoride and chlorhexidine varnishes on occlusal caries: a clinical study. *Quintessence Int* 2005;36:183-189.
21. Sköld-Larsson K, Fornell AC, Lussi A, Twetman S. Effect of topical applications of a chlorhexidine/thymol-containing varnish on fissure caries assessed by laser fluorescence. *Acta Odontol Scand* 2004;62:339-342.
22. Araujo AM, Naspitz GM, Chelotti A, Cai S. Effect of Cervitec on mutans streptococci in plaque and on caries formation on occlusal fissures of erupting permanent molars. *Caries Res* 2002;36:373-376.
23. Attin R, Tuna A, Attin T, Brunner E, Noack MJ. Efficacy of differently concentrated chlorhexidine varnishes in decreasing Mutans streptococci and lactobacilli counts. *Arch Oral Biol* 2003;48:503-509.
24. Petersson LG, Magnusson K, Andersson H, Almquist B, Twetman S. Effect of quarterly treatments with a chlorhexidine and a fluoride varnish on approximal caries in caries-susceptible teenagers: A 3-year clinical study. *Caries Res* 2000;34:140-143.
25. Twetman S, Petersson LG. Interdental caries incidence and progression in relation to mutans streptococci suppression after chlorhexidine-thymol varnish treatments in schoolchildren. *Acta Odontol Scand* 1999;57:144-148.
26. Petersson LG, Magnusson K, Andersson H. Effect of semi-annual applications of a chlorhexidine/fluoride varnish mixture on approximal caries incidence in schoolchildren. A three-year radiographic study. *Eur J Oral Sci* 1998;106:623-627.
27. Twetman S, Petersson LG. Comparison of the efficacy of three different chlorhexidine preparations in decreasing the levels of mutans streptococci in saliva and interdental plaque. *Caries Res* 1998;32:113-118.
28. Twetman S, Petersson LG. Efficacy of a chlorhexidine and a chlorhexidine-fluoride varnish mixture to decrease interdental levels of mutans streptococci. *Caries Res* 1997;31:361-365.
29. Twetman S, Petersson LG. Effect of different chlorhexidine varnish regimens on mutans streptococci levels in interdental plaque and saliva. *Caries Res* 1997;31:189-193.
30. Cosyn J, Wyn I, De Rouck T, Collys K, Bottenberg P, Matthijs S, Sabzevar MM. Short-term anti-plaque effect of two chlorhexidine varnishes. *J Clin Periodontol* 2005;32:899-904.
31. Heintze SD. Interdental mutans streptococci suppression in vivo: a comparison of different chlorhexidine regimens in relation to restorative material. *Am J Dent* 2002;15:103-108.
32. Øgaard B, Larsson E, Henriksson T, Birkhed D, Bishara SE. Effects of combined application of antimicrobial and fluoride varnishes in orthodontic patients. *Am J Orthod Dentofacial Orthop* 2001;120:28-35.
33. Yucel-Lindberg T, Twetman S, Sköld K, Modeer T. Effect on an antibacterial dental varnish on the levels of prostanooids, leukotriene B4, and interleukin-1 in gingival crevicular fluid. *Acta Odontol Scand* 1999;57:23-27.
34. Eronat C, Alpoz AR. Effect of Cervitec varnish on the salivary *Streptococcus mutans* levels in the patients with fixed orthodontic appliances. *J Marmara Univ Dent Fac* 1997;2:605-608.
35. Sköld K, Twetman S, Hallgren A, Yucel-Lindberg T, Modeer T. Effect of a chlorhexidine/thymol-containing varnish on prostaglandin E2 levels in gingival crevicular fluid. *Eur J Oral Sci* 1998;106:571-575.
36. Øgaard B, Larsson E, Glans R, Henriksson T, Birkhed D. Antimicrobial effect of a chlorhexidine-thymol varnish (Cervitec) in orthodontic patients. A prospective, randomized clinical trial. *J Orofac Orthop* 1997;58:206-213.
37. Øgaard B, Larsson E, Henriksson T, Birkhed D, Bishara SE. Effects of combined application of an antimicrobial varnish Cervitec and a fluoride varnish Fluor Protector on white spot lesions, gingivitis, plaque and mutans streptococci in orthodontic patients. A prospective, randomized clinical trial. *Swed Dent J* 2000;24:0.

38. Twetman S, Hallgren A, Petersson LG. Effect of an antibacterial varnish on Mutans Streptococci in plaque from enamel adjacent to orthodontic appliances. *Caries Res* 1995;29:188-191.
39. Lynch E, Beighton D. Short-Term Effects of Cervitec on the Microflora of Primary Root Caries Lesions Requiring Restoration. *Caries Res* 1993;27:236.
40. Wicht MJ, Haak R, Lummert D, Noack MJ. Treatment of root caries lesions with chlorhexidine-containing varnishes and dentin sealants. *Am J Dent* 2003;16:25A-30A.
41. Johnson G, Almqvist H. Non-invasive management of superficial root caries lesions in disabled and infirm patients. *Gerodontology* 2003;20:9-14.
42. Brailsford SR, Fiske J, Gilbert S, Clark D, Beighton D. The effects of the combination of chlorhexidine/thymol- and fluoride-containing varnishes on the severity of root caries lesions in frail institutionalised elderly people. *J Dent* 2002;30:319-324.
43. Ekenbäck SB, Linder LE, Lönnies H. Effect of four dental varnishes on the colonization of cariogenic bacteria on exposed sound root surfaces. *Caries Res* 2000;34:70-74.
44. Schmeiser R, Schiffner U, Gülzow HJ. Auswirkungen eines chlorhexidinhaltigen Lackes auf die Demineralisation unter experimenteller Plaque. *Dental Spiegel* 1994;49:834-836.
45. Marren PV, Lynch E, Heath MR. Cervitec's effect on plaque mutans streptococci covering healthy root surfaces. *J Dent Res* 1997;76:93.
46. Huizinga ED, Ruben J, Arends J. Effect of an antimicrobial-containing varnish on root demineralisation in situ. *Caries Res* 1990;24:130-132.
47. Benz C, Benz B, El-Mahdy KR, Hickel R. The effect of preventive regimens on dentate elderly after 3 years. *J Dent Res* 2002;81(Spec Iss A):Abstract #2734.
48. Wallman C, Birkhed D. Effect of chlorhexidine varnish and gel on mutans streptococci in margins of restorations in adults. *Caries Res* 2002;36:360-365.
49. Karaman AI, Uysal T. Effectiveness of a hydrophilic primer when different antimicrobial agents are mixed. *Angle Orthod* 2004;74:414-419.
50. Polat O, Uysal T, Karaman AI. Effects of a chlorhexidine varnish on shear bond strength in indirect bonding. *Angle Orthod* 2005;75:1036-1040.
51. Johansson LA, Ekfeldt A, Petersson LG, Edwardsson S. Antimicrobial effect of a chlorhexidine varnish, Cervitec in healing caps in osseointegrated prosthetic treatment. *Swed Dent J* 1994;16:255-256.
52. Besimo CHE, Guindy JS, Lewetag D, Meyer J. Prevention of bacterial leakage into and from prefabricated screw-retained crowns on implants in vitro. *Int J of Oral & Maxillofacial Implants* 1999;14:654-660.
53. Besimo CHE, Guindy JS, Lewetag D, Besimo RH, Meyer J. Marginale Passgenauigkeit und Bakteriendichtigkeit von verschraubten implantatgetragenen Suprastrukturen. *Parodontologie* 2000;3:217-229.
54. Duarte AR, Rossetti PH, Rossetti LM, Torres SA, Bonachela WC. In vitro sealing ability of two materials at five different implant-abutment surfaces. *J Periodontol* 2006;77:1828-1832.
55. Busch J. Final Report on the Safety Assessment of Vinyl Acetate/ Crotonic Acid Copolymer. *Journal of the American College of Toxicology* 1983;2:125-140.
56. Zondlo Fiume M. Final report on the safety assessment of Acrylates Copolymer and 33 related cosmetic ingredients. *Int J Toxicol* 2002;21 Suppl 3:1-50.
57. EMEA. Thymol - Summary Report Report No. EMEA/MRL/075/96-Final. 1996.
58. EMEA. Chlorhexidine. EMEA Report No. EMEA/MRL/107/96-FINAL. 1996.
59. Berufsgenossenschaft\_der\_chemischen\_Industrie. Thymol - Toxikologische Bewertung No 259. 2000.
60. Ullmann L. Primary skin irritation study. RCC Report No. 295007. 1991.

61. Meurer K. Cytotoxicity assay in vitro: Evaluation of test items in the XTT-Test. RCC-CCR Report No. 1023501. 2006.
62. Sokolowski A. Salmonella typhimurium and Escherichia coli reverse mutation assay. RCCCCR Report No. 1023502. 2006.
63. EPA. Thymol - Biopesticide Registration Action Document.
64. EPA. Chlorhexidine diacetate - Reregistration Eligibility Decision (RED). United States Environmental Protection Agency. 1996.
65. Ullmann L. Contact Hypersensitivity. RCC Report No. 294996. 1991.

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Ivoclar Vivadent AG