Description

Variolink® Esthetic is a colour-stable, adhesive luting system for the permanent cementation of ceramic and composite resin restorations. Variolink Esthetic is offered in a purely light-curing (VarioLink Esthetic LC) and a dual-curing version (VarioLink Esthetic DC).

The special filler composition gives Variolink Esthetic a very high radiopacity.

Shades

Variolink Esthetic is available in the following five shade gradations:

Shade	Effect
Light +	Translucency approx. 5% (opaque)
Light	Translucency approx. 10%
Neutral	Translucency approx. 17%
Warm	Translucency approx. 12%
Warm +	Translucency approx. 8.5% (opaque)

The gradation of Variolink Esthetic shades is based on the effect that a luting composite exerts on the brightness value of the final restoration. Variolink Esthetic Neutral does not affect the brightness value. At the same time, it demonstrates the highest translucency and is thus neutral in shade. Light and Light + make the restorations lighter, while Warm and Warm + create a darker overall shade impression. In order to evaluate the overall effect of the restoration in conjunction with the various Variolink Esthetic shades prior to permanent cementation, we recommend using the Variolink Esthetic try-in pastes.

Composition

The monomer matrix of Variolink Esthetic is composed of urethane dimethacrylate and further methacrylate monomers. The inorganic fillers are ytterbium trifluoride and spheroid mixed oxide. Initiators, stabilizers and pigments are additional ingredients.

The particle size is 0.04-0.2 µm. The mean particle size is 0.1 µm.

The total volume of inorganic fillers is approx. 38%.

Indication

- Permanent adhesive luting of glass-ceramic, lithium disilicate glass-ceramic and composite restorations (inlays, onlays and veneers)
- Only use Variolink Esthetic LC for restorations with a low thickness of <2mm that have sufficient translucency (e.g. restorations made of IPS e.max® HT).
- The following overview provides recommendations for the selection of the suitable Variolink Esthetic LC shade:

Indications			
Shade	Effect	Veneers	Inlays / Onlays
Light +	considerably lightens	✓	--*
Light	slightly lightens	✓	✓
Neutral	no shade effect	✓	✓
Warm	slightly darkens/ more yellowish	✓	✓
Warm +	considerably darkens/ more yellowish	✓	--*

* The use of the shades Light+ and Warm+ may lead to visible restoration margins when cementing inlays/onlays.

Contraindication

The application of Variolink Esthetic is contraindicated

- if a dry working field cannot be established or the stipulated working procedures cannot be applied;
- if a patient is known to be allergic to any of the ingredients of Variolink Esthetic.

Side effects

Systemic side effects are not known to date. In individual cases, allergic reactions to individual components have been reported.

Variolink® Esthetic LC

EN Instructions for Use	FI Käyttöohjeet	NO Bruksanvisning	NL Gebruiksaanwijzing	EL Οδηγίες Χρήσεως	TR Kullanma Talimatı
- Light-curing resin-based dental luting material	- Valkokovetteinen yhdistelmämäruuvimestintti	- Lysherdende, permanent sementeringskomposit	- Tandheelkundig, lichtuithardend bevestigingsmateriaal op basis van kunsthars	- Φωτοολυμέριζουμενη αδονυπατρική ρητινώδης κονία	- Isılaç seritleşen rezin esaslı dental yapıştırma materyali
DE Gebrauchsinformation	FR Mode d'emploi	IT Istruzioni d'uso	ES Instrucciones de uso	PT Instruções de Uso	
- Lichthärtendes zahnärztliches Befestigungscomposite	- Composite de collage photopolymérisable à usage dentaire	- Cemento compósito dental fotoinducente	- Material dental de cementado fotopolimerizable basado en resina	- Materiais de cimentação dental à base de resina fotoativaada	

Interactions

Phenolic substances (e.g. eugenol, wintergreen oil) inhibit polymerization. Consequently, the application of products containing these components, e.g. mouth rinses and temporary cements, must be avoided.

Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process.

Therefore, do not disinfect the preparation and the syringe using oxidative agents. Disinfection may be carried out, for instance, by wiping with medical alcohol.

Alkaline jet mediums applied on dentin (e.g. Airflow) may compromise the effect of self-etching adhesives.

Hemostatic substances may inhibit polymerization and/or lead to discolouration. Therefore, the directions for use of these substances must be observed.

Application

For more detailed information also refer to the Instructions for Use of the products used in conjunction with Variolink Esthetic.

1 Removal of temporary restoration and cleaning of the cavity

Remove possible residue of the temporary luting composite from the cavity or the preparation with a polishing brush and an oil- and fluoride-free cleaning paste (e.g. Proxyt® fluoride-free prophy paste). Rinse with water spray.

Subsequently, dry with oil- and moisture-free air. Avoid overdrying.

Note:

Cleaning with alcohol can lead to dehydration of dentin.

2 Try-in of the restoration and isolation

Seat the restoration using the desired Variolink Esthetic try-in paste and check the shade, fit and occlusion of the restoration.

Care should be taken when checking the occlusion of fragile and brittle ceramic objects before they are permanently cemented, as there is a risk of fracture.

If necessary, make adjustments with fine diamonds at medium speed and with slight pressure and adequate water cooling. Polish ground surfaces.

Adequate relative or absolute isolation using suitable auxiliaries, such as OptraGate® or OptraDam® Plus, is required.

If the dental hard tissue is contaminated with saliva or blood during try-in, it needs to be cleaned again as described in section 1.

3 Pre-treatment of the restoration

3.1 If a restoration which has been pre-treated in the dental laboratory is contaminated with saliva or blood during try-in, the contaminated bonding surface of the restoration needs to be cleaned as follows:

- Following try-in, thoroughly rinse the restoration with water spray and dry with oil-free air.
- Shake Ivoclean before use and completely coat the bonding surface of the restoration with the product using a microbrush or brush.
- Leave Ivoclean to react for 20 s, then thoroughly rinse with water spray and dry with oil-free air.

3.2 Restorations that have not been pre-treated by the dental laboratory should be rinsed with water spray and dried after the try-in procedure.

3.2.1 Subsequently, condition the bonding surfaces of the restoration as follows:

a) Glass-ceramic restorations (e.g. IPS Empress®)

- Etch with 5% hydrofluoric acid (e.g. IPS® Ceramic etching gel) for 60 s or according to the instructions of the manufacturer of the restorative material.
- Thoroughly rinse with water spray and dry with oil-free air.

b) Lithium disilicate glass-ceramic restorations (e.g. IPS e.max Press, IPS e.max CAD)

- Etch with 5% hydrofluoric acid (e.g. IPS Ceramic etching gel) for 20 s or according to the instructions of the manufacturer of the restorative material.
- Thoroughly rinse with water spray and dry with oil-free air.

3.2.2 Apply Monobond Plus with a brush or a microbrush to the pre-treated surfaces, let it react for 60 s and then disperse it with a strong stream of air.

3.3 Alternatively, restorations made of glass-ceramics and lithium disilicate glass-ceramics (e.g. IPS Empress, IPS e.max Press, IPS e.max CAD) may be conditioned with Monobond Etch & Prime® as follows, independent of their pre-treatment:

- Following try-in, thoroughly rinse the restoration with water spray and dry with oil- and moisture-free air.
- Apply Monobond Etch & Prime on the bonding surface using a microbrush and agitate it into the surface for 20 seconds. Allow to react for another 40 seconds.
- Thoroughly rinse off Monobond Etch & Prime with water spray until the green colour has disappeared. If any residue remains in the microporosities, the restoration may be cleaned with water in an ultrasonic unit for up to 5 minutes.
- Dry the restoration with a strong stream of oil- and moisture-free compressed air for approx. 10 seconds.

3.4 Composite and fiber-reinforced composite restorations should always be conditioned according to the instructions of the manufacturer of the restorative materials. With Tetric® CAD restorations, please observe the instructions stated in the Tetric CAD Instructions for Use.

4 Pre-treatment of the preparation and application of the adhesive (e.g. Adhese® Universal)

Please observe the Instructions for Use of the adhesive used.

4.1 Conditioning with phosphoric acid gel (optional)

The bond to enamel can be improved by selectively etching the enamel or by applying the "etch & rinse" technique. Unprepared enamel surfaces must be conditioned with phosphoric acid (e.g. Total Etch). Please observe the Instructions for Use for the phosphoric acid gel.

a) Selective enamel etching

Apply phosphoric acid gel (e.g. Total Etch) onto the enamel and allow it to react for 15-30 s. Then rinse thoroughly with a vigorous stream of water for at least 5 s and dry with compressed air until the etched enamel surfaces appear chalky white.

b) Etch & rinse technique

Apply phosphoric acid gel (e.g. Total Etch) onto the prepared enamel first, and then to the dentin. The etchant should be left to react on

the enamel for 15-30 s and on the dentin for 10-15 s. Then rinse thoroughly with a vigorous stream of water for at least 5 s and dry with compressed air until the etched enamel surfaces appear chalky white.

4.2 Application of Adhese Universal

- Starting with the enamel, thoroughly coat the tooth surfaces to be treated with Adhese Universal.
- The adhesive must be scrubbed into the tooth surface for at least 20 s. This time must not be shortened. Applying the adhesive on the tooth surface without scrubbing is inadequate.
- Disperse Adhese Universal with oil- and moisture-free compressed air until a glossy, immobile film layer results.

Important note: Avoid pooling as this may compromise the accuracy of fit of the definitive restoration.

- Light-cure Adhese Universal for 10 s using a light intensity of $\geq 500 \text{ mW/cm}^2$ (e.g. Bluephase® Style).

5 Seating of the restoration

- Variolink Esthetic LC is sensitive to any blue light (operative light, ambient light). Therefore, Variolink Esthetic LC should not be dispensed until right before use. Exposure to intensive light should be avoided during application.

- Apply Variolink Esthetic LC with the application tip or a brush/spatula directly to the internal surface of the restoration.

Note: Direct application of Variolink Esthetic LC to a preparation that has been pre-treated with adhesive may - depending on the adhesive used - lead to an acceleration of the curing process and affect the accuracy of fit of the restoration.

- Seat the restoration and **fix/hold it in place during excess removal**.

- Remove all excess luting material.

- a) Wiping technique

Remove excess material immediately with a micro-brush/brush/foam pellet/dental floss or scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal areas, gingival margins).

- b) Circular technique suitable for veneers / inlays / onlays – Light-curing of excess and subsequent removal

Light-cure excess material with a polymerization light (light intensity $\geq 500 \text{ mW/cm}^2$; e.g. Bluephase Style) for **2 s** at a distance of 10-15 mm by running the light probe along the entire cement line.

In the case of laminate **veneers**, choose a start and end point in the incisal area and light-cure the cement line by moving the polymerization light in a circle in clockwise direction.

In the case of **inlays / onlays**, choose a start and end point in the mesial or distal region and light-cure the cement line by performing a circular movement with the polymerization light.

Thereafter, excess cement is easy to remove with a scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal or gingival margins).

6 Polymerization


- As with all composite systems, Variolink Esthetic is subject to oxygen inhibition. This means that the surface layer (approx. 50 µm) does not polymerize during curing, as it comes in contact with atmospheric oxygen. To prevent this, cover the restoration margins with glycerine gel/air block (e.g. Liquid Strip) immediately after excess removal.
- Polymerize Variolink Esthetic in segments, starting with the proximal margins:

Light intensity	Exposure time per mm ceramic and segment	
500 - 1,000 mW/cm ²	20 s	
$\geq 1,000 \text{ mW/cm}^2$	10 s	e.g. Bluephase Style

- Rinse off glycerine gel/airblock (e.g. Liquid Strip).

7 Finishing of the completed restoration

- Check occlusion and functional movements and make adjustments if necessary.
- Finish the cement lines with finishing diamonds if necessary.
- Smooth out the cement lines using finishing and polishing strips and polish them with suitable polishing instruments (e.g. OptraPol®).
- If necessary, finish the restoration with suitable polishers (ceramic: e.g. OptraFine®; composite resin: e.g. OptraPol).

Warning

- Prevent any contact of uncured Variolink Esthetic with the skin / mucous membrane and eyes.
- Uncured Variolink Esthetic may have a slight irritating effect and cause sensitization to methacrylates.
- Customary medical gloves do not provide protection against the sensitizing effect of methacrylates.

Shelf life and storage

- Store Variolink Esthetic at 2-28 °C/36-82 °F.
- Do not use Variolink Esthetic after the expiry date.
- Do not disinfect syringes with oxidizing disinfectants.
- Close Variolink Esthetic LC syringes immediately after use. Exposure to light causes premature polymerization.
- Expiry date: see note on syringes and packages.

Keep out of the reach of children!

For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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www.ivoclarvivadent.com

Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive, P.O. Box 367, Noble Park, Vic. 3174, Australia

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www.ivoclarvivadent.us

Ivoclar Vivadent Shanghai Trading Co., Ltd.

2/F Building 1, 881 Wuding Road, Jing An District, 200040 Shanghai, China

Tel. +86 21 6032 1657, Fax +86 21 6176 0968

www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

Calle 134 No. 7-B-83, Of. 520, Bogotá, Colombia

Tel. +57 1 627 3399, Fax +57 1 633 1663

www.ivoclarvivadent.co

Ivoclar Vivadent SAS

B.P. 118, 74410 Saint-Jorioz, France

Tel. +33 4 50 88 64 00, Fax +33 4 50 68 91 52

www.ivoclarvivadent.fr

Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2, 73479 Ellwangen, Jagst, Germany

Tel. +49 7961 889 0, Fax +49 7961 6326

www.ivoclarvivadent.de

Ivoclar Vivadent Marketing (India) Pvt. Ltd.

503/504 Raheja Plaza, 15 B Shah Industrial Estate

Veera Desai Road, Andheri (West), Mumbai, 400 053, India

Tel. +91 22 2673 0302, Fax +91 22 2673 0301

www.ivoclarvivadent.in

Ivoclar Vivadent Marketing Ltd.

The Icon, Horizon Broadway BSD, Block M5 No. 1

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www.ivoclarvivadent.com

Ivoclar Vivadent Ltd.

12 Omega St, Rosedale, PO Box 303011 North Harbour

Auckland 0751, New Zealand

Tel. +64 9 914 9999, Fax +64 9 914 9990

www.ivoclarvivadent.co.nz

Ivoclar Vivadent Polska Sp. z o.o.

Al. Jana Pawła II 78, 00-175 Warszawa, Poland

Tel. +48 22 635 5496, Fax +48 22 635 5469

www.ivoclarvivadent.pl

Ivoclar Vivadent LLC

Prospekt Andropova 18 korp. 6/, office 10-06, 115432 Moscow, Russia

Tel. +7 499 418 0300, Fax +7 499 418 0310

www.ivoclarvivadent.ru

Ivoclar Vivadent Marketing Ltd.

Qlaya Main St., Siricon Building No.14, 2nd Floor, Office No. 204

P.O. Box 300146, Riyadh 11372, Saudi Arabia

Tel. +966 11 293 8345, Fax +966 11 293 8344

www.ivoclarvivadent.com

Ivoclar Vivadent S.L.U.

Carretera de Fuencarral nº24, Portal 1 – Planta Baja
28108-Alcobendas (Madrid), Spain

Tel. +34 91 375 78 20, Fax +34 91 375 78 38

www.ivoclarvivadent.es

Ivoclar Vivadent AB

Dalvägen 14, 169 56 Solna, Sweden

Tel. +46 8 514 939 30, Fax +46 8 514 939 40

www.ivoclarvivadent.se

Ivoclar Vivadent Liaison Office

: Tesvikiye Mahallesi, Sakayik Sokak, Nisantas' Plaza No:38/2

Kat:5 Daire:24, 34021 Sisli – Istanbul, Turkey

Tel. +90 212 343 0802, Fax +90 212 343 0842

www.ivoclarvivadent.com

Ivoclar Vivadent Limited

Compass Building, Feldspar Close, Warrens Business Park

Enderby, Leicester LE19 4SD, United Kingdom

Tel. +44 116 284 7880, Fax +44 116 284 7881

www.ivoclarvivadent.co.uk

Ivoclar Vivadent, Inc.

175 Pineview Drive, Amherst, N.Y. 14228, USA

Tel. +1 800 533 6825, Fax +1 716 691 2285

www.ivoclarvivadent.us