# Vivaglass<sup>®</sup> Liner

English

**Brugsanvisning**  
**Käyttöohjeet**  
**Bruksanvisning**  
**Productinformatie**  
**Οδηγίες Χρήσεως**  
**Kullanma Talimatı**  
**Инструкция по применению**  
**Instrukcja stosowania**

**Instructions for Use**  
**Gebrauchsinformation**  
**Mode d'emploi**  
**Istruzioni d'uso**  
**Instrucciones de uso**  
**Instruções de Uso**  
**Bruksanvisning**

- Lyspolymeriserende glasionomercement
- Valokovetteinen lasi-ionomeerisementti
- Lysherdende glasionomercement
- Lichtuithardend glasionomercement
- Φωτοπολυμεριζόμενη υαλοϊονομερής κονία
- Işıklı sertleşen cam iyonomer siman
- Светоотверждаемый стеклоиономерный цемент
- Światłoutwardzalny cement szklano-jonomerowy

- Light-curing glass ionomer cement
- Lichthärtender Glasionomercement
- Ciment Verre Ionomère photopolymérisable
- Cemento vetroionomericofotoindurente
- Cemento de ionomero de vidrio fotopolimerizable
- Cimento de ionômero de vidro fotopolimerizável
- Ljushärdande glasjonomercement

For dental use only.  
Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein  
Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

CE 0123

556281/0310/WE3/G

ivoclar  
vivadent<sup>®</sup>  
clinical

## Instructions for Use

### Description

Vivaglass Liner is a light-curing glass ionomer cement, which is especially suitable for cavity lining.

### Composition

#### Powder

Aluminosilicate glass, initiator and pigments

#### Liquid

HEMA, polyacrylic acid, water, dimethacrylate and initiator

### Indication

Light-curing glass ionomer cement for cavity lining

### Contraindication

Vivaglass Liner is contraindicated

- if the stipulated technique cannot be applied.
- if the patient is known to be allergic to any of the ingredients of Vivaglass Liner.
- if pulp is near the dentin surface or exposed.

### Side effects

In rare cases, ingredients of Vivaglass Liner may cause sensitization. In such cases, the product should cease to be used.

### Interactions

Materials containing eugenol may inhibit the polymerization of Vivaglass Liner.

### Application procedure

#### 1. Cavity preparation

Remove caries. Prepare the cavity in the usual manner.

## 2. Isolation

Isolate teeth using a rubber dam (e.g. OptraDam® Plus) or cotton rolls.

## 3. Cleaning and drying the cavity

Clean the cavity in the customary manner using a water jet and dry it (do not dry out).

## 4. Pulpal protection

When working on deep preparations, cover the areas close to the pulp with a small amount of calcium hydroxide material.

## 5. Mixing

### Powder/Liquid Ratio:

The standard powder/liquid ratio of 1.4 g/1.0 g can be achieved with a level Vivaglass measuring spoon of powder and one drop of liquid. Hold bottle vertically when dosing the liquid to avoid the formation of bubbles. Dispense the required amount of powder and liquid onto the mixing pad. Divide the powder into two equal parts. Mix the first half with the liquid for approx. 5–10 seconds. Add the second half of the powder and mix for another 10–15 seconds. The total mixing time should not exceed 20 seconds. Add more powder to obtain a thicker consistency and to shorten working time.

## 6. Application

Apply Vivaglass Liner into the cavity with a suitable applicator.

## 7. Light curing

Light-cure Vivaglass Liner for 20 seconds at a light intensity of more than 500 mW/cm<sup>2</sup> (e.g. bluephase in the LOP mode). Use a layering technique if lining depth exceeds 2 mm.

## 8. Restoration

Place the composite (e.g. Tetric EvoCeram) or amalgam in the usual manner.

## Note

1. When dosing the liquid, hold the bottle vertically (upside down). To avoid the formation of bubbles, hold bottle 5 cm from the mixing pad with a steady hand and squeeze gently.
2. If the nozzle becomes clogged, wipe it clean with a damp cloth.
3. High temperatures decrease working time, while low ones increase it.
4. Do not mix the components with conventional glass ionomer cements.
5. If etching gel is used to condition the enamel/dentin after having lined the cavity with Vivaglass Liner, wetting of the liner surface with etching gel should be avoided.

## Warning

The contact of uncured Vivaglass Liner with the mucous membrane, skin or eyes may cause irritations and should be avoided. If Vivaglass Liner accidentally comes in contact with the eyes, immediately rinse with copious amounts of water and see a physician. Commercial medical gloves do not provide protection against the sensitizing effects of methacrylates.

## Storage

- Do not use Vivaglass Liner after the expiration date.
- To protect the powder from humidity and prevent the liquid from evaporating, close the respective bottles tightly after use.
- Store packages in use at room temperature (2–28 °C / 36–82 F°).

**Keep material out of children's reach.**

**For use in dentistry only!**

**Date information prepared 03/2010, Rev. 1**

#### Manufacturer

Ivoclar Vivadent AG  
FL-9494 Schaan/Liechtenstein

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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### Ivoclar Vivadent Limited
### Ivoclar Vivadent, Inc.