# Bluephase® Style 20i

The high-performance curing light

The powerful LED  
for every use

**EN Instructions for Use**

Page 2

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Pagina 56

**Appendix**

Page 110

**DE Gebrauchsinformation**

Seite 20

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Página 74

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Page 38

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Página 92

CE 0123

**Rx ONLY**

For dental use only!  
Made in Austria

**Manufacturer**

Ivoclar Vivadent AG  
FL-9494 Schaan/Liechtenstein  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

In USA distributed by  
Ivoclar Vivadent Inc.,  
175 Pineview Drive,  
Amherst, NY 14228, USA

ivoclar  
vivadent®  
clinical

*Dear Customer*

*Optimum polymerization is an important requirement for all light cured materials in order to consistently produce high quality restorations. The polymerization light selected also plays a decisive role in this respect. Therefore, we would like to thank you for having purchased Bluephase® Style 20i.*

*Bluephase Style 20i is a high-quality medical device which has been designed according to the latest standard of science and technology in compliance with the relevant industry standards.*

*These Instructions for Use will help you safely start up the device, make full use of its capabilities and ensure a long service life.*

*Should you have any further questions, please do not hesitate to contact us (see addresses on the reverse page).*


*Your Ivoclar Vivadent Team*

# Table of Contents

<b>1. Product Overview</b>	<b>4</b>
1.1 List of parts	
1.2 Indicators on the charging base	
1.3 Indicators on the handpiece	
1.4 Operating the light	
<b>2. Safety</b>	<b>7</b>
2.1 Intended use	
2.2 Indication	
2.3 Signs & symbols	
2.4 Safety notes	
2.5 Contraindication	
<b>3. Start-Up</b>	<b>10</b>
<b>4. Operation</b>	<b>13</b>
<b>5. Maintenance and Cleaning</b>	<b>16</b>
<b>6. What if ... ?</b>	<b>18</b>
<b>7. Warranty /Procedure in Case of Repair</b>	<b>19</b>
<b>8. Product Specifications</b>	<b>19</b>

# 1. Product Overview

## 1.1 List of parts


## 1.2 Indicators on the charging base

- Indicator is black = Battery is charged
- Indicator lights up in blue with different light intensity = battery is being charged

### 1.3 Indicators on the handpiece


#### Curing time / light intensity

- 20 seconds / 1,200 mW/cm<sup>2</sup>
  - 15 seconds / 1,200 mW/cm<sup>2</sup>
  - 10 seconds / 1,200 mW/cm<sup>2</sup>
  - 5 seconds / 2,000 mW/cm<sup>2</sup>
- } HIGH POWER (LED green)
- } TURBO (LED turquoise)

### 1.4 Operating the light

#### Intuitive two-button operation

With the handpiece switched on, the current charging status is shown on the handpiece as follows:

• **Turquoise (5s) / Green (10s, 15s, 20s):**

**Battery fully charged**

Curing capacity of approximately 20 minutes in the High Power program and approximately 10 minutes in the Turbo program.

• **Orange:**


**Battery weak**

Time/intensity can still be set and a polymerization time of approximately 3 minutes in the High Power program is left. Place the light into the charging base as soon as possible.

• **Red:**


**Battery completely discharged**

The light can no longer be called up and the curing time can no longer be set. However, the handpiece can be used in the Click & Cure corded operation.

**Acoustic signal – Setting the signal volume**

With the light switched off, press the time/program selection button for 5 seconds until a blue light appears. The volume can be adjusted using the time/program selection button. There is a selection of four volume levels. The volume is saved by pressing the Start button.

## 2. Safety

### 2.1 Intended use

Bluephase Style 20i is an LED polymerization light that produces energy-rich blue light. It is used for the polymerization of light-curing dental materials immediately at the dental unit. The intended place of application is in the dental practice, in the medical practice or in the hospital. The intended use also includes the observation of the notes and regulations in these Instructions for Use.

### 2.2 Indication

With its "Polywave<sup>®</sup>" broadband spectrum, Bluephase Style 20i is suitable for the polymerization of all light-curing dental materials curing in the wavelength range of 385–515 nm. These materials include restoratives, bonding agents/adhesives, bases, liners, fissure sealants, temporaries, as well as luting materials for brackets and indirect restorations, such as ceramic inlays.

### 2.3 Signs and symbols

Contraindication

*Symbols on the curing light*

Double insulation  
(device complies with safety class II)

Protection against electrical shock  
(BF type apparatus)

Observe Instructions for Use

Observe Instructions for Use

Caution

The curing light must not be disposed of in the normal domestic waste. Information regarding disposal of the light can be found on the respective national Ivoclar Vivadent homepage.

Recyclable

AC voltage

DC voltage

### 2.4 Safety notes

Bluephase Style 20i is an electronic device and a medical product which is subject to IEC 60601-1 (EN 60601-1) and EMC directives IEC 60601-1-2 (EN60601-1-2) Edition 3.0, as well as the 93/42/EEC Medical Device Directive. The curing light complies with the relevant EU regulations.

CE 0123

The curing light has been shipped from the manufacturer in a safe and technically sound condition. In order to maintain this condition and to ensure risk-free operation, the notes and regulations in these Instructions for Use have to be observed. To prevent damage to equipment and risks for patients, users and third parties, the following safety instructions have to be observed.

### 2.5 Contraindications

Materials, the polymerization of which is activated outside the wavelength range of 385–515 nm (no materials known to date). If you are not sure about certain products, please ask the manufacturer of the corresponding material.

Do not charge or use the appliance near flammable or combustible substances.

Use without light probe.

Use of a light probe other than the one provided in the delivery form.

The use of this device close to other equipment or stacked with it should be avoided because the correct function can be disrupted. If such use is unavoidable, the devices need to be monitored and checked for correct function.

Portable and mobile high-frequency communication devices may interfere with medical equipment. The use of mobile phones during operation is not allowed.

Caution – The use of controls or adjustment devices or performing procedures other than those specified herein may result in hazardous radiation exposure.

### Warning

This unit should not be used near flammable anaesthetics or mixtures of flammable anaesthetics with air, oxygen or nitric oxide.

### Usage and liability

- Bluephase Style 20i must only be employed for the intended use. Any other uses are contraindicated. Do not touch defective, open devices. Liability cannot be accepted for damage resulting from misuse or failure to observe the Operating Instructions.
- The user is responsible for testing Bluephase Style 20i for its use and suitability for the intended purposes. This is particularly important if other equipment is used in the immediate vicinity of the curing light at the same time.
- Use only original spare parts and accessories from Ivoclar Vivadent (see Accessories). The manufacturer does not accept any liability for damage resulting from the use of other spare parts or accessories.
- The light probe is an applied part and may warm up to a maximum of 45 °C (113 °F) at the interface to the handpiece during operation.

### Operating voltage

Before switching on, make sure that

- a) the voltage indicated on the rating plate complies with the local power supply and
- b) the unit has acquired the ambient temperature.

If the battery or power pack are used separately, e.g. during start-up or Click & Cure corded operation, contact with patients or third parties must be prevented. Do not touch the exposed contacts of the connection plug (power pack).

### Rating plate on the charging base

### Assumption of impaired safety

If it has to be assumed that safe operation is no longer possible, the power must be disconnected and the battery removed to avoid accidental operation. This may be the case, for example, if the device is visibly damaged or no longer works correctly. A complete disconnection from the power supply is only ensured when the power cord is disconnected from the power source.

### Eye protection

Direct or indirect exposure of the eyes must be prevented. Prolonged exposure to the light is unpleasant to the eyes and may result in injury. Therefore, using the supplied anti-glare cones is recommended. Individuals who are sensitive to light, who take photosensitizing drugs, have undergone eye surgery, or people who work with the apparatus or in its vicinity for long periods of time should not be exposed to its light and wear protective orange goggles that absorb light below a wavelength of 515 nm. The same is true for patients.

## Battery

Caution: Use only original spare parts, particularly Ivoclar Vivadent batteries and charging bases. Do not short circuit battery. Do not store at temperatures above 40 °C / 104 °F (or 60 °C / 140 °F for a short period). Always store batteries charged. The storage period must not exceed 6 months. May explode if disposed of in fire.

Please note that lithium-polymer batteries may react with explosion, fire and smoke development if handled improperly or mechanically damaged. Damaged lithium-polymer batteries must no longer be used.

The electrolytes and electrolyte fumes released during explosion, fire and smoke development are toxic and corrosive. In case of accidental contact with the eyes or skin, immediately wash with copious amounts of water. Avoid inhalation of fumes. In case of indisposition, see a physician immediately.

## Heat development

As it is the case with all high-performance lights, the high light intensity results in a certain heat development. Prolonged exposure of areas near the pulp and soft tissues may result in irreversible damage. Therefore, this high-performance curing light must only be operated by trained professionals.

The recommended curing times have to be observed. Uninterrupted curing times of more than 5 seconds in the Turbo program on the same tooth surface, as well as direct contact with the gingiva, oral mucous membrane or skin have to be prevented. If the irradiation of soft tissue cannot be avoided, work with reduced light intensity (use High Power program). Polymerize indirect restorations at intermittent intervals of 10 seconds in the High Power program and 1 x 5 seconds each in the Turbo program or use external cooling with an air stream. The instructions regarding curing programs and curing times must be observed (see Selecting the curing program and the curing time). In addition, the light emission window must be placed exactly on the material to be cured at all times (e.g. by holding it in place using a finger).

**After several curing cycles on the same tooth, there is a risk that the pulpa suffers damage caused by the increased temperature!**

### 3. Start-Up

Check the delivery for completeness and any possible transportation damage (see delivery form). If parts are damaged or missing, contact your Ivoclar Vivadent representative.

#### Charging base



Before you switch on the device, make sure that the voltage mentioned on the rating plate complies with your local power supply.

Connect the power cord with the power supply and the power pack. Make sure that the power cord is easily accessible at all times and can easily be disconnected from the power supply. The Power On indicator on the power pack lights up in green and the illuminated ring on the charging base briefly lights up in blue.

Slide the connection plug of the power pack into the socket on the underside of the charging base. Tilt it slightly and apply slight pressure until you hear and feel it snap into place. Place the charging base on a suitable, flat table top.

#### Handpiece

Remove the handpiece from its packaging and clean the light probe and the handpiece (see chapter Maintenance and Cleaning). The light probe can be removed and reattached by slightly rotating it.

After that, mount the anti-glare cone on the light probe.

Gently place the handpiece in the corresponding rest in the charging base without using any force. If a hygiene sleeve is used, please remove it before you charge the battery. If possible, use the light always with a fully charged battery. This will prolong the service life. It is therefore recommended to place the handpiece into the charging base after each patient. If the battery is fully discharged, the charging time is 2 hours.

### Battery


We recommend fully charging the battery before the first use! If the battery is fully charged, it features a curing capacity of approximately 20 minutes in the High Power program and approximately 10 minutes in the Turbo program.

Slide the battery straight into the handpiece until you hear and feel it click into place.

Since the battery is an expendable part, it has to be replaced after its typical life cycle has expired after approximately 2.5 years. See battery label for the age of the battery.

1106000944  
#637 692  
+ -  
dd/mm/yy  
30/01/13

### Battery charging status

The respective charging status is indicated on the charging base as described on Page 6.


### Click & Cure corded operation

Bluephase Style 20i can be used in corded operation at any time, but particularly when the battery is completely empty.


For this purpose, remove the battery from the handpiece by pressing the battery release button.

Then remove the power pack from the underside of the charging base. Do not pull on the power cord.

Insert the connection plug straight into the handpiece until you hear and feel it click into place.

During corded operation, the charging base cannot charge the battery, since it is not connected to a power source.

A complete disconnection from the power supply is only ensured when the power cord is disconnected from the power socket.

## 4. Operation

Disinfect contaminated surfaces of the curing light as well as light probes and anti-glare cones before each use. Additionally, the light probe can be sterilized using the autoclaves intended for this purpose. Furthermore, make sure that the stipulated light intensity permits adequate polymerization. For that purpose, check the light probe for contamination and damage as well as the light intensity at regular intervals.

### Selecting the curing program and time

Bluephase Style 20i is equipped with the following 4 selectable curing times and 2 curing programs for the different indications. Use the time/program selection button to adjust the desired curing time and hence the specified light intensity.

#### HIGH POWER Program, 1,200 mW/cm<sup>2</sup> ± 10%:


By selecting a curing time of 10, 15 or 20 seconds, the High Power program with a light intensity of 1,200 mW/cm<sup>2</sup> is automatically used.

#### TURBO Program, 2,000 mW/cm<sup>2</sup> ± 10%:

By selecting a curing time of 5 seconds, the Turbo program with a light intensity of 2,000 mW/cm<sup>2</sup> is automatically used.

Observe the Instructions for Use of the material applied when selecting the curing time. The curing recommendations for composite materials apply to all shades and, if not mentioned otherwise in the Instructions for Use, to a maximum layer thickness of 2 mm. Generally, these recommendations apply to situations where the emission window of the light probe is placed directly over the material to be polymerized. Increasing the distance between the light source and the material will require the curing time to be extended accordingly. For instance, if the distance to the material is 8 mm, the effective light output is reduced by approx. 50%. In this case, the recommended curing time has to be doubled.

- 1) The information provided herein applies to the 10>8 mm light probe supplied in the delivery form.
- 2) The information regarding heat development and burn hazards must be taken into consideration (see Safety notes)

Curing times		HIGH POWER 1,200 mW/cm <sup>2</sup> ± 10%	TURBO 2,000 mW/cm <sup>2</sup> ± 10%
Restorative materials	<b>Composites</b> • 2 mm <sup>1)</sup> IPS Empress Direct / IPS Empress Direct Flow Tetric EvoCeram / Tetric EvoFlow Tetric / Tetric Basic White	10 seconds	5 seconds
	• 4 mm <sup>2)</sup> Tetric EvoCeram Bulk Fill / Tetric EvoFlow Bulk Fill Tetric Basic White	10 seconds	5 seconds
	Heliomolar / Heliomolar HB / Heliomolar Flow All conventional composites <sup>1)</sup>	15 seconds	2 x 5 seconds
	<b>Compomers</b> <sup>3)</sup> Compoglass F / Compoglass Flow	20 seconds	2 x 5 seconds
Indirect restorations / cementation materials	Variolink Esthetic LC <sup>4)</sup> / Variolink Esthetic DC <sup>5)</sup>	per mm ceramic: 10 seconds per surface	per mm ceramic: 5 seconds per surface
	Multilink Automix <sup>5)</sup> / SpeedCEM Plus <sup>5)</sup>	per mm ceramic: 20 seconds per surface	per mm ceramic: 2 x 5 seconds per surface
Adhesives	Adhese / Adhese Universal ExcITE F / ExcITE F DSC Heliobond / Syntac	10 seconds	Use HIGH POWER program
Temporary materials	Telio CS Link Telio CS Inlay/Onlay Telio Stains Telio Add-On Flow	10 seconds per surface 10 seconds 10 seconds 15 seconds	5 seconds per surface 5 seconds 5 seconds 2 x 5 seconds
Miscellaneous	Helioseal / Helioseal F / Helioseal Clear	20 seconds	Use HIGH POWER program
	Heliosit Orthodontic	10 seconds	2 x 5 seconds
	IPS Empress Direct Color IPS Empress Direct Opaque MultiCore Flow / Multicore HB	20 seconds 20 seconds 20 seconds	5 seconds 2 x 5 seconds 2 x 5 seconds

1) Applies to a maximum layer thickness of 2 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case e.g. with dentin shades)

2) Applies to a maximum layer thickness of 4 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case e.g. with dentin shades)

3) Applies to a maximum layer thickness of 3 mm

4) Applies to light-curing

5) Applies to dual-curing

### **Cure Memory function**

The last settings used, together with the combination of curing program and curing time, are automatically saved.

### **Start**

The light is switched on by means of the start button. It is recommended that the emission window of the light probe is placed directly on the material to be polymerized. Once the selected curing time has elapsed, the curing program is automatically terminated. If desired, the light can be switched off before the set curing time has elapsed by pressing the start button again.

### **Acoustic signals**

Acoustic signals can be heard for the following functions:

- Start (Stopp)
- Every 10 seconds
- Curing time change, you hear the acoustic signal twice when changing to the Turbo program
- Inserting battery
- Error message

### **Light intensity**

The light intensity is maintained at a consistent level during operation. If the supplied 10>8 mm light probe is used, the light intensity has been calibrated to  $2,000 \text{ mW/cm}^2 \pm 10\%$


If another light probe than the one supplied is used, it directly influences the emitted light intensity.

In a light probe with parallel walls (10 mm), the diameter of the light entry and the light emission window is the same. When using focussing light probes (10>8 mm light probe, Pin-Point light probe 6>2 mm), the diameter of the light entry is larger than that of the light emission window. The incident blue light is thus bundled to a smaller area. In this way, the emitted light intensity is increased.

Pin-Point light probes are suitable for spot-on polymerization, e.g. to fix veneers before the removal of excess. For complete polymerization, the light probe must be changed.

## 5. Maintenance and Cleaning

For reasons of hygiene, we recommend using a disposable protective sleeve for each patient. Make sure to fit the protective sleeve closely to the light probe. Disinfect contaminated surfaces of the device and anti-glare cones (FD 366/Dürr Dental, Incidin Liquid/Ecolab). Sterilize the light probe before each use if disposable protective sleeves are not used. Make sure that no liquids or other foreign substances enter the handpiece, charging base and particularly the power pack during cleaning (risk of electrical shock). Disconnect the charging base from the power source when cleaning it.

### Housing

Wipe the handpiece and handpiece holder with a customary aldehyde-free disinfecting solution. Do not clean with highly aggressive disinfecting solutions (e.g. solutions based on orange oil or with an ethanol content of more than 40%), solvents (e.g. acetone), or pointed instruments, which may damage or scratch the plastic. Clean dirty plastic parts with a soapy solution.

### Light probe

Before cleaning and/or disinfecting the light probe, pretreat it. This applies to both automated and manual cleaning and disinfection:

- Remove substantial contamination immediately after use or 2 hours after that at the latest. For this purpose, thoroughly rinse the light probe under running water (for at least 10 seconds). Alternatively, use a suitable aldehyde-free disinfecting solution to prevent the fixation of blood.
- To remove contamination manually, use a soft brush or soft cloth. Partially polymerized composite can be removed with alcohol and a plastic spatula, if necessary. Do not use sharp or pointed objects, as they may scratch the surface.

### Cleaning and disinfection

For cleaning, immerse the light probe in a cleaning solution and ensure that it is sufficiently covered with liquid (ultrasound or careful brushing with a soft brush may support the effect). A neutral-enzymatic cleaning agent is recommended.

When cleaning and disinfecting, please make sure that the agents used are free of:

- organic, mineral and oxidizing acids (the minimum admissible pH value is 5.5)
- alkaline solution (the maximum admissible pH value is 8.5)
- oxidizing agent (e.g. hydrogen peroxide)

Afterwards, remove the light probe from the solution and thoroughly rinse it under running water (for at least 10 seconds). Cleaning in a thermal disinfectant is an effective alternative.

### **Sterilization**

Thorough cleaning and disinfecting is imperative to ensure that the subsequent sterilization is effective. Use only autoclave sterilization for this purpose. The sterilization time (exposure time at sterilization temperature) is 4 minutes at 134 °C (273 °F); pressure should be 2 bar (29 psi). Dry the sterilized light probe using either the special drying program of your steam autoclave or hot air. The light probe has been tested for up to 200 sterilization cycles.

After that, check the light probe for damage. Hold it against light. If individual segments appear black, glass fibres are broken. If this is the case, replace the light probe with a new one.

### **Disposal**

The curing light must not be disposed of as urban waste. Dispose unserviceable batteries and polymerization lights according to the corresponding legal requirements in your country. Batteries must not be incinerated.

## 6. What if ...?

Indicator	Causes	Error rectification
<p>All LED's orange</p> 	<p>The device is overheated.</p>	<p>Allow the device to cool down and try again after a certain time. If the error persists, please contact your dealer or your local Service Centre.</p>
<p>All LED's red</p> 	<p>Electronic component of the handpiece is defective.</p>	<p>Remove and reinsert the battery. If the error persists, please contact your dealer or your local Service Centre.</p>
<p>The charging base is not illuminated during charging</p>	<ul style="list-style-type: none"> <li>– Power pack not connected or defective</li> <li>– Battery fully charged</li> </ul>	<p>Check if the power pack is correctly positioned in the charging base or if the power pack is connected to the power supply by means of the power cord (LED on the power pack lights up in green if it works correctly).</p>
<p>No indicator activity of the device with the battery in place</p>	<p>Battery empty</p>	<p>Place the device in the charging base and charge for at least 2 hours.</p>
	<p>Battery contacts dirty</p>	<p>Remove battery and clean the battery contacts.</p>

## 7. Warranty / Procedure in Case of Repair

The warranty period for Bluephase Style 20i is 3 years from the date of purchase (battery: 1 year). Malfunctions resulting from faulty material or manufacturing errors are repaired free of charge during the warranty period. The warranty does not provide the right to recover any material or non-material damage other than the ones mentioned. The apparatus must only be used for the intended purposes. Any other uses are contraindicated. The manufacturer does not accept any liability resulting from misuse and warranty claims cannot be accepted in such cases. This is particularly valid for:

- Damage resulting from improper handling, especially incorrectly stored batteries (see Product Specifications: Transportation and storage conditions).
- Damage to components resulting from wear under standard operating conditions (e.g. battery).
- Damage resulting from external influences, e.g. blows, drop to the floor.
- Damage resulting from incorrect set-up or installation.
- Damage resulting from connecting the unit to a power supply, the voltage and frequency of which do not comply with the ones stated on the rating plate.
- Damage resulting from improper repairs or modifications that have not been carried out by certified Service Centres.

In case of a claim under warranty, the complete apparatus (handpiece, charging base, power cord and power pack) must be returned, carriage paid, to the dealer or directly to Ivoclar Vivadent, together with the purchase document. Use the original packaging with the corresponding cardboard inserts for transportation. Repair work may only be carried out by a certified Ivoclar Vivadent Service Centre. In case of a defect that cannot be rectified, please contact your dealer or your local Service Centre (see addresses on the reverse side). A clear description of the defect or the conditions under which the defect occurred will facilitate locating the problem. Please enclose this description when returning the apparatus.

## 8. Product Specifications

Light source	Ivoclar Vivadent Polywave® LED
Wavelength range	385–515 nm
Light intensity	TURBO program: 2,000 mW/cm <sup>2</sup> ± 10 % HIGH POWER program: 1,200 mW/cm <sup>2</sup> ± 10 %
Operation	3 min on / 7 min off (intermittent)
Light probe	10>8 mm black, autoclavable
Anti-glare cones	autoclavable
Signal transmitter	acoustic signal after 5 or every 10 seconds and every time the start button or the time/ program selection button is pressed
Dimensions of the handpiece (without light probe)	L = 180 mm, B = 30 mm, H = 30 mm
Weight of the handpiece	120 g (including battery and light probe)
Operating voltage handpiece	3.7 VDC with battery 5 VDC with power pack
Operating voltage Charging Base	5 VDC
Power supply	Input: 100–240 VAC, 50–60 Hz 310 mA Output: 5 VDC / 2 A Manufacturer Friwo Type: FW7401M/05
Operating conditions	Temperature +10 °C to +30 °C (+ 50 °F to + 86 °F) Relative humidity 30% to 75% Ambient pressure 700 hPa to 1060 hPa
Dimensions of the charging base	D = 125 mm, H = 70 mm
Weight of the charging base	195 g
Charging time	Approx. 2 hours (with the battery empty)
Power supply of the handpiece	Li-Po battery (approx. 20 min. with a new, fully charged battery in the High Power program)
Transportation and storage conditions	Temperature –20 °C to +60 °C (+4 °F to +140 °F) Relative humidity 10% to 75% Ambient pressure 500 hPa to 1060 hPa The curing light has to be stored in closed, roofed rooms and must not be exposed to severe jarring.  Battery: – Do not store at temperatures above 40 °C / 104 °F (or 60 °C / 140 °F for a short period). Recommended storage temperature 15–30 °C (59–86 °F). – Keep the battery charged and store no longer than 6 months.
Delivery form	1 Charging base with power cord and power pack 1 Handpiece 1 Light probe 10>8 mm, black 1 Anti-glare shield 3 Anti-glare cones 1 Pack of sleeves (1 x 50 pcs.) 1 Instructions for Use

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# Appendix

Bluephase Style 20i is tested according to IEC 60601-1-2, Edition 3.0.

Medical electrical devices are subject to particular preventive actions according to EMC rules and must be installed and operated according to the EMC guidelines in the accompanying documents.

## Guidance and manufacturer's declaration – electromagnetic emission

The following tables are guidelines according to the 3rd edition of the medical standard IEC 60601-1-2.

Bluephase Style 20i is intended for use in the electromagnetic environment specified below. The customer or the user of Bluephase Style 20i should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	Bluephase Style 20i uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	Bluephase Style 20i is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	N/A
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	N/A

Table: According to IEC 60601-1-2, Edition 3.0

## Guidance and manufacturer's declaration – electromagnetic immunity

Bluephase Style 20i is intended for use in the electromagnetic environment specified below. The customer or the user of Bluephase Style 20i should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of typical commercial or dental environment.
Surge IEC 61000-4-5	± 1 kV line - line ± 2 kV line - earth	± 1 kV line - line no prot. earth	Mains power quality should be that of typical commercial or dental environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	Mains power quality should be that of typical commercial or dental environment. If the user of Bluephase Style 20i requires continued operation during power mains interruptions, it is recommended that Bluephase Style 20i be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic Field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or dental environment.

Table: According to IEC 60601-1-2, Edition 3.0

NOTE:  $U_T$  is the a.c. mains voltage prior to application of the test level.

## Guidance and manufacturer's declaration - electromagnetic immunity

Bluephase Style 20i is intended for use in the electromagnetic environment specified below. The customer or the user of Bluephase Style 20i should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should not be used closer to any part of Bluephase Style 20i, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			<b>Recommended separation distance</b>
Conducted RF IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	10 V	$d = 0.35 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m	$d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz $d = 0.70 \sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol: 

Table: According to IEC60601-1-2, Edition 3.0

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered.

If the measured field strength in the location in which Bluephase Style 20i is used exceeds the applicable RF compliance level above, Bluephase Style 20i should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating Bluephase Style 20i.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strength should be less than 10 V/m.

## Recommended separation distances between portable and mobile RF communications equipment and the "Bluephase Style 20i"

Bluephase Style 20i is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Bluephase Style 20i can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Bluephase Style 20i as recommended below, according to the maximum output power of the communication equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 0.35 \sqrt{P}$	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
0.01	0.035	0.12	0.23
0.1	0.11	0.38	0.73
1	0.35	1.2	2.3
10	1.1	3.8	7.3
100	3.5	12	23

Table: According to IEC 60601-1-2, Edition 3.0

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 4: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

# Ivoclar Vivadent – worldwide

## **Ivoclar Vivadent AG**

Bendererstrasse 2  
9494 Schaan  
Liechtenstein  
Tel. +423 235 35 35  
Fax +423 235 33 60  
www.ivoclarvivadent.com

## **Ivoclar Vivadent Pty. Ltd.**

1 – 5 Overseas Drive  
P.O. Box 367  
Noble Park, Vic. 3174  
Australia  
Tel. +61 3 9795 9599  
Fax +61 3 9795 9645  
www.ivoclarvivadent.com.au

## **Ivoclar Vivadent GmbH**

Tech Gate Vienna  
Donau-City-Strasse 1  
1220 Wien  
Austria  
Tel. +43 1 263 191 10  
Fax: +43 1 263 191 111  
www.ivoclarvivadent.at

## **Ivoclar Vivadent Ltda.**

Alameda Caiapós, 723  
Centro Empresarial Tamboré  
CEP 06460-110 Barueri – SP  
Brazil  
Tel. +55 11 2424 7400  
Fax +55 11 3466 0840  
www.ivoclarvivadent.com.br

## **Ivoclar Vivadent Inc.**

1-6600 Dixie Road  
Mississauga, Ontario  
L5T 2Y2  
Canada  
Tel. +1 905 670 8499  
Fax +1 905 670 3102  
www.ivoclarvivadent.us

## **Ivoclar Vivadent Shanghai Trading Co., Ltd.**

2/F Building 1, 881 Wuding Road,  
Jing An District  
200040 Shanghai  
China  
Tel. +86 21 6032 1657  
Fax +86 21 6176 0968  
www.ivoclarvivadent.com

## **Ivoclar Vivadent Marketing Ltd.**

Calle 134 No. 7-B-83, Of. 520  
Bogotá  
Colombia  
Tel. +57 1 627 3399  
Fax +57 1 633 1663  
www.ivoclarvivadent.co

## **Ivoclar Vivadent SAS**

B.P. 118  
F-74410 Saint-Jorioz  
France  
Tel. +33 4 50 88 64 00  
Fax +33 4 50 68 91 52  
www.ivoclarvivadent.fr

## **Ivoclar Vivadent GmbH**

Dr. Adolf-Schneider-Str. 2  
D-73479 Ellwangen, Jagst  
Germany  
Tel. +49 7961 889 0  
Fax +49 7961 6326  
www.ivoclarvivadent.de

## **Wieland Dental + Technik**

**GmbH & Co. KG**  
Lindenstrasse 2  
75175 Pforzheim  
Germany  
Tel. +49 7231 3705 0  
Fax +49 7231 3579 59  
www.wieland-dental.com

## **Ivoclar Vivadent Marketing**

**(India) Pvt. Ltd.**  
503/504 Raheja Plaza  
15 B Shah Industrial Estate  
Veera Desai Road, Andheri (West)  
Mumbai, 400 053  
India  
Tel. +91 22 2673 0302  
Fax +91 22 2673 0301  
www.ivoclarvivadent.in

## **Ivoclar Vivadent s.r.l.**

Via Isonzo 67/69  
40033 Casalecchio di Reno (BO)  
Italy  
Tel. +39 051 6113555  
Fax +39 051 6113565  
www.ivoclarvivadent.it

## **Ivoclar Vivadent K.K.**

1-28-24-4F Hongo  
Bunkyo-ku  
Tokyo 113-0033  
Japan  
Tel. +81 3 6903 3535  
Fax +81 3 5844 3657  
www.ivoclarvivadent.jp

## **Ivoclar Vivadent Ltd.**

12F W-Tower  
54 Seocho-daero 77-gil, Seocho-gu  
Seoul, 06611  
Republic of Korea  
Tel. +82 2 536 0714  
Fax +82 2 596 0155  
www.ivoclarvivadent.co.kr

## **Ivoclar Vivadent S.A. de C.V.**

Calzada de Tlalpan 564,  
Col Moderna, Del Benito Juárez  
03810 México, D.F.  
México  
Tel. +52 55 50 62 10 00  
Fax +52 55 50 62 10 29  
www.ivoclarvivadent.com.mx

## **Ivoclar Vivadent BV**

De Fruittuinen 32  
2132 NZ Hoofddorp  
Netherlands  
Tel. +31 23 529 3791  
Fax +31 23 555 4504  
www.ivoclarvivadent.com

## **Ivoclar Vivadent Ltd.**

12 Omega St, Rosedale  
PO Box 303011 North Harbour  
Auckland 0751  
New Zealand  
Tel. +64 9 914 9999  
Fax +64 9 914 9990  
www.ivoclarvivadent.co.nz

## **Ivoclar Vivadent Polska Sp. z o.o.**

Al. Jana Pawla II 78  
00-175 Warszawa  
Poland  
Tel. +48 22 635 5496  
Fax +48 22 635 5469  
www.ivoclarvivadent.pl

## **Ivoclar Vivadent Marketing Ltd.**

Prospekt Andropova 18 korp. 6/  
office 10-06  
115432 Moscow  
Russia  
Tel. +7 499 418 0300  
Fax +7 499 418 0310  
www.ivoclarvivadent.ru

## **Ivoclar Vivadent Marketing Ltd.**

Qlaya Main St.  
Siricon Building No.14, 2<sup>nd</sup> Floor  
Office No. 204  
P.O. Box 300146  
Riyadh 11372  
Saudi Arabia  
Tel. +966 11 293 8345  
Fax +966 11 293 8344  
www.ivoclarvivadent.com

## **Ivoclar Vivadent S.L.U.**

Carretera de Fuencarral n°24  
Portal 1 – Planta Baja  
28108-Alcobendas (Madrid)  
Spain  
Tel. +34 91 375 78 20  
Fax: +34 91 375 78 38  
www.ivoclarvivadent.es

## **Ivoclar Vivadent AB**

Dalvägen 14  
S-169 56 Solna  
Sweden  
Tel. +46 8 514 939 30  
Fax +46 8 514 939 40  
www.ivoclarvivadent.se

## **Ivoclar Vivadent Liaison Office**

: Tesvikiye Mahallesi  
Sakayik Sokak  
Nisantas' Plaza No:38/2  
Kat:5 Daire:24  
34021 Sisli – Istanbul  
Turkey  
Tel. +90 212 343 0802  
Fax +90 212 343 0842  
www.ivoclarvivadent.com

## **Ivoclar Vivadent Limited**

Ground Floor Compass Building  
Feldspar Close  
Warrens Business Park  
Enderby  
Leicester LE19 4SE  
United Kingdom  
Tel. +44 116 284 7880  
Fax +44 116 284 7881  
www.ivoclarvivadent.co.uk

## **Ivoclar Vivadent, Inc.**

175 Pineview Drive  
Amherst, N.Y. 14228  
USA  
Tel. +1 800 533 6825  
Fax +1 716 691 2285  
www.ivoclarvivadent.us