# Colado® CAD Ti5

## EN Instructions for Use

- Titanium disc for the CAD/CAM technology

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### Product Description

Colado® CAD Ti5 are milling discs made from titanium alloy, titanium grade 5 (Ti Al6 V4), for the production of dental restorations using CAD/CAM technology. Colado CAD Ti5 fulfills the requirements according to ISO 22674, type 4.

The milling discs can be processed using conventional milling machines (please observe the manufacturer's instructions).

**Note:** When using Colado CAD Ti5 8 mm, check the insertion holder (distance ring).

### Indications

Colado CAD Ti5 is suitable for frameworks with thin cross-sections which are exposed to very high loads, e.g. clasps, veneered crowns, large spanned bridges or bridges with thin cross-sections, bars, fixtures, implant-supported superstructures according to ISO 22674, type 4.

In the anterior region, the bridge frameworks can be made with up to three consecutive pontics. The connector cross-section must not be less than 6 mm². In the posterior region, the bridge frameworks can be made with up to three consecutive pontics. The connector cross-section must not be less than 9 mm².

### Contraindications

If the patient is known to be allergic or sensitive to any components of this alloy, the material must not be used to fabricate restorations. Different types of alloy in the same oral environment can lead to a galvanic reaction.

### Processing Colado® CAD Ti5 after milling

#### 1. Separating the frameworks from the disc

The milled restorations may be separated from the milling disc using cross-cut tungsten carbide burs or separating discs suitable for processing titanium.

**Note:** Do not use any residue material from the disc as a casting alloy!

#### 2. Finishing / Cleaning

Thoroughly clean the separated restorations with hot steam and subsequently degrease with ethyl alcohol in an ultrasonic bath. The restorations can be finished using clean tungsten carbide burs suitable for processing titanium so that no sharp angles or edges remain. Using the bur, grind the metal in one direction only in order to avoid overlapping. Please observe the manufacturer's recommended maximum rpm for these instruments. Next, sandblast the surface using pure aluminum oxide (approx. 180 µm) at 2–3 bar pressure. Rinse off the jet medium under running water and then thoroughly clean the restoration with hot steam. After that, contamination of the object must be prevented (e.g. do not touch it with your fingers). When grinding the metal, do not inhale the grinding dust! Never use hydrofluoric acid!

#### 3. Veneering

The metal wall thickness (single crowns 0.3 mm/ abutment crowns 0.5 mm) should always be observed. Colado CAD Ti5 frameworks can be veneered using the laboratory composite SR Nexo®/SR® Link (Ivoclar Vivadent).

#### Note:

The alloy must not be exposed to temperatures above 800°C (1472 °F) since this can result in a transformation of the lattice structure in the alloy micro-structure.

#### Caution:

Metal vapours and metal dust are harmful if inhaled; therefore, always use extraction equipment and/or suitable protective masks! When processing titanium and titanium alloys, there is a fire hazard. Therefore always use a cooling liquid when processing Colado CAD Ti5 in the milling machine.

#### Storage

No specific measures required.

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

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