# IPS® PressVEST Premium

**EN Instructions for Use**  
Page 2

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**NO Bruksanvisning**  
Side 38

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For dental use only.

## Description

IPS® PressVest Premium is a universal, phosphate bonded investment material (type 1, class 1 and 2) for use with Ivoclar Vivadent press ceramics. It can be used with the speed heating procedure as well as with conventional heating methods.

## Definition

### *Speed heating procedure*

- The set investment ring is placed directly into a furnace which has been preheated to the final temperature and then transferred to the press furnace after a defined time.

### *Conventional heating procedure*

- The set investment ring is placed into a cold preheating furnace (room temperature) and heated to the end temperature according to a designated heating process. After a defined time the investment ring is then transferred into the press furnace.

## Indication

IPS PressVest Premium can be used with the following press ceramics in the Ivoclar Vivadent press furnaces:

- IPS e.max® Press
- IPS e.max® ZirPress
- IPS InLine® PoM (Press-on Metal)
- IPS Empress® Esthetic

## Contraindications

Any application not included within the indications

## Storage

- Store at 12–28 °C / 54–82 °F
- Store powder in a dry place
- The liquid must not be exposed to temperatures below +5 °C / 41 °F (frost susceptible!). **If the liquid should freeze, it will be unusable (crystallization).**

## Working temperature

- Room temperature 18 °C – max. 23 °C / 64 °F – max. 73 °F
- Any other processing temperature decisively affects the setting behaviour.
- A lower processing temperature delays the chemical process.
- A higher processing temperature accelerates the chemical process. (With temperatures from 23 – 28 °C / 73 °F – 82 °F the liquid and the water can be cooled to 12 °C / 54 °F).

## Mixing ratio

By diluting the investment liquids the expansion of the investment material can be controlled and therefore influence the fit of the pressed restoration on the die or preparation. The liquid concentration is shown in %.

**Important: The amount of mixing liquid (investment liquid + distilled water) is related to the amount of investment powder and must be strictly adhered to.**

- 100 g powder: 26 ml mixing liquid (investment material liquid + distilled water)
- 200 g powder: 52 ml mixing liquid (investment material liquid + distilled water)
- 300 g powder: 78 ml mixing liquid (investment material liquid + distilled water)

The following information on liquid concentrations are guidelines and may be adjusted depending on the object size and processing conditions. Do not dilute the liquid to below 50%!

Material	Indication	Concentration Mixing liquid approx.	100 g powder	200 g powder	300 g powder
			Liquid : dist. water	Liquid : dist. water	Liquid : dist. water
<b>IPS e.max® Press</b>	Crown, veneer	70%	18 ml : 8 ml	36 ml : 16 ml	–
<b>IPS Empress® Esthetic</b>	Inlay / Onlay	60%	16 ml : 10 ml	31 ml : 21 ml	–
<b>IPS e.max® Press</b>	3-unit bridge	70%	–	36 ml : 16 ml	–
	Hybrid abutment, hybrid abutment crown	85%	22 ml : 4 ml	44 ml : 8 ml	–
<b>IPS e.max® ZirPress</b>	Overpressing zirconium oxide frameworks	70%	18 ml : 8 ml	36 ml : 16 ml	54 ml : 24 ml
<b>IPS InLine® PoM</b>	Overpressing metal frameworks	70%	18 ml : 8 ml	36 ml : 16 ml	54 ml : 24 ml

### Mixing

Premixing by hand	20 – 30 seconds
Vacuum mixing time	90 seconds
Speed	approx. 350 rpm
Processing time	approx. 5 – 6 minutes

Dependable results can only be achieved by following the processing instructions and maintaining consistent conditions. Please observe the following instructions:

- Wax the restoration(s) onto the selected investment ring base.
- Remove excess wax separator on the restorations before investing. Wax isolation can react with the investment material and have a negative effect on the pressed results.
- Do not use a wax wetting agent (surfactant) and remove excess separating liquid before investing. Reaction with the investment material is possible.
- Mix max. 400 g of investment material in one go.
- Observe the ratio of powder to liquid exactly.
- Set the concentration of the mixing liquid – according to the material and indication – by adjusting the measurements of investment liquid and distilled water. The higher the water content in the mixing liquid, the lower the setting expansion. In addition, the investment material will be weakened in strength. Do not dilute the mixing liquid to below 50%.
- Reseal the liquid bottle well after use, risk of crystal formation.
- Pour the mixing liquid into a clean mixing beaker. Only process the investment material using clean, dry instruments and do not use plaster mixing beakers.
- Carefully cut the powder bag open and pour the contents into the liquid.
- Mix the investment material well using a spatula for 20 – 30 seconds until all the powder is mixed with the mixing liquid.
- Machine mix with a vacuum mixing device. Check the function of the vacuum mixing device regularly.
- The processing time depends on the material temperature, the amount, the stirring time and the intensity of mixing. A higher material temperature or longer mixing time shortens the processing time.
- It is essential to adhere to the times indicated (setting time / holding time).

### Investing / preheating

- Carry out the investing on a vibrator unit under gentle vibration.
- Allow the invested ring to set without vibration.
- Do not touch the handle of the ring whilst the investment material is setting.
- Do not use the pressure investment method as it changes the product properties.
- When using the conventional method, preheat the investment ring at the latest after 12 hours. Do not invest over the weekend (risk of efflorescence, cracks or press flash).

## Setting time and transfer into the preheating furnace

Speed heating procedure	Conventional heating procedure
30 minutes (min. 30 minutes to max. 45 minutes)	min. 30 minutes to max. 12 hrs

### Factors which influence the processing time and setting expansion of the investment material:

- Temperature of the mixing beaker, investment material and liquid
- Mixing duration and intensity
- Mixing amount

### Preheating / holding times

	Speed heating procedure	Conventional heating procedure
Starting temperature	850 °C / 1562 °F	Room temperature
Temperature increase	–	3 °C / 5.4 °F / min
Holding temperature / time	–	280 °C / 536 °F, 45 min
Holding time after the final temperature of 850 °C / 1562 °F is reached	100 g investment ring – min. 45 minutes	100 g investment ring – min. 45 minutes
	200 g investment ring – min. 60 minutes	200 g investment ring – min. 60 minutes
	300 g investment ring – min. 90 minutes	300 g investment ring – min. 90 minutes

- When placing speed investment rings in the preheating furnace, make sure that the furnace temperature does not drop substantially. The stipulated holding time counts from the point when the preheating temperature of 850 °C / 1562 °F has been reached again.
- If several speed investments are to be conducted (more than two 200 g rings), they should be invested consecutively and placed into the preheating furnace at intervals of at least 20 minutes.
- If more than 2 speed investment rings are preheated at the same time, the holding time at the final temperature should be extended by 15 minutes for each additional ring.

### Important information

- Always remove pressed investment rings from the press furnace immediately after the end of the program and place them on a grid for cooling.
- When divesting, ensure that the coarse sand blasting procedure (at 4 bar pressure) is only used to remove the excess investment material and is not used on the ceramic object itself.

### To ensure thorough heating of the investment rings in the preheating furnace, the following points must be observed:

- Load the preheating furnace to maximally half of the available utility space, always place the investment rings in the rear part of the preheating furnace.
- Always place the investment rings into the furnace with the opening facing down, if possible at a 45° angle.
- Do not place the investment rings in direct contact with each other. Allow air circulation.
- Place the investment rings which are to be preheated using the speed heating procedure into a furnace preheated to 850 °C / 1562 °F.
- Place the investment rings as quickly as possible in the preheating furnace when using the speed method. Make sure that the furnace temperature does not drop substantially. If several speed investments are to be conducted (more than two 200 g rings), they should be invested consecutively and placed into the preheating furnace at intervals of at least 20 minutes.
- If more than 2 speed investment rings are preheated at the same time, the holding time at the final temperature should be extended by 15 minutes for each additional ring.
- Place investment rings which are to be heated using the conventional preheating procedure into a cold preheating furnace and begin heating from room temperature. Ensure that the investment rings do not dry out before the preheating process starts.

## Safety notes

### DANGER!


- This material contains silica that causes damage to lungs through prolonged or repeated exposure via inhalation.
- Do not breathe dust.
- In case of inadequate ventilation, wear respiratory protection.
- Dispose of contents/containers in accordance with local regulations.
- Do not open the furnace during the first 20 minutes of preheating when processing the material using the speed method. Burn hazard caused by wax fumes that might ignite in the air!

### Physical data

According to ISO/FDIS 15912:2015

- |                             |            |
|-----------------------------|------------|
| - Flow characteristics:     | 115 mm     |
| - Initial set:              | 10 Minutes |
| - Compressive strength:     | 6 MPa      |
| - Linear thermal expansion: | 0.8 %      |

Curve of typical linear thermal expansion:

### Composition

Quartz, cristobalite, aluminium oxide, binding agent: ammonium phosphate, magnesium oxide

This material has been developed solely for use in dentistry and must be processed according to the Instructions. Liability cannot be accepted for damages resulting from misuse or failure to observe the Instructions. The user is solely responsible for testing the material for its suitability for any purpose not explicitly stated in the Instructions.

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**Ivoclar Vivadent Limited**
**Ivoclar Vivadent, Inc.**