Instructions for Use EN Gebrauchsinformation DE Mode d'emploi FR Istruzioni d'uso IT Instrucciones de uso ES Instruções de Uso PT Bruksanvisning SV Bruksanvisning DA Käyttöohjeet FI

■ **PRODUCT DESCRIPTION**
Au-based dental metal-ceramic alloy, Type 4

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■ **WAXING/MODELLATION**
Design the framework in a reduced anatomic shape taking the planned veneer into consideration. Single wax thickness of 0.3 mm. Abutment corners require a minimum thickness of 0.5 mm. Ensure the framework design provides adequate support for the veneering material. Avoid sharp angles. Connectors must have the required dimensions to provide resistance to deformation. Create large surface areas for alignment and soldering, with a gap of 0.05–0.2 mm.

■ **SPRUNG**
Provide the modeled single-tooth restorator or bridge framework with sprues of a suitable size. In general the reservoir, sprue leads, and connector sprues, whether per shaped or traditional, must be sized according to the specific technique used. When using the direct or indirect technique be sure that the reservoir is positioned in the heat center. The connector sprues between the reservoir and the casting should be a maximum of 2.5–3.0 mm in length and width. The wax pattern including the sprues must be weighed in grams in order to determine the needed amount of alloy. Wax conversion formula: wax weight (gram) x alloy density = grams of alloy required.

■ **INVESTING**
Use a phosphate-bonded investment material. Follow the manufacturer's instructions.

■ **PREHEATING/ BURN-OUT**
Recommended burn-out temperature: 850 °C

■ **MELTING AND CASTING**
Torch: Propane 0.35 bar Oxygen 0.7 bar
Cast into the well of the casting machine. It is recommended to use a separate and clean carbon/ceramic crucible for each alloy. Preheat the ceramic crucible in the burnout furnace. The recommended ratio of used material to new material is 1:1. Do not use flux.
Casting Temperature: 1260 °C

■ **FRAMEWORK FINISHING**
After bench cooling, carefully divert and clean the casting with aluminum oxide (Al₂O₃). Do not use a hammer for diverting. Finish the casting with carbide burrs and/or with carbide-bonded grinding instruments. Blast the surface with 50–110 micron aluminum oxide (Al₂O₃) at 2.0 bar. Subsequently, steam clean or ultrasonic clean with distilled water or ethanol and dry the framework.

■ **OXIDATION**
Place the framework on the firing tray providing adequate support. To achieve a uniform result follow the oxidation cycle.
Temperature: 930 °C. Holding time: 5 min; Vacuum: No
If the oxide layer is stained, grind and blast the surface again. Repeat the oxide firing. In the case of Zn-containing alloys, the additional picking of the oxide layer in 10% sulfuric acid or similar pickling solutions is recommended. Note: Always observe proper safety procedures when handling acid. Before application clean framework carefully with distilled water or ethanol and dry the framework.

■ **HIGHEST RECOMMENDED FIRING TEMPERATURE:** 950 °C

■ **HEAT TREATMENT**
Hardening: 450 °C for 15 min; bench cool

■ **SOLDERING AND LASER SWEISSING**
The soldering gap should not be wider than the thickness of the soldering material. Allow the soldered casting to cool slowly. Use flux sparingly.
Pre Solder: BioPorta Solder 1020 Flux: High Fusing Bondal Flux
Post Soldering: BioPorta Solder W-2 Flux: Bondal Flux
Laser Soldering Wire: BioPorta G

■ **POLISHING**
Carefully remove any oxide and flux residue. Smooth the metal surfaces with rubber polishers. Polish to a high gloss finish using polishing paste. Subsequently, clean using ultrasonic cleaning equipment or careful steam cleaning.

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