## Instructions for Use EN

**PRODUCT DESCRIPTION**  
Au-based dental ceramic / metal-ceramic alloy, Type 4

**INDICATIONS\***  
Inlays, Onlays, 3/4 Crowns, Crowns, Telescope Crowns, Conus Crowns, Bridges, Wide Bridges, Cast Posts /Cores, Bars, Attachments, Implant Retained Superstructures, Partial Dentures

**WAXING/MODELLATION**  
Full contour or in a reduced anatomical shape taking the planned veneer into consideration. For Lab Composite veneering material use mechanical retention. Single crowns require a minimum thickness of 0.3 mm. Abutment crowns require a minimum thickness of 0.5 mm. Ensure the wax is of adequate strength to support the veneering material. Avoid sharp angles. Connectors must have the required dimensions to provide resistance to deformation. Create large surface areas for planned soldering, with a gap of 0.05 – 0.2 mm.

**SPRUNG**  
Use the modeled single-tooth restoration or bridge framework with sprues of a suitable size. In general the reservoir, sprue leads, and connector sprues, whether per shaped or traditional, must be sized according to the specific technique used. When using the direct or indirect technique be sure that the reservoir is positioned in the best center. The reservoir and the casting should be a maximum of 2.5 – 3.0 mm in length and width. The wax pattern including the sprues must be weighed in grams in order to determine the needed amount of wax. Wax conversion formula: wax weight (gram) x alloy density = grams of alloy required.

**INVESTING**  
Use a phosphate-bonded investment material. Follow the manufacturer's instructions.

**PREHEATING / BURN-OUT**  
Recommended burn-out temperature: 800 °C

**MELTING AND CASTING**  
Torch: Propano 0.15 bar, Oxigeno 0.35 bar  
Other specifics may be required by the type of casting machine. It is recommended to use a separate and clean investment crucible for each alloy and the casting should be a maximum of 2.5 – 3.0 mm in length and width. The wax pattern including the sprues must be weighed in grams in order to determine the required amount of wax. Wax conversion formula: wax weight (gram) x alloy density = grams of alloy required.

**CASTING TEMPERATURE:** 1200 °C

**FRAMWORK FINISHING**  
After being casted, carefully clean the casted with aluminum oxide (AlO<sub>3</sub>). Do not use a hammer for finishing. Finish the casting with carbide burrs and/or with ceramic-bonded grinding instruments. Blast the surface with 50 – 110 micron aluminum oxide (AlO<sub>3</sub>) at 2.0 bar. For veneering with Lab Composite material the framework must be finished and polished. Subsequently, steam clean or ultrasonic clean with distilled water to ethanol and dry the framework.

**OXIDATION**  
Place the framework on the firing tray providing adequate support. To achieve a uniform result follow the oxidation cycle.

**TEMPERATURE:** 800 °C. **Holding time:** 5 min; **Vacuum:** No  
If the oxide layer is stained, grind and blast the surface again. Repeat the oxide firing. In the case of Zn-containing alloys, use a 10% sulfuric acid or similar pickling solutions if recommended. Note: Always observe proper safety procedures when handling acid. Before applying application clean framework carefully with water. Use the appropriate ceramic veneering material, following the manufacturer's instructions. High temperature firing temperature: 860 °C

**HEAT TREATMENT**  
**HARDENING:** 500 °C for 15 min; bench cool

**SOLDERING AND LASER WELDING**  
The soldering gap should not be wider than the thickness of the soldering material. Allow the soldered casting to cool slowly.

**Post Solder:** Porta Optimum Solder 940 Flux: High Fusing Bondal Flux  
**Pre Solder / Flux:** Porta Optimum Solder 710 Flux: Bondal Flux  
**Laser Welding Wire:** Portadur IN

**POLISHING**  
Carefully remove any oxide and flux residue. Smooth the metal surfaces with rubber polishes. Polish to a high gloss finish using polishing paste. Subsequently, clean using ultrasonic cleaning equipment or careful steam finishing.

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**SMALTING GOCH GUTTING**  
After being casted, carefully clean the casted with aluminum oxide (AlO<sub>3</sub>). Do not use a hammer for finishing. Finish the casting with carbide burrs and/or with ceramic-bonded grinding instruments. Blast the surface with 50 – 110 micron aluminum oxide (AlO<sub>3</sub>) at 2.0 bar. For veneering with Lab Composite material the framework must be finished and polished. Subsequently, steam clean or ultrasonic clean with distilled water to ethanol and dry the framework.

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