Gnathometer CAD

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Seite 3

Operating Instructions

Page 15

Mode d'emploi

Page 43

Istruzioni d'uso

Pagina 63

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Página 83

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Página 83

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


Table of Contents

1 Introduction / Signs and Symbols	17
1.1 Preface	17
1.2 Introduction	17
1.3 Notes regarding the Operating Instructions	17
2 Safety First	17
2.1 Indications	17
2.2 Health and safety instructions	17
3 Product Description	18
3.1 Construction and functional description	18
3.2 Indications, Contraindications	18
4 Installation and Initial Start-Up	18
4.1 Unpacking and checking the contents	18
4.2 Installation and initial start-up	18
5 Handling, Operation	19
5.1 Attaching the bite rim supports	19
5.2 Attaching the registration components	19
5.3 Adjusting the height of the registration pin	19
5.4 Mounting the registration plate with the fixing plates	19
6 Practical Use on the Patient	20
6.1 Preparing the registration	20
6.2 Performing the registration	20
7 Maintenance, Cleaning, Diagnosis	23
7.1 Monitoring and maintenance	23
7.2 Cleaning	23
7.3 Maintenance	23
8 What If ...	24
8.1 Technical malfunctions	24
8.2 Repair	24
9 Product Specifications	25
9.1 Delivery form	25
9.2 Technical data	25
10 Miscellaneous	25

Views of Gnathometer CAD / List of Parts

List of Parts

- 1 Mandibular registration plate
- 2 Maxillary registration plate
- 3 Bite rim supports
- 4 Registration pin holder
- 5 Retentions
- 6 Registration pin
- 7 Fixing plate
- 8 Snap fit
- 9 Fixation screw
- 10 Funnel hole
- 11 Stylus
- 12 Washer

1 Introduction/Signs and Symbols

1.1 Preface

Dear Customer,

Thank you for having purchased the "Gnathometer CAD". It is a highly-technical quality instrument. This product is distinguished for its good quality and high precision.

The Gnathometer CAD has been designed according to the latest industry standards. Inappropriate use may damage the equipment and be harmful to personnel. Please observe the relevant notes and read the Operating Instructions.

Enjoy working with Gnathometer CAD.

1.2 Introduction

The signs and symbols in these Operating Instructions facilitate the finding of important points and have the following meanings:

 Important information

 Risks and dangers

 Contraindication

1.3 Notes regarding the Operating Instructions

Device concerned: Gnathometer CAD
Target group: Dentists, dental professionals

These Operating Instructions facilitate the correct, safe and economical use of the Gnathometer CAD.

Should you lose the Operating Instructions, extra copies can be ordered from the respective Service Centre or downloaded as PDF from the internet.

2 Safety First

We recommend this chapter to all individuals working with the device or perform cleaning and maintenance work.

2.1 Indications

When manipulating the Gnathometer CAD while it is attached to the patient, make sure that the tongue of the patient is not injured by the stylus of the registration pin.

Risks and dangers:

2.1.1 Risk of swallowing and choking hazard:

Warning:

Since there are small components that may fall off if not attached correctly, there is a risk of swallowing these parts.

2.1.2 Risk of injury

Warning:

To prevent injury by the registration pin, the patient must be instructed to slide the tongue below the registration plate.

2.2 Health and safety instructions

Gnathometer CAD must only be used for the purpose described in Chapter 3. Further instructions to assure proper use of Gnathometer CAD:

- The instructions, regulations, and notes in these Operating Instructions must be observed.
- Gnathometer CAD must be properly maintained (see Chapter 7).

3 Product Description

3.1 Construction and functional description

Gnathometer CAD is an intraoral support pin registration device or the individualized determination of the centric jaw position in edentulous patients.

The basic principle of the intraoral support pin method is based on the idea of McGrane, Gysi, Gerber etc. of the three-point rest of the two jaw joints and the registration pin, which is placed in the balance point of the mandible. The support pin draws the typical image of the Gothic Arch on a coloured registration plate if the patient performs movements towards the right lateral and the left lateral out of the dorsal position. The tip of the Gothic Arch marks the retral contact position, while the centric occlusal position is marked by the crossing point. The desired centric position can then be fixed using the fixing plate. Subsequently, the entire registration is three-dimensionally locked with a suitable material (e.g. registration silicone or stone). After that, the registration is mounted on a special scan holder and placed into the scanner for digitization.

The special feature of this registration device is that the registration plates are equipped with a snap-action mechanism. With this snap-action mechanism, the bite rim supports and the registration plates can be quickly and easily attached and removed in order to conduct the different working steps. In a first step, the bite rim supports are used for the functional impressions. After that, the bite rim supports are removed and replaced with the registration plates. Then, needle point tracing is performed to determine the centric function.

3.2 Indications, Contraindications

Indications:

Determining the centric TMJ position in edentulous patients.

Contraindication:

The device cannot be used for dentulous or partially dentulous cases.

4 Installation and Initial Start-Up

4.1 Unpacking and checking the contents

Please make sure that the contents are complete.

- 1 Mandibular registration plate
- 1 Maxillary registration plate
- 2 Bite rim supports
- 1 Registration pin
- 1 Fixing plate
- 1 Fixation screw
- 1 Washer
- 1 Operating Instructions


4.2 Assembly and initial start-up

Important information:

Clean the components before they are used for the first time (see Cleaning in Chapter 7.2).

The Gnathometer CAD is already partly assembled. This general drawing depicts the various components and individual parts.


Please refer to the List of Parts on Page 16 for the part designations.

5 Handling, Operation

5.1 Attaching the bite rim supports


The bite rim supports (3) can be easily secured on and removed from the impression trays with the snap-action mechanism by means of the retentions (5).

5.2 Attaching the registration components

The registration components (1 & 2) can be easily secured on and removed from the impression trays with the snap-action mechanism.

Before needle point tracing is started, the registration pin should be secured with Pattern Resin or wax.

5.3 Adjusting the height of the registration pin


The height of the registration pin (6) can be adjusted by means of a thread. The screw thread provides more free space for increasing the vertical dimensions than for lowering the vertical dimensions. Therefore, it is recommended to perform the Centric Tray registration rather at a lower occlusal position and then lift it to the correct height with the registration.

Also, the fact that the functional impression results in the slight increase of the occlusal height must be taken into consideration.

The height of the registration pin (6) can be increased with a right turn.

The height of the registration pin (6) can be reduced with a left turn.

If necessary, the registration pin can be secured with Pattern Resin or wax.


5.4 Mounting the mandibular registration plate with fixing plates

The fixing plate (7) can be secured on the mandibular registration plate (1) by means of the fixation screw (9).

The the fixing plate (7) can be tightened with a right turn.

The the fixing plate (7) can be loosened with a left turn.


6 Practical Use on the Patient


6.1 Taking the functional impression

6.1.1 Placing the bite rim supports


Secure the bite rim supports (3) with the snap-action mechanism and the retentions (5) on the mandibular and maxillary impression trays. The bite rim supports (3) are identical in construction and fit both the maxillary and mandibular impression tray.

6.1.2 Applying the impression material into the impression tray and taking functional impressions


6.1.3 Removing the bite rim supports from the impression trays

6.2 Performing the registration

6.2.1 Attaching the registration components to the impression trays with the snap-action mechanism

Maxillary impression tray = secure the maxillary registration plate (2) with registration pin (6) in the registration pin holder (4).


Secure the mandibular impression tray = mandibular registration plate (1).


6.2.2 Verifying the vertical height

Adjust height, if necessary. Depending on the method, adequate auxiliaries (e.g. ruler, sliding caliper) or techniques (speech-sounds) can be used.

6.2.3 Coating the mandibular registration plate with colourant

Coat the mandibular registration plate (1) with a suitable colourant (crayon, grease pencil etc.).

6.2.4 Conducting the needle point tracing

Both functional impressions together with the registration components are replaced in the mouth of the patient. Now instruct the patient to execute the following movements:

It is recommended to practice these movements beforehand!


1. Move the mandible forward (protrusion) – and back, repeat if necessary
2. Move the mandible to the right (lateral movement) – and back, repeat if necessary
3. Move the mandible to the left (lateral movement) – and back, repeat if necessary
4. Move the mandible backward (retrusion) – and back, repeat if necessary

Important information:

The patient must sit upright and hold the head straight (no leaning).

The typical result of needle point tracing may look like the one shown in the picture. Please refer to the corresponding specialist literature for the interpretation of the various Gothic Arch tracings and their meaning.


Warning:

To prevent injury by the registration pin, the patient must be instructed to slide the tongue below the registration plate.

6.2.5 Aligning the fixing plate with the Gothic Arch tip and securing with the fixation screw

Now align the funnel hole (10) of the fixing plate (7) with the tip of the Gothic Arch and secure the fixing plate with the fixation screw (9).


6.2.6 Placing the registration in the mouth and verifying the centric position

Now, have the patient carefully close the mouth and verify that the stylus (11) of the registration pin (6) fits into the funnel hole (10) of the fixing plate (7) and the secured Gothic Arch registration. If this is not the case, the registration must be repeated. If the registration is correct, step 6.2.7 may be performed.

6.2.7 Locking the registration


The key is best fabricated with a registration silicone (e.g. CADBite), which is directly injected between the registration components from the syringe using the corresponding tip.

Moreover, esthetic lines (e.g. centre line, canine line, laugh line, lip closure line) can also be noted on the registration with a waterproof pen.

6.2.8 Removing the registration from the mouth

The entire locked registration must now be removed from the mouth of the patient.

6.2.9 Disinfecting the impression and registration components

Before passing on the components, they must be cleaned and disinfected using suitable disinfectants.

7 Maintenance, Cleaning, Diagnosis

This chapter describes the user maintenance and cleaning procedures for Gnathometer CAD. Only those tasks are listed that may be performed by dental professionals. All other tasks must be performed by qualified service personnel at a certified Ivoclar Vivadent Service Centre.

7.1 Monitoring and maintenance

The time for these maintenance procedures depends on the frequency of use and the working habits of the user. For that reason, the recommended times are only approximations.

7.2 Cleaning

Before disinfection, the Gnathometer must be completely disassembled.

Item	Frequency	Cleaning material
Registration components and small parts (metal)	After every use	Sterilization in the autoclave according to the valid standards or regulations
Fixing plates (plastic)	After every use	Wipe disinfection using a spray disinfectant according to the valid standards
Bite rim supports (plastic)	After every use	Wipe disinfection using a spray disinfectant according to the valid standards

7.3 Maintenance

Avoid any contact with strong acids and solvents (e.g. MMA, acetone) to prevent damage of the surfaces.

Sterilize and/or disinfect the components before the first use and after every use and/or after contact with saliva or blood.

- Use only disinfectants suitable for stainless steel and light metals.
- Avoid any contact with strong acids, brines or solvents.
- Clean in an ultrasonic bath, with steam jets or sterilization devices.

8 What If...

This chapter will help you to recognize malfunctions and take appropriate measures or, if possible, to perform some minor repairs.

The expendable parts are as follows:

- Fixing plates
- Bite rims

8.1 Technical malfunctions

Description	Instructions for users	Measures
The surface of the device is scratched and damaged	Was the surface cleaned or manipulated with a sharp object?	Do not scratch the surface with a sharp object
The surface of the device is damaged or partially dissolved	Was the surface cleaned with acid or a solvent?	Do not use any acids or solvents for cleaning
Fixing plate is dissolved and destroyed	Was the surface cleaned with an unsuitable disinfectant?	Do not use any unsuitable disinfectants. Order a new part (expendable part)
Bite rim support is dissolved and destroyed	Was the surface cleaned with an unsuitable disinfectant?	Do not use any unsuitable disinfectants. Order a new part (expendable part)

8.2 Repair

Repairs may only be carried out by a certified Ivoclar Vivadent Service Centre. If repairs during the warranty period are not carried out by a certified Service Centre, the warranty will expire immediately. Please also refer to the corresponding warranty regulations.

9 Product Specifications

9.1 Delivery form

Delivery forms may vary from country to country.

9.2 Technical data

Weight (approximately):

- Mandibular registration plate with registration pin: 13.1 g
- Maxillary registration plate with fixing plate, screw and washer: 12.8 g
- Bite rim support 1.6 g
- Stylus 1.0 g

Dimensions:

- Registration plate (into which the stylus is screwed)
Width: 29 mm / Length: 55.8 mm /
Material thickness: 2 mm
- Registration plate (to which the fixing plate is attached)
Width: 22 mm / Length: 55.8 mm /
Material thickness: 2 mm
- Bite rim support
Width: 32.5 mm / Length: 54.3 mm /
Material thickness: 2 mm

Material:

- Registration plates: Stainless steel
- Bite rim supports: POM white
- Fixing plates: PMMA, crystal clear

Adjustment options:

- Height adjustment of the stylus: 5 mm

10 Miscellaneous

Tips for Gothic Arch analysis and interpretation

Please also refer to the specialized literature available on this subject matter.

The device has been developed solely for use in dentistry. Start-up and operation should be carried out strictly according to the Operating Instructions. Liability cannot be accepted for damages resulting from misuse or failure to observe the Instructions. The user is solely responsible for testing the apparatus for its suitability for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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© Ivoclar Vivadent AG, Schaan/Liechtenstein
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.