Centric Tray

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For dental use only.

Description

Recording of the occlusal position

- Edentulous cases
- Partially edentulous cases

Indications

Determination of the occlusal position in edentulous patients

Contraindications

The use of Centric Tray is contraindicated if the patient is known to be allergic to any of its ingredients.

Side effects

There are no known side effects to date.


Composition

- Centric Tray: Polyamide
- Regular support, supporting wings: Polypropylene

Application

The Centric Tray and its components

- 1 Centric Tray
- 2 UTS Transferbow connection
- 3 Lip support
- 4 Retaining track
- 5 Regular support
- 6 Supporting wings

Clinical application

a) Edentulous cases

1. Take separate impressions, as normal.
2. Slip the support wings (6) distally onto the retaining track (4). Try-in the Ivoclar Vivadent Centric Tray, and adjust by trimming as necessary.
3. To record the occlusal height, first make marks on the patient's nose and chin, using a suitable marker.

Then proceed by selecting one of the options below:

- Replace the patient's dentures and measure the present vertical dimension.
 - Position the empty Centric Tray in the mouth. Ensure that the lips are moist and seated on the lip support (3) and that they gently support the tray. Leave the tray in-situ for 1–2 minutes, then measure and record the vertical dimension. Remove tray from mouth.
4. Mix silicone impression material, alginate or similar material and load both sides of the tray.
 5. Insert the loaded tray into the mouth sideways, rotating into position over the lower ridge. Place the moist lips on the lip support (3), whilst slowly closing the mouth until the marked occlusal height is reached. Allow the silicone, alginate or other material used to set completely. Check the vertical measure and remove the tray.
 6. The UTS Transferbow conveys patient-related information to the articulator. Align the Transferbow to the correct coordinates and attach it to the Ivoclar Vivadent Centric Tray at the UTS connection point (2). Fix the registration joint and remove the Transferbow.
 7. The occlusal record, with the UTS Transferbow or the registration joint attached, is now handed to the technician, who will mount it in the articulator based on the patient's coordinates.

b) Partially edentulous cases

1. Take separate impressions, as normal.
2. Check the closed occlusal position. If, due to the loss of contact between the opposing antagonists this is not possible, proceed as in Edentulous Cases (3).
3. Slip the regular supports (5) distally onto the retaining track (4), placing them into the gaps. Try-in the Ivoclar Vivadent Centric Tray and adjust by trimming as necessary.
4. Mix silicone impression material and load both sides of the tray.
5. Insert the loaded tray into the mouth sideways, rotating into position over the lower ridge. Place the moist lips on the lip support (3), whilst slowly closing the mouth until the marked occlusal height is reached. Allow the silicone, alginate or other material used to set completely. Check the vertical measure and remove the tray.
6. The UTS Transferbow conveys patient-related information to the articulator. Align the Transferbow to the correct coordinates and attach it to the Ivoclar Vivadent Centric Tray at the UTS connection point (2). Fix the registration joint and remove the Transferbow.
7. The occlusal record, with the UTS Transferbow or registration joint attached, is now handed to the technician, who will set-up the teeth according to the patient's coordinates.

Technical instructions

1. Trim back the vestibular edges of the silicone impression. Try to identify three reference points to the models.
2. Seat the models in the silicone impression and check the fit.
3. Mount the models in the articulator, positioning them in relation to the occlusal record.
4. If a UTS Transferbow has been used, attach it to the articulator, applying the philosophy relative to that determined by the Transferbow record.
5. Secure the Ivoclar Vivadent Centric Tray in place. Place and fix the maxillary model. Remove lower jig and locate and fix the mandibular model.
6. Fabricate individual trays, centre pin recording devices, wax occlusal rims or other prescription requests.

Please observe the following:

- Keep a record of the measured occlusal height.
- Use a viscous silicone, alginate or similar material.
- Use sufficient impression material.
- Check the seating of the models.
- Consider the articulator philosophy.

Reprocessing

- Autoclavable tray for multiple use (recommended for 10 autoclave cycles. If the tray is subjected to further autoclave cycles, the inherent stability of the tray may be compromised and thus the exact dimensions are no longer ensured).
- Manual disinfection:
Disinfection procedures:
Disassemble the tray before disinfecting it. Place the Centric Tray into a mesh basket. Immerse the mesh basket with the individual parts of the tray in disinfectant solution, e.g. MD 520 (Dürr Dental), at 20 ± 2 °C/ 68 ± 3.6 °F for 5 min. Make sure that the surfaces of the components are completely covered with disinfectant solution and that they are free of bubbles during the reaction time. Rinse Centric Tray with cold demineralized water to remove any disinfectant residues.
- Sterilization: For steam sterilization, use the following parameters: 3x pre-vacuum procedure, 134 °C/273.2 °F, sterilization time 5.5 min.

Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclarvivadent.com, and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclarvivadent.com).

Disposal information

- Remaining stocks must be disposed of according to the corresponding national legal requirements.

Shelf life and storage

This product does not require any special storage conditions.

Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions.

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.