
IPS e.max[®] Abutment Solutions

Cem Kit

EN Instructions for Use

Page 2

Set for the clinical try-in and permanent extraoral cementation of hybrid abutment solutions

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

Hybrid abutment crown

Description

The IPS e.max® Abutment Solutions Cem Kit is a professional kit containing quality materials for the permanent cementation of ceramic structures made of lithium disilicate glass-ceramic (LS₂) or zirconium oxide on bases (e.g. abutment or adhesive basis) made of titanium/titanium alloy or zirconium oxide in the fabrication of hybrid abutments or hybrid abutment crowns. Furthermore, the kit contains materials for the clinical try-in procedure.

Indication

Extraoral, permanent cementation of ceramic structures made of lithium disilicate glass-ceramic (LS₂) or zirconium oxide on titanium/titanium alloy or zirconium oxide bases. Provisional securing of ceramic structures on the base as well as securing of the crown during try-in.

Contraindication

The materials are contraindicated

- for intraoral luting;
- if the stipulated working technique cannot be applied;
- if a patient is known to be allergic to any of the ingredients of Multilink Hybrid Abutment, Monobond Plus and Virtual Extra Light Body Fast Set.

Composition

The IPS e.max Abutment Solution Cem Kit includes the following components:

- **Multilink® Hybrid Abutment** – self-curing luting composite for the permanent cementation of ceramic structures made of lithium disilicate glass-ceramic (LS₂) or zirconium oxide on bases (e.g. abutment or adhesive basis) made of titanium/titanium alloy or zirconium oxide. The monomer matrix is composed of dimethacrylate and HEMA. The inorganic fillers include barium glass, ytterbium trifluoride, spheroidal mixed oxide and titanium dioxide. The particle size is between 0.15 and 3.0 µm. The mean particle size is 0.9 µm. The total volume of inorganic fillers is approx. 36%.

Shade

- H00 (High Opacity)

The working and curing times depend on the ambient temperature. Once Multilink Hybrid Abutment has been dispensed from the automix syringe, the following times apply:

	At room temperature 23 °C ± 1 °C/73 °F ± 2 °F
Working time	approx. 2 minutes
Curing time (including working time)	approx. 7 minutes

The working and curing times become shorter at higher temperatures, while they become longer at lower temperatures.

- **Monobond® Plus** – universal primer promoting an adhesive bond between luting composites and all indirect restorative materials

Alcohol solution of silane methacrylate, phosphoric acid methacrylate and sulphide methacrylate

- **Liquid Strip** – oxygen-impermeable glycerine gel for the prevention of an oxygen-inhibited layer on composite surfaces

Glycerine gel

- **Virtual® Extra Light Body Fast Set** – vinylpolysiloxane based impression material (addition-curing) for provisionally securing components in place during try-in

Addition-curing silicone containing vinylpolysiloxane, methylhydrogensiloxane, organoplatinic complex, silica and food dyes

Side effects

- **Multilink Hybrid Abutment, Monobond Plus**
Systemic side effects are not known. In rare cases, allergic reactions to individual components have been reported.

Interactions

- **Multilink Hybrid Abutment:**
Phenolic substances (e.g. wintergreen oil, eugenol) inhibit polymerization. Therefore, materials containing these components must not be used. Oxidative disinfectants (e.g. hydrogen peroxide) may interact with the initiator system and impair the setting process. As a result, the automix syringe must not be disinfected with oxidative agents. The syringe can be disinfected by wiping it with medical alcohol for example.
- **Monobond Plus:**
IMPORTANT! In order to create a strong bond, zirconium oxide and metal surfaces must not be cleaned with phosphoric acid.
Do not insert the brush into the Monobond Plus bottle. This may lead to contamination.
- **Virtual Extra Light Body Fast Set:**
The setting of vinylpolysiloxanes is inhibited by latex gloves. The surfaces to which the silicone will be applied (ceramic structure, titanium base, etc.) must not be touched with latex gloves.

e.max Abutment Solutions

Optional: Clinical Try-in

Provisional securing of the ceramic structure on the base

Before the ceramic structure is permanently luted to the base, a clinical try-in can take place. To facilitate intraoral handling, the components are temporarily attached to one other with silicone material, e.g. Virtual Extra Light Body Fast Set.

The following procedure should be observed in the temporary attachment of the components:

- The untreated base and the ceramic structure are cleaned with steam and subsequently dried with blown air.
- The ceramic structure is placed on the base (which is screwed on the model analog) and the relative position of the components is marked with a waterproof pen. This step makes it easier to attain the correct position when the parts are temporarily assembled.
- The screw channel is sealed with a foam pellet.
- The Virtual cartridge is inserted in the dispenser and the protective cap is removed.
- The mixing tip is screwed on and the Oral Tip is attached to the mixing tip.
- Virtual Extra Light Body Fast Set is applied to the base and directly into the ceramic structure.
- The base is introduced into the ceramic structure. The alignment of the two components must be checked (rotation lock/markings).
- The components should be held firmly in the correct position for 2:30 minutes until Virtual Extra Light Body Fast Set has set.
- Any excess that has been displaced must be carefully removed with a suitable instrument, e.g. a scalpel.

Cleaned, untreated ceramic structures

The ceramic structure is placed on the base and the relative position is marked.

The screw channel is sealed with a foam pellet.


The Virtual cartridge is inserted in the dispenser. The mixing tip is screwed on and the Oral Tip is attached.


Virtual Extra Light Body Fast Set is applied to the base and...

...directly on the ceramic structure.

The base is introduced into the ceramic structure. In doing so, the alignment of the two components is checked (rotation lock/markings). The components are firmly held in place for approx. 2:30 minutes until Virtual Extra Light Body Fast Set has set.

Displaced excess is carefully removed with a suitable instrument, e.g. a scalpel.

Excess Virtual Extra Light Body Fast Set is removed from the screw channel with an instrument.

Clinical try-in

Hybrid abutment with dedicated crown

Important note: Any intraoral inspection of the occlusion/articulation and any necessary grinding adjustments may only be carried out if the components have been attached to one other with Virtual Extra Light Body Fast Set. Virtual has a cushioning effect during the try-in procedure, in particular, if any grinding adjustments have to be made. Therefore, it prevents chipping in the transition area between the hybrid abutment and the crown.

The following procedure should be observed during the clinical try-in:


- The prepared hybrid abutment (provisionally secured in place) and the clean corresponding crown are laid out.
- The provisional restoration is removed.
- The hybrid abutment is screwed in manually with the dedicated screw.
- The geometry of the hybrid abutment is checked (e.g. fit, gingival anaemia) in relation to the gingival margin.
- If desired, the screw channel on the hybrid abutment can be sealed with a foam pellet.
- **Tip:** Isolate the inner aspect of the crown, e.g. with Try-In-Paste, Liquid Strip
- **The crown is placed on the hybrid abutment intraorally to check and if necessary adjust the proximal contacts. Note: No occlusal functional inspection must be performed at this stage.**
- **For the functional inspection, the crown has to be secured on the hybrid abutment with Virtual Extra Light Body Fast Set. Try-in paste must not be used for this purpose, as this material is not sufficiently resistant to compressive force.**
- The Virtual cartridge is inserted in the dispenser and the protective cap is removed.
- The mixing tip is screwed on and the Oral Tip is attached to the mixing tip.
- Virtual Extra Light Body Fast Set is applied to the inner aspect of the crown.
- The crown is pressed onto the hybrid abutments using the fingers until the final position is reached. The crown is held in the final position until the Virtual material has set.
- Excess Virtual material is removed.
- The occlusion/articulation is checked and if necessary adjustments are made with suitable grinding instruments (see separate IPS e.max recommended grinding instruments for ceramics – use in the dental practice). If adjustments have been made by grinding, the restoration has to be re-polished or glaze fired again.
- The crown is carefully removed from the hybrid abutment, and the hybrid abutment (including the base).
- The temporary restoration is placed.

The hybrid abutment is manually screwed in place with the dedicated screw. The geometry of the hybrid abutment (e.g. fit, gingival anaemia) is checked in relation to the gingival margin.


If desired, the screw channel of the hybrid abutment can be sealed with a foam pellet.

Tip: The inner aspect of the crown can be isolated with glycerine gel.

The crown is placed on the hybrid abutment intraorally to check and if necessary adjust the proximal contacts. **Note: No occlusal functional inspection must be performed at this stage.**


Virtual Extra Light Body Fast Set is applied to the inner aspects of the crown.

The crown is pressed onto the hybrid abutment using the fingers until the final position is reached. The crown is held in the final position until the Virtual material has set.

Excess Virtual material is removed.

The occlusion/articulation is checked and if necessary adjustments are made with suitable grinding instruments.

The crown is carefully lifted from the hybrid abutment and the Virtual Extra Light Body Fast Set is removed.

The hybrid abutment is unscrewed.

Hybrid abutment crown

The following procedure should be observed during the clinical try-in:

- The prepared and cleaned hybrid abutment crown (provisionally secured in place with Virtual Extra Light Body Fast Set) is laid out.
- The provisional restoration is removed.
- The hybrid abutment crown is placed on the implant intraorally in order to check and if necessary adjust the proximal contacts. **Note: No occlusal functional inspection must be performed at this stage.**
- The hybrid abutment crown is screwed in manually with the dedicated screw.
- The geometry of the hybrid abutment crown is checked (e.g. fit, gingival anaemia) in relation to the gingiva.
- The occlusion/articulation is checked and if necessary adjustments are made with suitable grinding instruments (see separate IPS e.max recommended grinding instruments for ceramics – use in the dental practice). If adjustments have been made by grinding, the restoration has to be re-polished or glaze fired again.
- The crown is carefully removed from the hybrid abutment crown (including the base).
- The implant site is rinsed, e.g. with Cervitec® Liquid (antibacterial mouth rinse containing chlorhexidine) to clean and disinfect it.
- The temporary restoration is placed.

The hybrid abutment crown is placed on the implant intraorally in order to check and if necessary adjust the proximal contacts. **Note: No occlusal functional inspection must be performed at this stage.**

The hybrid abutment crown is screwed in with the dedicated screw.

The geometry of the hybrid abutment crown (e.g. fit, gingival anaemia) is checked in relation to the gingiva.

The occlusion/articulation is checked and if necessary adjusted with suitable grinding instruments.

The hybrid abutment crown (including base) is carefully removed.

IPS e.max® Abutment Solutions

Permanent cementation of base / ceramic structure

Careful preparation of the bonding surface is a requisite for the successful adhesive cementation of the base and the ceramic structure. The following paragraphs outline the required procedures. The procedure is the same for hybrid abutments and hybrid abutment crowns.

Hybrid abutments and hybrid abutment crowns made of IPS e.max lithium disilicate glass-ceramic (LS₂)

	Ceramic structure (LS ₂)	Base
Blasting	–	According to the instructions of the manufacturer
Etching	The bonding area with IPS® Ceramic Etching Gel for 20 s	–
Conditioning	The bonding area with Monobond Plus for 60 s	
Adhesive cementation	Multilink Hybrid Abutment	
Covering the cementation joint	Glycerine gel, e.g. Liquid Strip	
Curing	Approx. 7 min auto-polymerization	
Polishing the cementation joint	Conventional polishers for ceramic/composite resin	

Hybrid abutments made of zirconium oxide (ZrO₂)

	Ceramic structure (ZrO ₂)	Base
Blasting	According to the instructions of the manufacturer	According to the instructions of the manufacturer
Conditioning	The bonding area with Monobond Plus for 60 s	
Adhesive cementation	Multilink Hybrid Abutment	
Covering the cementation joint	Glycerine gel, e.g. Liquid Strip	
Curing	Approx. 7 min auto-polymerization	
Polishing the cementation joint	Conventional polishers for ceramic/composite resin	

Preparation of the base

The following procedure should be observed in the preparation of the base for the cementation with the ceramic structure:

- The base should be prepared according to the instructions of the manufacturer.
- The base is cleaned in an ultrasonic bath or with a steam cleaner and then dried with blown air.
- The base is screwed on the model analog.
- The ceramic structure is placed on the base and the relative position of the components is marked with a water-proof pen. This facilitates locating the correct position when the parts are assembled at a later stage.

- The emergence profile of the base must not be blasted or modified in any way!
- **If the manufacturer recommends that the bonding surfaces be blasted, the following procedure should be observed:**
 - Hard modelling wax is applied to protect the emergence profile, since this type of material is easy to remove later on.
 - The screw channel is also sealed with wax.
 - The bonding area is carefully blasted according to the instructions of the manufacturer.
 - An instrument and a steam cleaner are used for cleaning. It is important to remove all the wax meticulously.
- After the bonding surface has been cleaned, it must not be contaminated under any circumstances, as this would impair the bond.
- Monobond Plus is applied on the clean bonding surface and allowed to react for 60 s. After the reaction time, any residue is dispersed with air that is free of water and oil.
- The screw channel is sealed with a foam pellet or wax. The bonding surface must not be contaminated in the process.


If the manufacturer recommends that the bonding surfaces be blasted, the following procedure should be observed:

The base is screwed on the model analog. The relative position to the ceramic structure is marked with a waterproof pen.

Optional: Wax is applied in order to protect the emergence profile. In addition, the screw channel is sealed with wax.

Optional: The bonding surface can be carefully blasted according to the instructions of the manufacturer.

Optional: Cleaning with an instrument and a steam cleaner. It is important to remove all the wax meticulously.

Monobond Plus is applied to the clean bonding surface and allowed to react for 60 s. After the reaction time, any remaining residue is dried with blown air that is free of water and oil.


The screw channel is sealed with a foam pellet or wax. The bonding surface must not be contaminated in the process.

Preparing ceramic structures made of IPS e.max lithium disilicate glass-ceramic (LS₂)

The following procedure must be observed in the preparation of the ceramic structure for cementation on the base:


- The ceramic structure **must not** be blasted in preparation for cementation.
- The ceramic structure is cleaned in an ultrasonic bath or with a steam cleaner and subsequently blown dry.
- After cleaning, the bonding must not be contaminated under any circumstances, as this would impair the bond.
- Wax can be applied to protect the outer surfaces or the glazed areas.
- The bonding surface is etched with 5% hydrofluoric acid gel (IPS Ceramic Etching Gel) for 20 s.
- Subsequently, the bonding surface is thoroughly rinsed under running water and dried with air that is free of water and oil.
- Monobond Plus is applied to the clean bonding surface and allowed to react for 60 s. After this reaction time, any residue is dispersed with air that is free of water and oil.

Ceramic structures made of IPS e.max lithium disilicate glass-ceramic (LS₂) **must not** be blasted.

Etching with IPS Ceramic Etching Gel for 20 s and subsequently rinsing it off with water.

Monobond Plus is allowed to react for 60 s and then it is dried with blown air.

Preparing ceramic structures made of zirconium oxide (ZrO₂)

The following procedure must be observed in the preparation of the ceramic structure for cementation on the base:

- In order to facilitate its handling, the ceramic structure is attached to a holder (e.g. Microbrush) before it is blasted.
- The outer surfaces of the ceramic structure are protected with a wax coating.
- In order to better control the blasting procedure, the bonding surface can be coloured with a permanent marker.
- The bonding surface of the ceramic structure is cleaned by means of blasting with Al₂O₃ 100 µm at max. 1 bar and at a distance of max. 1–2 cm from the blasting nozzle (**The instructions of the manufacturer must be observed**).
- The ceramic structure is cleaned with steam or in an ultrasonic bath.
- After cleaning, the bonding surface must not be contaminated under any circumstances, since this may negatively influence the bonding results.
- Monobond Plus is applied to the clean bonding surface and allowed to react for 60 s. After this reaction time, any remaining residue is dried with air that is free of water and oil.

Optional: Wax is applied to protect the outer surfaces.

Optional: The bonding surface is cleaned by blasting it (**The instructions of the manufacturer must be observed**).

Monobond Plus is allowed to react for 60 s and then dried with blown air.

Cementation with Multilink Hybrid Abutment

The following instructions must be observed in the cementation procedure:

- The cleaned and conditioned components (ceramic structure, base) are laid out ready for cementation.
- **The subsequent cementation procedure must be carried out quickly and without interruption.**
- **The working time of Multilink Hybrid Abutment is approx. 2 min at 23 °C (± 1 °C)/73 °F (± 2 °F).**
- As a general rule, a new mixing tip is attached to the Multilink Hybrid Abutment prior to each use.
- A thin layer of Multilink Hybrid Abutment is directly applied from the mixing syringe to the bonding surface of the base and **the bonding surface of the ceramic structure.**
- The mixing tip is left on the Multilink Hybrid Abutment syringe until the next use. The remaining cement polymerizes in the tip and functions as a seal.
- The ceramic structure is placed on the base in such a way that the position markings are aligned.
- The parts are lightly and evenly pressed together and the correct relative position of the components is checked (transition base/ceramic structure).
- Subsequently, the parts are tightly pressed together for 5 s.
- Excess in the screw channel is carefully removed, e.g. with a Microbrush or brush, using rotary movements.
- **Important: Excess must not be removed before curing has started, i.e. 2–3 minutes after mixing. For the purpose, a suitable dental lab instrument (e.g. Le Cron) is used and the components are held in place with light pressure in the process.**
- Glycerine gel is applied (e.g. Liquid Strip) on the cementation joint to prevent the formation of an inhibition layer. The glycerine gel must be applied cautiously to avoid blending it with or displacing the composite. The gel must be left on the cementation joint until polymerization is complete.
- Next, the luting composite is completely auto-polymerized within 7 min.
- **Important: The components must not be moved until Multilink Hybrid Abutment has completely cured. They can be held immobile with, e.g. diamond-coated tweezers.**
- After the completion of auto-polymerization, the glycerine gel is rinsed off with water.
- **The cementation joint should be cautiously polished with rubber polishers at a low speed (< 5,000 rpm) to avoid overheating.**
- Any cement residue left in the screw channel is removed with suitable rotating instruments.
- The restoration is steam cleaned.

The cleaned and conditioned components are laid out ready for cementation.

A new mixing tip is attached to the Multilink Hybrid Abutment prior to each use.


A thin layer of Multilink Hybrid Abutment is directly applied from the mixing syringe to the bonding surface of the base.

A thin layer of Multilink Hybrid Abutment is directly applied from the mixing syringe on the bonding surface of the ceramic structure.

The ceramic structure is placed on the base in such a way that the position markings are aligned.

The components are joined using even and light pressure. Subsequently, the components are tightly pressed together for 5 s.

Excess in the screw channel is carefully removed, e.g. with a Microbrush or brush, using rotary movements.

Important: Excess must not be removed before curing has started, i.e. 2-3 minutes after mixing. For the purpose, a suitable dental lab instrument (e.g. Le Cron) is used and the components are held in place with light pressure in the process.

Glycerine gel (e.g. Liquid Strip) is applied on the cementation joint to prevent the formation of an inhibition layer.


The luting composite auto-polymerizes within 7 min.
Important: The components must not be moved until auto-polymerization is completed. The components must be immobilized during this time.

After the auto-polymerization process, the glycerine gel is rinsed off with water.

The cementation joint is cautiously polished with rubber polishers at a low speed (< 5,000 rpm) to avoid overheating.

Any remaining cement residue in the screw channel is removed with suitable rotating instruments. The base must not be damaged.

Completed hybrid abutment and hybrid abutment crown after cementation.

Additional notes

- Multilink Hybrid Abutment should be processed at room temperature. Low(er) temperatures may cause difficulty in dispensing and mixing the material and may prolong the working and curing time.
- Unpolymerized material may stain clothing.

Warning

- Contact of uncured Multilink Hybrid Abutment with skin/mucous membrane and eyes must be prevented. In an unpolymerized state, Multilink Hybrid Abutment may cause slight irritation and lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.
- Monobond Plus:
 - Highly flammable
 - Irritating
 - The product contains ethanol.
 - The product must be kept away from sources of ignition.
- Virtual Extra Light Body Fast Set:
If uncured material comes in contact with the eyes, they should immediately be rinsed with plenty of water. If irritation persists, a physician should be consulted. If the material comes in contact with the skin, it should be rinsed off with water and soap.

Shelf life and storage


- The individual components of the IPS e.max Abutment Solution Cem Kit should not be used after the indicated date of expiration.
- Shelf life: see note on syringes/bottles/cartridges
- Storage temperature of Multilink Hybrid Abutment, Monobond Plus, Liquid Strip and Virtual Extra Light Body Fast Set: 2–28 °C/36–82 °F.
- Monobond Plus and Liquid Strip: Close immediately after use.
- Multilink Hybrid Abutment and Virtual Extra Light Body Fast Set:
The mixing tip is left on the syringe/cartridge after use. The material contained within will polymerize and allow the mixing tip to function as a seal.

Keep out of the reach of children!


For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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© Ivoclar Vivadent AG, Schaan, Liechtenstein
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.