AcrylCAD

Instructions for Use Gebrauchsinformation Mode d'emploi Istruzioni d'uso Instruções de uso

615212/0712/REV1

English

Instructions for Use

Product description

IPS AcrylCAD Blocks are burn-out acrylic polymer blocks. They can be used instead of modelling wax for the automated fabrication of elements which are subsequently invested and cast, pressed or pressed over a substrate material.

Main indication

- Milling of cast restorations which can be invested for
 - the fabrication of primary crowns using the telescope technique
 - the fabrication of crown and bridge frameworks (metal-ceramics)
 - the fabrication of crowns and bridges (full-cast technique)
- Milling of investable restorations for the press technique for
 - the fabrication of crown and bridge frameworks, e.g. made of IPS e.max® Press
 - the fabrication of full-contour restorations, e.g. made of IPS e.max Press, IPS Empress® Esthetic
- Press-on technique: Milling of investable veneering elements for the press-on technique
 - the fabrication of metal-supported crowns and bridges, made e.g. of IPS InLine® PoM
 - the fabrication of crowns and bridges supported by a zirconium oxide framework, made e.g. of IPS ZirPress

Secondary indication

Try-in of milled elements before further processing

Limitations of use

- IPS AcrylCAD must not be used for the fabrication of temporary restorations or cemented (whether temporarily or permanently) in the patient's oral cavity.
- Any use not listed in the indications.

Composition

Polymethyl methacrylate (PMMA)

Note

- Do not use in case of known hypersensitivity to methacrylates.
- Failure to observe the stipulated limitations of use and processing instructions may lead to failure.

Processing with CEREC® MC-L / MC-XL and inLab® MC-L / MC-XL (Sirona)

In comparison with ceramic materials, a reduced quantity of Dentatec fluid per tank filling is required. Please refer also to the Ivoclar Vivadent abrasives overview.

The CAD-Waxx Starter Kit from Sirona with a modified tank and reinforced filter system (specialized dental dealers) is required to process IPS AcrylCAD in MC-L machines.

1. Scan the plaster model. When fabricating elements for the press-on technique, the double-scan technique or the "Quadrant" function (virtual placement) can be used.
2. Select IPS AcrylCAD from the material selection and construct the element using the software. Observe the minimum thickness for the intended application. The minimum thicknesses programmed in the software refer to alloy frameworks for the casting technique. If frameworks for the press technique, resin over-structures for the press-on technique or full-contour (pressed) restorations are designed, the minimum layer thickness must be at least 0.8 mm in order to ensure complete pressing. In the case of pressed bridges, the connector size must be at least 16 mm². Please observe the specifications of the respective manufacturer.
3. Clamp in the IPS AcrylCAD block and mill with the appropriate grinding tools.
4. After milling, cut the element from the attachment points. Use cross-cut tungsten carbide burs to remove the attachment points or to make shape adjustments.
5. If desired, the marginal fit and the shape can be adapted and optimized using burn-out modelling wax or resin. Due to potential deformations resulting from increased temperatures (>50°C) it is not recommended to use light-curing modelling waxes.
6. The further processing, e.g. casting with metal or pressing with ceramic material, is carried out according to the instructions for use of the respective manufacturer.
7. The Ivoclar Vivadent investment materials IPS® PressVEST (conventional) and IPS PressVEST Speed may be used to press invested IPS AcrylCAD elements with IPS e.max Press, IPS e.max ZirPress and IPS InLine PoM.

Date of information prepared:

07/2012, REV1

Manufacturer:

Ivoclar Vivadent AG
FL-9494 Schaan / Liechtenstein

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damage resulting from other uses. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.