# UTS 3D

## Universal Transferbow System

### **Operating Instructions**

Page 3

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Seite 19

### **Mode d'emploi**

Page 35

### **Istruzioni d'uso**

Pagina 51

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Página 67

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Pagina 83

# UTS 3D

	Page
<b>List of Parts</b>	<b>4</b>
<b>1. Introduction / Signs and Symbols</b>	<b>7</b>
1.1 Preface	
1.2 Introduction	
1.3. Notes regarding the Instructions for Use	
<b>2. Safety First</b>	<b>7</b>
2.1 Stipulated use	
2.2 Health and safety instructions	
<b>3. Product Description</b>	<b>8</b>
3.1 Design and functional description	
3.2 Indication, contraindication	
<b>4. Assembly and Initial Start-Up</b>	<b>9</b>
4.1 Unpacking and checking the contents	
4.2 Assembly and initial start-up	
<b>5. Handling</b>	<b>10</b>
5.1 Earpieces	
5.2 Nosepiece	
5.3 Reference plane indicator	
5.4 3D-registration joint	
5.5 Adjusting the width	
5.6 Additional height of the nosepiece	
<b>6. Practical Use on the Patient</b>	<b>12</b>
6.1 Preparing the bite fork	
6.2 Reference plane	
6.3 Positioning the transferbow	
6.4 Positioning the nosepiece	
6.5 Aligning the transferbow	
6.6 Other versions	
6.7 Securing the 3D-registration joint	
6.8 Removing the transferbow	
6.9 Removing the 3D-registration joint	
<b>7. Maintenance, Cleaning, Diagnosis</b>	<b>15</b>
7.1 Inspection and maintenance	
7.2 Cleaning	
7.3 Cleaning instructions for the 3D-registration joint	
<b>8. What If .....</b>	<b>16</b>
8.1 Technical malfunctions	
8.2 Repair	
<b>9. Product Specifications</b>	<b>17</b>
9.1 Delivery form	
9.2 Technical data	
<b>10. Miscellaneous</b>	<b>17</b>

## List of Parts

### **1.0 Universal transferbow (complete)**

- 1.1. Earpiece mounting
- 1.2. Earpiece rod screw
- 1.3. Side arms
- 1.4. Side arm support
- 1.5. Scale
- 1.6. Multifunctional mounting
- 1.7. Width setting screw
- 1.8. Registration joint connection
- 1.9. Connection for nosepiece
- 1.10. Nosepiece holder screw
- 1.11. Registration joint connection screw
- 1.12. Thread for screw 3.1

### **2.0 Individually adjustable nosepiece (complete)**

- 2.1. Nosepiece holder
- 2.2. Nosepiece rod
- 2.3. Nosepad holder
- 2.4. Nosepiece rod screw

### **3.0 Reference indicator**

- 3.1. Reference indicator screw

### **4.0 Axis pins**

- 4.1. Thread
- 4.2. Adapter pin for articulator

### **5.0 Intercondylar distance (S,M,L)**

### **6.0 Earpieces (complete)**

- 6.1. Earpiece rod
- 6.2. Threaded openings for axis pins
- 6.3. Earpiece
- 6.4. Threaded pin for earpieces
- 6.5. Earpiece support

### **7.0 3D-Registration joint**

- 7.1. Connector
- 7.2. Tommy screw
- 7.3. Bite fork connector
- 7.4. Bite fork connector screw
- 7.5. 3D-bite fork connector

### **8.0 3D-Bite fork**

- 8.1. Connector

### **9.0 Earpiece, black, 'large'**



- 9.1. Threaded connector

### **10.0 Nosepads**

- 10.1. Nosepad, green (medium)
- 10.2. Nosepad, blue (soft)

### **11.0 3D-Bite fork "S" (small)**

- 11.1. 3D-bite fork "S" connector


3D Registration joint (7.0)


3D Bite fork (8.0)

3D Bite fork "S" (11.0)

Individually adjustable nosepiece (2.0)

soft (10.2) medium (10.1)

Nosepads, 'M' and 'S'

UTS, complete (1.0)

Earpieces, complete (6.0)

Reference indicator (3.0)

Earpiece, 'large' (9.0)

# 1. Introduction / Signs and Symbols

## 1.1 Preface

Dear Customer,  
Thank you for having purchased the UTS 3D transferbow. This instrument is a highly technical product that is distinguished for its excellent quality and high precision.

This transferbow has been designed according to the latest industry standards. Inappropriate use may result in certain risks. Please observe the corresponding notes and become familiar with the Instructions for Use.

We wish you much success with the UTS 3D transferbow.

## 1.2 Introduction

The signs and symbols in these Instructions for Use facilitate the finding of important points and have the following meanings:

Risks and dangers

Important information

Contraindication

## 1.3 Notes regarding the Instructions for Use

Device concerned: UTS 3D System  
Target group: Dentists, dental technologists

These Instructions for Use facilitate the correct, safe, and economic use of the UTS 3D.

Should you lose these Instructions for Use, extra copies can be ordered at a nominal fee from your local Ivoclar Vivadent Service Center.

# 2. Safety First

## 2.1 Stipulated use

This chapter is especially important for personnel who work with the UTS 3D or who have to carry out maintenance or repair work.


Risks and dangers:

When manipulating the transferbow on the patient, make sure that the reference indicator cannot injure the patient's eyes.

Risk and dangers:

Do not apply too much pressure when positioning the ear-pieces. This may result in injury to the auditory meatus.

## 2.2 Health and safety instructions


The UTS 3D must exclusively be used for the indications listed in Chapter 3. Further instructions to ensure proper use of the transferbow:

- The instructions, regulations, and notes in these Instructions for Use must be observed.
- The UTS 3D must be properly maintained (see Chapter 7).

### 3. Product Description


#### 3.1 Design and functional description

The Universal Transferbow System is a transfer system developed on a scientific basis. It is a fact that there are larger and smaller jaws. The Bonwill triangle of various patients may differ up to 50 mm, depending on the size of the skull.


Individual, spatial model orientation by means of a transferbow takes these parameters into account and helps avoid occlusion errors.

The UTS 3D permits the spatial transfer of a patient's individual Bonwill triangle according to the joint elements of the articulator.

With its adjustable nosepiece, the UTS 3D is also suitable for model orientation according to Camper's plane (CP) and the Frankfort horizontal (FH). The adjustable ear-pieces also make the UTS 3D suitable for model transfer according to specific coordinates.

The Stratos 200 permits working according to these reference planes (CP and FH). For this purpose, two different registration joint holders are available (CP and FH). The registration joint holders are used according to the reference plane selected. Thanks to the registration joint holder of the Stratos 200, the transferbow no longer has to be transported to the dental laboratory. Only the 3D-registration joint and the bite fork have to be sent to the lab.

e.g. CP version

The UTS 3D may also be directly connected to the articulator by means of the vertically adjustable type 2 support pins. For this purpose, however, the entire transferbow including the 3D-registration joint have to be sent to the dental laboratory.

#### 3.2 Indication, contraindication

##### Indication:

- *CP version*  
Arbitrary skull/joint-related model orientation according to Camper's plane
- *FH version*  
Arbitrary skull-joint related model orientation according to the Frankfort horizontal
- *Other*  
Model orientation according to specific coordinates

##### Contraindication:

No contraindications are known to date, provided the UTS 3D is used strictly according to the Instructions for Use.

## 4. Assembly and Initial Start-Up

### 4.1 Checking the contents


Check the delivery for completeness:

- 1 Universal transferbow, complete
- 1 Nosepiece, complete
- 1 Reference indicator
- 1 3D-Registration joint
- 1 3D-Bite fork
- 2 Earpieces with support, complete
- 5 Nosepads, green (medium)
- 5 Nosepads, blue (soft)
- 1 Instructions for Use

### 4.2 Assembly and initial start-up


#### Mount earpieces

Mount the desired earpieces (6.3) (M or L) on the threaded pins for earpieces (6.4). Loosen the earpiece rod screw (1.2) and mount the entire assembly on the ear-piece mounting (1.1). Make sure the earpieces are correctly positioned.


#### Attach nosepiece with nosepads

Place the desired nosepad (10.1 or 10.2) in the nosepad holder (2.3). Insert the nosepiece rod (2.2) into the nosepiece holder (2.1).


The nosepiece holder (2.1) may now be attached to the connection for the nosepiece (1.9) using the nosepiece holder screw (1.10) of the multifunctional mounting (1.6).


#### Reference indicator

Mount the reference indicator (3.0) on the thread for screw 3.1 (1.12) by means of the reference indicator screw (3.1).

#### 3D-registration joint with bite fork

Insert the connector (8.1) of the bite fork (8.0) into the bite fork connector (7.5) and tighten the bite fork connector screw (7.4). Insert the connector (7.1) into the registration joint connection (1.8) and secure the 3D-registration joint (7.0) with the registration joint connection screw (1.11).

#### Axis pins

The axis pins (4.0) are mounted in the laboratory if the UTS 3D is directly connected to the articulator.


When using the UTS 3D on the patient, the axis pins (4.0) must not be used in order to prevent injury.


## 5. Handling

### 5.1 Earpieces


The earpiece rod (6.1) may be adjusted by loosening the earpiece rod screw (1.2).


**i** Please note that the earpiece support (6.5) is touching the earpiece mounting (1.1) if the CP or FH versions are used.

### 5.2 Nosepiece


The height of the nosepiece holder may be adjusted to the patient's individual reference point by loosening the nosepiece holder screw (1.10).


By loosening the nosepiece rod screw (2.4) the nosepiece rod (2.2) may be adjusted horizontally to fit the patient.

### 5.3 Reference indicator


The reference indicator (3.0) may be adjusted to the patient's individual reference point by loosening the reference indicator screw (3.1).


#### 5.4 3D-registration joint


The 3D-bite fork (8.0) is attached to the 3D-bite fork connection (7.5) by means of the 3D-bite fork connector screw (7.4).

The bite fork (older version) can be attached to the connector (7.3), as the connector is equipped with the corresponding adapter on the opposite side. For this purpose, the connector screw (7.4) has to be changed to the opposite side.


The 3D-registration joint (7.0) is secured on the multifunctional mounting (1.6) with the registration joint connection screw (1.11).


The tommy screw (7.2) secures the 3D-registration joint.

#### 5.5 Adjusting the width


The width of the transferbow can be adjusted by loosening the 'width setting' screw. If the minimum width is not sufficient, the black earpieces 'L' (9.0) can be used. In this way, an additional 20 mm of space is available.

#### 5.6 Additional height in the area of the nosepiece

If the nosepiece rod (2.2) is rotated by 180°, an additional 10 mm is gained.

## 6. Practical Use on the Patient

### 6.1 Preparing the bite fork

Coat the bite fork with sufficient hard wax (warmed) or silicone. Press wax or silicone on the maxillary tooth row. Make sure that clear occlusal support on both sides is provided.

### 6.2 Reference plane

#### 6.2.1 CP version

Camper's plane runs from the spina nasalis to the lower edge of the external auditory meatus. Mark the spina nasalis.

#### 6.2.2 FH version

The Frankfort horizontal runs from the suborbital point to the upper outer edge of the auditory meatus. Mark the orbital point.

### 6.3 Positioning the transferbow

The patient should be in an upright position. Loosen the width setting screw (1.7) and completely open the transferbow. Introduce the earpieces (6.3) into the external auditory meatus by slowly pushing the side arms together. To relieve some of the weight, the patient may support the UTS 3D by holding the side arms with both hands.

Do not apply too much pressure when positioning the earpieces.


Next, tighten the width setting screw (1.7) to set the width.

## 6.4 Positioning the nosepiece


With the 'nosepiece' screw open (2.4), exert slight pressure with the thumb to position the nosepiece rod (2.2) on the glabella so that the bow sits firmly.


## 6.5 Aligning the transferbow

Align the transferbow by vertical adjustment of the nosepiece, pointing with the reference indicator either at the spina nasalis (for CP) or at the orbital point (for FH).

**CP version = spina nasalis**

**FH version = orbital point**

When manipulating the reference indicator (3.0), make sure to observe a certain distance to the eye of the patient in order to prevent injury.

## 6.6 Other versions

By means of the vertically adjustable earpieces, the UTS 3D-transferbow is also suitable for model transfer according to specific coordinates. Please refer to the literature for further details on this subject.

## 6.7 Securing the 3D-registration joint

Check again if the transferbow is tightly positioned with regard to the relevant reference points. Attach the 3D-bite fork (8.0) to the open 3D-registration joint (7.0) and position the 3D-bite fork on the mandibular tooth row. Additional support for the 3D-bite fork may be provided by cotton rolls. The 3D-registration joint (7.0) is now secured at the 3D-registration joint connection (1.8).

Tighten the tommy screw (7.2) for the rotary joint. The registration procedure is now complete and the corresponding parameters are recorded in the 3D-registration joint.


## 6.8 Removing the transferbow

Loosen the width setting screw (1.7) for the nosepiece (1.9). Pull the side arms outward and instruct the patient to open the mouth. Now the whole apparatus including the registration can be removed.

## 6.9 Removing the 3D-registration joint

Loosen the registration joint connection screw (1.11) and remove the 3D-registration joint (7.0).



The transferbow may be used in combination with the Gnathometer 'M' or the Centric Tray in the manner described above.

## Important note on the UTS 3D registration joint

Please make sure that the tommy screw at the UTS 3D registration joint is always attached and fixed on the right hand side. Otherwise, problems may occur with the registration joint holder (CP or FH version) upon transfer of the data in the articulator.

## 7. Maintenance, Cleaning, and Diagnosis

This chapter describes the user maintenance and cleaning procedure for the UTS 3D-transferbow. This list contains only tasks that may be performed by dental personnel. All other tasks must be performed by qualified service personnel at a certified Ivoclar Vivadent Service Center.

### 7.1 Inspection and maintenance

The time for these maintenance procedures depends on the frequency of use and the working habits of the users. For that reason, the recommended times are only approximate values.

### 7.2 Cleaning

What	When	Cleaning material
Screws that are contaminated with wax or plaster.	If required	Rinse with hot water and dry with a dry cloth
Screws that are contaminated with silicone.	If required	Clean with a cloth.
3D-bite fork that has come into contact with saliva or blood.	After every use	see chapter 7.3
3D-registration joint	After every use	Clean with a spray disinfectant
UTS 3D-transferbow	After every use	Clean with a spray disinfectant
Nosepads	After every use	Use new nosepad (nosepads are disposable)
Earpieces	After every use	Autoclaving device. Please observe the corresponding instructions.

Avoid contact with strong acids or solvents (e.g. MMA). In this way, damage to the surfaces can be prevented.

### 7.3 Cleaning instructions for the 3D-registration joint

- After every use or contact with saliva or blood, wipe with a disinfectant.
- Only use disinfectants that are suitable for stainless steel and light metals.
- Avoid any contact with strong acids, bases or solvents.
- Do not clean with ultrasonic bath, water, steam jet or sterilization units.

## 8. What if....

This chapter will help you to recognize errors and take appropriate measures or, if possible, to perform some repairs.

### 8.1 Technical errors

Description	Double check	Corrective action
Patient has pressure sores and complains about pain of the nose	Is the nosepad in place?	Use nosepad
Patient complains about pain in the ear	Is the width of the transferbow too narrow?	Do not close the transferbow too tightly.
The reference to the corresponding plane (CP or FH) has not been correctly recorded	Were the earpieces correctly mounted so that they came into contact with the earpiece support?	When recording according to the CP and FH method, the earpieces must always come into contact with the earpiece support.
Surface of the transferbow is damaged	Was the surface cleaned with acid or a solvent?	Do not use acids or solvents to clean the transferbow.
Surface of the transferbow is scratched	Was a sharp instrument used?	Do not use sharp instruments in conjunction with the transferbow.

### 8.2 Repair

Repairs may only be carried out by a certified Ivoclar Vivadent Service Center. Please refer to the addresses on the last page of these Instructions for Use.

If repairs during the warranty period are not carried out by a certified Ivoclar Vivadent Service Center, the warranty will expire immediately. Please refer to the warranty regulations for more details.

## 9. Product Specifications

### 9.1 Delivery form

- 1 Universal transferbow with set screws
- 1 3D-registration joint with tommy screws
- 1 3D-Bite fork
- 1 Nosepiece
- 5 Nosepads, green (medium)
- 5 Nosepads, blue (soft)
- 2 Earpieces, white
- 2 Earpieces, black (large)
- 1 Reference indicator
- 2 Axis pins
- 1 Instructions for Use

**The delivery form may vary from country to country.**

#### Recommended accessories

- UTS type 2 support pins for direct model transfer to the articulator

#### Additional accessories and spare parts

- 3D-Bite fork (8.0)
- 3D-Registration joint (7.0)
- Nosepads 'medium' (10.1)
- Nosepads 'soft' (10.2)
- Earpieces (6.3)
- Earpieces 'L' (9.0)

### 9.2 Technical data

#### Weights:

Nosepiece:	25 g
3D-Registration joint:	94 g
UTS 3D-transferbow:	160 g

#### Materials:

Metal parts:	aluminium, steel
Side arms:	carbon composite material
Earpieces:	DELTRIN resin

#### Adjustment possibilities

Adjustable width with	
- earpieces 'L':	85 mm to 155 mm
- earpieces "S":	105 mm to 175 mm
Adjustable width:	87 mm to 160 mm
Adjustable nosepiece height:	85 mm
Adjustable nosepiece length:	22 mm
Adjustable earpiece height:	22 mm
Adjustable reference indicator height:	66 mm

## 10. Miscellaneous

### Tips for model transfer to the Stratos articulator

Please observe the Operating Instructions of the Stratos articulator or the corresponding literature.

This device has been developed solely for use in dentistry. Start-up and operating should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from misuse or failure to observe the Instructions. The user is solely responsible for testing the apparatus for its suitability for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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© Ivoclar Vivadent AG, FL-9494 Schaan
## **Ivoclar Vivadent Limited**
## **Ivoclar Vivadent, Inc.**