# Pretreatment and cementation

Firing						Polishing							
<p><b>IPS e.max® CAD</b></p>  <p>Lithium disilicate glass-ceramics</p>						<p><b>IPS e.max® ZirCAD</b></p>  <p>Zirconium oxide ceramics</p>		<p><b>IPS Empress® CAD</b></p>  <p>Leucite glass-ceramics</p>		<p><b>Tetric® CAD</b></p>  <p>Composite</p>		<p><b>Telio® CAD</b></p>  <p>PMMA</p>	
<b>Flexural strength</b>	530 MPa <sup>1)</sup>					MT Multi: 850 MPa <sup>2)</sup> LT: 1 200 MPa <sup>2)</sup>		185 MPa <sup>1)</sup>		272 MPa <sup>2)</sup>		135 MPa <sup>2)</sup>	
<b>Types of restorations</b>	Veneers, inlays, onlays, minimally invasive crowns (min. 1.0 mm)		Crowns (min. 1.5 mm), 3-unit bridges up to the 2 <sup>nd</sup> premolar			Crowns, 3-unit bridges		Veneers, inlays, onlays, crowns		Occlusal veneers, veneers, inlays, onlays, crowns		Temporary crowns and bridges	Long-term crowns and bridges
<b>Cementation method</b>	adhesive		self-adhesive	conventional		adhesive	self-adhesive/ conventional	adhesive		adhesive		temporary	adhesive
<b>Blasting</b>	–					Al <sub>2</sub> O <sub>3</sub> , 25–70 µm, 1 bar or Al <sub>2</sub> O <sub>3</sub> , 70–110 µm, 1,5 bar		–		Al <sub>2</sub> O <sub>3</sub> , 25–70 µm, 1 bar or Al <sub>2</sub> O <sub>3</sub> , 70–110 µm, 1,5 bar		Al <sub>2</sub> O <sub>3</sub> , 25–70 µm, 1 bar or Al <sub>2</sub> O <sub>3</sub> , 70–110 µm, 1,5 bar	
<b>Etching</b>	<b>Option 1:</b> Agitate Monobond Etch & Prime® for 20 s and allow it to react for another 40 s	<b>Option 2:</b> 20 s with IPS® Ceramic Etching Gel	<b>Option 1:</b> Agitate Monobond Etch & Prime® for 20 s and allow it to react for another 40 s	<b>Option 2:</b> 20 s with IPS® Ceramic Etching Gel	20 s with IPS® Ceramic Etching Gel	–		<b>Option 1:</b> Agitate Monobond Etch & Prime® for 20 s and allow it to react for another 40 s	<b>Option 2:</b> 60 s with IPS® Ceramic Etching Gel	–		–	
<b>Conditioning</b>	60 s with Monobond® Plus	–	60 s with Monobond® Plus	–	–	60 s with Monobond® Plus	–	60 s with Monobond® Plus	–	20 s with Adhese® Universal	30 s with Multilink® Primer A+B	–	2–3 min with SR® Connect
<b>Cementation system</b>	Variolink® Esthetic, Multilink® Automix <sup>3)</sup>		SpeedCEM® Plus		Vivaglass® CEM	Multilink® Automix	SpeedCEM® Plus, Vivaglass® CEM	Variolink® Esthetic, Multilink® Automix <sup>3)</sup>		Variolink® Esthetic	Multilink® Automix <sup>3)</sup>	Telio® Link	Variolink® Esthetic, Multilink® Automix

<sup>1)</sup> Average biaxial flexural strength, over a period of 10 years, R&D Ivoclar Vivadent, Schaan, Liechtenstein

<sup>2)</sup> Typical mean value of biaxial flexural strength, R&D Ivoclar Vivadent, Schaan, Liechtenstein

<sup>3)</sup> Not recommended for veneers.

CE 0123

Manufacturer:  
Ivoclar Vivadent AG, 9494 Schaan/Liechtenstein  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

see instructions

Date information prepared: 2019-11-08, Rev. 1  
702850

more information  
[www.cementation-navigation.com](http://www.cementation-navigation.com)

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