# ProArt CAD Try-In

## EN Instructions for Use

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For dental use only.

**Description**

- Functional try-in mock-ups for removable denture prosthetics
- Functional try-in mock-ups for fixed denture prosthetics
- Impression and functional trays for denture prosthetics

**Indication**

Suitable for use in the treatment of partially and completely edentulous patients

**Contraindications**

- The use of ProArt CAD Try-In is contraindicated if the patient is known to be allergic to any of the product's ingredients.

**Limitations of use**

- Permanent restorations
- Temporary restorations
- Abutment mock-ups

**Side effects**

Systemic side effects are not known to date. In individual cases local allergic reactions to methyl methacrylate materials have been reported.

**Composition**

Polymethyl methacrylate

Contains: Methyl methacrylate, ethylene glycol dimethacrylate, dibenzoylperoxide

**Application****Processing with CAD/CAM equipment**

When designing the restoration, the specified minimum wall thicknesses and parameters of the final restoration must be complied with. If the minimum values are not observed, this can lead to failure (e.g. fracture).

ProArt CAD Try-In discs can be processed in CAD/CAM milling systems with a standard holder of 98.5 mm. Only use milling tools recommended for this milling system when processing the material. Please note that the material can be subject to overheating if milling tools, polishing brushes or steam cleaners are handled incorrectly. This may result in the material being damaged.

**Processing after the milling procedure****Finishing**

After milling, separate the completed object from the disc with a tungsten carbide bur. Use a cross-cut tungsten carbide bur for finishing or shape adjustments, and use PMMA polishing equipment if necessary.

The impression tray made of ProArt CAD Try-In can be used immediately for impression taking. High gloss polishing is not required due to the short time during which the material remains in the mouth.

**Adjustment options**

If necessary, impression trays made of ProArt CAD Try-In can be adjusted using auto-polymerizing tray materials (e.g. SR Ivolen®). To do so, roughen or sandblast the surface. Then, clean it and adapt the material.

**Impression materials**

If the fit of the impression tray is checked on the patient before impression taking, the tray must be cleaned under running water after the try-in. After cleaning and drying the impression tray, it is ready for impression taking. All impression materials which do not require moist storage can be used for impression taking with the respective adhesive.

**Cleaning and disinfection**

Clean the impression tray together with the impression material under running water. Use a suitable commercial disinfecting agent for disinfection.

**Important processing restrictions**

The following points have to be observed for the successful working with ProArt CAD Try-In:

- Observe the required minimum thickness.
- Only mill the discs in a compatible CAD/CAM system.
- Failure to observe the stipulated limitations of use and processing instructions may lead to failure.

**Safety information**

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: [www.ivoclarvivadent.com](http://www.ivoclarvivadent.com), and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website ([www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)).

**Warnings**

Do not inhale grinding dust during finishing – use extraction equipment and wear a mask.

**Disposal information**

Remaining stocks must be disposed of according to the corresponding national legal requirements.

**Shelf life and storage**

No special storage conditions are required.

**Additional information**

Keep material out of the reach of children!

The material has been developed solely for use in dentistry.

Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions.

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### Ivoclar Vivadent Limited
### Ivoclar Vivadent, Inc.