EN Instructions for Use

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Description

AdheSE is a self-etching light-curing filled two-component adhesive system for enamel and dentin. It comprises

1. AdheSE Primer
2. AdheSE Bond

Composition

AdheSE Primer: dimethacrylate, phosphonic acid acrylate, initiators and stabilizers in an aqueous solution

AdheSE Bond: HEMA, dimethacrylate, silicon dioxide, initiators and stabilizers

Indications

- Adhesive for direct light-curing composite and compomer restorations
- After light-curing of the bonding component for build-ups with dual or self-curing composites

Contraindication

If the patient is allergic to any of the components of AdheSE or if the stipulated working technique cannot be used.

Direct pulp capping. Not indicated for adhesive cementation of indirect restorations.

Side effects

Components of AdheSE may cause a sensitizing reaction in patients with a hypersensitivity to any of the material's ingredients. In these cases, the material should not be used.

Interactions

Substances containing eugenol/clove oil inhibit the polymerization of AdheSE.

Rx ONLY

Date information prepared

2014-02/Rev. 3
591948/WE3

CE 0123

Manufacturer
Ivoclar Vivadent AG
FL-9494 Schaan
Liechtenstein

ivoclar
vivadent
clinical

Application

A: Direct restorations – light-curing materials

1. A dry operating field must be ensured. Ideally, a rubber dam should be used.
2. Apply an adequate amount of AdheSE Primer with a brush. All the cavity surfaces, beginning at the enamel, should be adequately wetted. If necessary, fresh primer should be applied. Once the cavity surfaces have been completely coated, the primer must be brushed into the entire surface for another 15 s. The total reaction time should not be shorter than 30 s.
3. Disperse excess AdheSE Primer with a strong stream of air until the mobile liquid film is no longer visible.
4. Apply AdheSE Bond, beginning at the dentin.
5. Disperse AdheSE Bond with a very weak stream of air. Avoid pooling! Care should be taken not to blow away the bonding agent. All the prepared surfaces must be adequately coated with AdheSE Bond.
Note: AdheSE Bond does not contain a solvent that has to be evaporated!
6. Polymerize AdheSE Bond for 10 s using a light intensity of more than 500 mW/cm².
If other lights are used, the recommendations of the manufacturer of the device used must be followed to obtain fully polymerized results. The light-curing time depends on the intensity and the wavelength range of the light emitted by the different curing units.
Note: Light-curing times of less than 5 s are not recommended.
7. The restorative material must be applied according to the instructions of the manufacturer.

Warning

AdheSE is irritating. Avoid contact of the material with the skin, mucous membrane and eyes. If the product comes in contact with the skin, wash it off with plenty of water. In case of accidental contact with the eyes, wash immediately with plenty of water and seek medical attention. Commercial medical gloves do not protect against sensitivity towards methacrylate.

Storage

- See label / packaging for the expiry date
- Storage temperature 2–28 °C / 36–82 °F

Keep out of the reach of children.

For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.