Tetric EvoFlow®

Tetric EvoFlow®

- EN Instructions for Use**
– Light-curing, flowable resin-based dental restorative material
- DE Gebrauchsinformation**
– Lichthärtendes, fließfähiges, zahnrärztliches Füllungscomposito
- FR Mode d'emploi**
– Matériau composite fluide à base de résine photopolymérisable pour restauration dentaire
- IT Istruzioni d'uso**
– Composito da restauro diretto fluido, ad uso odontoiatrico
- ES Instrucciones de uso**
– Material de restauración dental basado en resina fluida fotopolimerizable
- PT Instruções de Uso**
– Material de restauração dental, baseado em resina, fluido e fotopolimerizável
- DA Brugsanvisning**
– Lysbærende, flydende, plastfyldningsmateriale til tænder
- FI Käyttöohjeet**
– Valokovetteinen, juokseva, yhdistelmämuovipohjainen hammastäyttemateriaali
- NO Bruksanvisning**
– Lysbærende, flytende odontologisk fyllingskompositt
- NL Gebruiksaanwijzing**
– Lichtuithardend, vloeibaar tandheelkundig restauratiemateriaal op basis van kunststof
- EL Οδηγίες Χρήσεως**
– Φωτοπολυμεριζόμενο λεπτόρρευστο οδοντιατρικό υλικό αποκαταστάσεων με βάση ρητίνη
- TR Kullanma Talimatı**
– Işıklı sertleşen, akıcı, rezin esaslı dental restoratif materyal

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Description

Tetric EvoFlow is a flowable, light-curing radiopaque nano-hybrid composite for the restorative therapy and the cementation of ceramic and composite restorations. Tetric EvoFlow cures with light in the wavelength range of 400–500 nm.

Shades

Tetric EvoFlow is available in the following 14 shades:

- Enamel shades: A1, A2, A3, A3.5, A4, B1, B3, T
- Dentin shades: A2 Dentin, A 3.5 Dentin, B2 Dentin
- Bleach shades: Bleach XL, Bleach L, Bleach M

Composition

The monomer matrix is composed of dimethacrylates (38 wt.%). The fillers contain barium glass, ytterbium trifluoride, highly dispersed silicon dioxide, mixed oxide and copolymer (62 wt.%). Additives, initiators, stabilizers and pigments are additional ingredients (<1.0 wt.%). The total content of inorganic fillers is 57.5 wt.% or 30.7 vol.%. The particle sizes of the inorganic fillers range between 40 nm and 3 µm.

Indication

- As initial layer under Class I to V restorations
- Anterior restorations (Classes III, IV)
- Class V restorations (cervical caries, root erosion, wedge-shaped defects)
- Small restorations of all types
- Extended fissure sealing
- Splinting of mobile teeth
- Blocking out of undercuts
- Adhesive cementation of light-transmissive indirect composite and ceramic restorations
- Repair of composite/ceramic veneers

Contraindication

The use of Tetric EvoFlow is contraindicated:

- If a dry working field cannot be established, or the stipulated working technique cannot be applied.
- If a patient is known to be allergic to any of the ingredients of Tetric EvoFlow.

Side effects

In individual cases, components of Tetric EvoFlow may lead to sensitization. Tetric EvoFlow should not be used in such cases. To avoid possible irritation of the pulp, areas close to the pulp should be protected with a suitable pulp/dentin protector (selectively apply a calcium hydroxide-based preparation in areas close to the pulp).

Interactions

Substances containing eugenol/clove oil may inhibit the polymerization of composite materials. Consequently, the application of such materials together with Tetric EvoFlow must be avoided. Cationic mouth rinses and plaque disclosing agents as well as chlorhexidine may cause discolouration.

Application

1. Shade determination

Clean the teeth prior to shade determination. The shade is selected with the tooth still moist.

2. Isolation

Adequate relative or absolute isolation using suitable auxiliaries, such as OptraGate® or OptraDam® Plus, is required.

3. Cavity preparation

The cavity is prepared according to the principles of the adhesive technique, i.e., by preserving as much of the tooth structure as possible. Do not prepare sharp internal edges and angles or additional undercuts in caries-free areas. The dimensions of the cavity are generally determined by the extent of the caries or the size of the old restoration. Bevel enamel margins of anterior teeth. In the posterior region, only the sharp enamel edges should be rounded (finishing diamonds, 25–40 µm). Caries-free cervical defects are not prepared, only cleaned with pumice or other suitable cleaning pastes using rubber cups or rotary brushes. Subsequently, remove all residue in the cavity with water spray and dry with water- and oil-free air.

4. Pulp protection / Base

Do not apply a base material when using an enamel/dentin bonding agent. Selectively coat only very deep areas close to the pulp with a calcium hydroxide material (e.g. ApexCal®) and subsequently cover with a pressure-resistant cement (e.g. a glass-ionomer cement, such as Vivaglass® Liner). Do not cover the other cavity walls so that they can be used to support the bond with an enamel/dentin adhesive.

5. Apply matrix / interdental wedge

Use a wrap around matrix for cavities affecting the proximal area or a sectional matrix and wedge it.

6. Conditioning / Application of the bonding agent

Condition and apply the bonding agent according to the Instructions for Use of the product in use. Ivoclar Vivadent recommends using the universal bonding agent Adhese® Universal, a single-component, light-cured adhesive for direct and indirect bonding procedures that features compatibility with all etching techniques.

7. Application of Tetric EvoFlow

- For an optimum result, apply Tetric EvoFlow in layers of max. 2 mm (or 1.5 mm Dentin shades, Bleach XL) and adapt with a suitable instrument (e.g. OptraSculpt).
- Sufficient exposure to the curing light prevents incomplete polymerization. For the recommendations regarding exposure time and light intensity see table 1.
- When using a metal matrix, additionally polymerize the composite material from the buccal or the lingual/palatal aspect after removing the matrix.

8. Finishing / Checking the occlusion / Polishing

Remove excess material with diamond or tungsten carbide finishers after polymerization. Check the occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restorations. Use silicone polishers (e.g. OptraPol®) as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Cementation of ceramic and composite restorations

Tetric EvoFlow is a purely light-curing material. Consequently, the composite is only suitable for the cementation of restorations that are sufficiently translucent and have an adequate thickness to permit enough polymerization light to pass through the restoration and reach Tetric EvoFlow so that the composite can cure completely. Before being cemented in place, the restoration has to be pre-treated according to the instructions of the manufacturer.

Additional notes

1. Tetric EvoFlow can be used in combination with Tetric EvoFlow Bulk Fill, Tetric EvoCeram Bulk Fill and Tetric EvoCeram.
2. In the case of repairs, additional Tetric EvoFlow can be directly applied to the polymerized material. If the Tetric EvoFlow restoration has already been polished, it must first be roughened and wetted with Adhese Universal before fresh Tetric EvoFlow is applied.
3. Tetric EvoFlow should have ambient temperature when applied. Cool temperatures render the material difficult to extrude.
4.  For single use only. If Tetric EvoFlow is applied from the Cavifil directly in the mouth of the patient, the Cavifil or the application cannula should only be used for one patient due to hygienic reasons (prevention of cross-contamination between patients).
5. Syringes or Cavifils should not be disinfected with oxidizing disinfection agents.

Warning

Avoid contact of unpolymerized Tetric EvoFlow with skin, mucous membrane or eyes. Unpolymerized Tetric EvoFlow may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)
- Close syringes/Cavifils immediately after usage. Exposure to light leads to premature polymerization.
- Do not use Tetric EvoFlow after the indicated date of expiration.
- Date of expiration: see information on Cavifils, syringes and packages.

Keep out of the reach of children!

For use in dentistry only!

The material has been developed solely for use in dentistry and should be used strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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