Description

Tetric® is a light-curing, radiopaque composite for restorative treatment. Tetric cures with light in the wavelength range of 400–500 nm.

Composition of Tetric

The monomer matrix is composed of Bis-GMA, urethane dimethacrylate and TEGDMA (total monomer content 18.0 wt%). The inorganic fillers are barium glass, ytterbium trifluoride, silicon dioxide and mixed oxide. Additives, initiators, stabilizers and pigments are additional contents (0.2 wt%). The total content of fillers is 82.0 wt% or 65.0 vol%. The particle size of inorganic fillers is between 40 nm and 7 µm.

Indication

- Class I–V restorations
- Restoration of deciduous teeth
- Extended fissure sealing
- Direct veneering
- Splinting of mobile teeth
- Blocking out of undercuts
- Repair of composite/ceramic veneers
- Inlays/onlays with extraoral post-tempering

Contraindication

The placement of Tetric restorations is contraindicated

- if a dry working field cannot be established or the prescribed application technique cannot be applied;
- if a patient is known to be allergic to any of the ingredients of Tetric

Side effects

In rare cases, components of Tetric may lead to sensitization. The product must not be used in such cases. To avoid possible irritations of the pulp, areas close to the pulp should be protected with a suitable pulp/dentin protector (apply a calcium hydroxide-based preparation to areas in close proximity of the pulp and cover it with an adequate lining).

Interactions

Phenolic substances such as eugenol/clove oil inhibit the polymerization of materials based on methacrylate. Consequently, the application of such materials together with Tetric must be avoided. Discolouration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.

Application**1. Shade selection**

Clean the teeth prior to shade determination. The shade is selected with the tooth still moist.

2. Isolation

Adequate relative or absolute isolation using suitable auxiliaries, such as OptraGate® or OptraDam® Plus, is required.

3. Cavity preparation

The cavity is prepared according to the principles of the adhesive technique, i.e. by preserving as much of the tooth structure as possible. Do not prepare sharp internal angles or additional undercuts in caries-free areas. The geometry of the cavity is predetermined by the dimensions of the caries lesion or the old filling. Bevel enamel edges of anterior teeth. In the posterior region, only the sharp enamel edges should be lightly broken or rounded (finishing diamonds, 25–40 µm). Caries-free cervical defects are not prepared, only cleaned with pumice or other suitable cleaning pastes with the help of rubber polishers or rotary brushes. Subsequently, rinse the cavity with water spray to remove all the residue and dry with water- and oil-free air.

4. Pulp protection / Base

Do not apply a base material when using an enamel/dentin bonding agent. Only cover very deep areas close to the pulp with a calcium hydroxide material (e.g. ApexCal®) and subsequently use a pressure-resistant cement (e.g. a glass ionomer cement, such as Vivaglass® Liner). Do not cover other cavity walls, since they can be used to support the bond with an enamel/dentin adhesive.

5. Placement of matrix / interdental wedge

Use a wrap around matrix for cavities affecting the proximal area or a sectional matrix and wedge it.

6. Conditioning / Application of the bonding agent

Condition and apply the bonding agent according to the Instructions for Use of the product in use. Ivoclar Vivadent recommends using the bonding agent Syntac® or ExciTE® F (both in conjunction with phosphoric acid etching), or Adhese® Universal (in the self-etching mode or in conjunction with phosphoric acid etching).

7. Application of Tetric

- In order to achieve optimum results, Tetric should be applied in increments of max. 2 mm or 1.5 mm for dentin shades and max. 4 mm for Tetric Basic White and adapted to the cavity walls with a suitable instrument (eg. OptraSculpt®).
- Sufficient exposure to the curing light prevents incomplete polymerization. For the recommendations regarding exposure time and light intensity see table 1.
- When using a metal matrix, additionally polymerize the composite material from the vestibular or lingual/palata aspect after removing the matrix, if no Bluephase polymerization light is used or the light probe cannot be ideally positioned, e.g. distant to composite or diverging scattering angle.

8. Finishing / Checking the occlusion / Polishing

After polymerization, remove excess material with suitable finishers or fine diamonds. Remove proximal excess with diamond carbide finishers, finishing strips or flexible finishing discs. Check the occlusion and articulation and grind in the necessary adjustments to prevent premature contacts or undesired articulation paths on the surface of the restoration.

Tetric®

EN Instructions for Use Light-curing resin-based dental restorative material	DA Brugsanvisning Lyshærdende plastfyldningsmateriale til tænder
DE Gebrauchsinformation Lichthärtendes zahnärztliches Füllungscomposit	FI Käyttöohjeet Valokovetteinen yhdistelmämuovityttemateriaali
FR Mode d'emploi Composite de restauration dentaire photopolymérisable	NO Bruksanvisning Lysherdende odontologisk fyllingskomposit
IT Istruzioni d'uso Composito fotoindurente a base resinosa per restauri dentali	NL Gebruiksaanwijzing Lichtuithardende composit voor restauratieve tandheelkunde
ES Instrucciones de uso Material restaurativo dental fotopolimerizable en base a resina	EL Οδηγίες Χρήσεως Φωτοπολυμερίζόμενο οδοντιατρικό αποκαταστατικό υλικό με βάση την ρητίνη

Use silicone polishers (e.g. OptraPol®) as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Additional notes

1. In the case of repairs, additional Tetric can be directly applied to the polymerized material. If the Tetric restoration has already been polished, it has to be roughened and wetted with Heliobond before a new layer of Tetric can be applied.
2. Use Tetric at room temperature. Cold material may be difficult to dispense.
3.  For single use only. If Tetric is directly applied from the Cavifil in the mouth of the patient, the Cavifil should only be used for one patient due to hygienic reasons (prevention of cross-contamination between patients).
4. Do not use oxidizing disinfectants to disinfect syringes and Cavifils.
5. The recommended thickness of the increment of 4 mm with Tetric Basic White is based on hardness measurements.

Warnings

Unpolymerized Tetric should not come in contact with skin, mucous membrane and eyes. Unpolymerized Tetric may have a slight irritating effect and may lead to a sensitization against methacrylates. Customary medical gloves do not provide protection against the sensitizing effects of methacrylates.

Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)
- Close syringes/Cavifils immediately after usage. Exposure to light causes premature polymerization.
- Do not use Tetric after the indicated date of expiration.
- Expiry date: see note on Cavifils, syringes and packages.

Keep out of the reach of children!

For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent Marketing Ltd.
Ivoclar Vivadent AB
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.