


# Ivoclar Vivadent All-Ceramics | Vollkeramik



# Ivoclar Vivadent All-Ceramics | Vollkeramik

## Indications | Indikationen

	IPS e.max®					IPS Empress®	
	IPS e.max® Press	IPS e.max® CAD	IPS e.max® ZirCAD	IPS e.max® ZirPress	IPS e.max® Ceram	IPS Empress® Esthetic	IPS Empress® CAD
Thin Veneers I Dünne Veneers	✓	✓					
Veneers	✓	✓		✓	✓	✓	✓
Occlusal veneers I Okklusale Veneers	✓	✓					
Inlays / Onlays	✓	✓				✓	✓
Partial crowns I Teilkronen	✓	✓			✓ <sup>2)</sup>	✓	✓
Minimally invasive crowns <sup>4</sup> I Minimal-invasive Kronen <sup>4</sup>	✓	✓	✓				
Crowns I Kronen	✓	✓/✓ <sup>3</sup>	✓	✓ <sup>1)</sup>	✓ <sup>2)</sup>	✓	✓
3-unit anterior bridges I 3-gliedrige Frontzahnbrücken	✓	✓/✓ <sup>3</sup>	✓	✓ <sup>1)</sup>	✓ <sup>2)</sup>		
3-unit premolar bridges I 3-gliedrige Prämolarenbrücken	✓	✓/✓ <sup>3</sup>	✓	✓ <sup>1)</sup>	✓ <sup>2)</sup>		
3-unit posterior bridges I 3-gliedrige Seitenzahnbrücken		✓ <sup>3)</sup>	✓	✓ <sup>1)</sup>	✓ <sup>2)</sup>		
3- and multiple-unit bridges I 3- und vielgliedrige Brücken		✓ <sup>3)</sup>	✓	✓ <sup>1)</sup>	✓ <sup>2)</sup>		
Hybrid abutments I Hybrid-Abutments	✓	✓					
Hybrid abutment crowns I Hybrid-Abutment-Kronen	✓	✓					

1) In combination with IPS e.max ZirCAD I in Kombination mit IPS e.max ZirCAD

2) One layering ceramic for all IPS e.max components I eine Schichtkeramik für alle IPS e.max-Komponenten

3) IPS e.max Veneering Solutions (CAD-on technique) in combination with IPS e.max ZirCAD I IPS e.max Veneering Solutions (CAD-on-Technik) in Kombination mit IPS e.max ZirCAD

4) Material thicknesses vary depending on the material I Schichtstärken variieren in Abhängigkeit vom Material



Please read the corresponding Instructions for Use  
Bitte entsprechende Gebrauchsinformation beachten

Date information prepared: 2017-05-30  
627459