



Furnace accessories ideally complement your furnace

Programat Accessories Set 2

610 164

The Programat Accessories Set is an ideal supplement to the Programat furnaces and comprises the following accessories:

- 1 Programat Firing Tray
- 1 Silicon Nitride Firing Tray «K»
- 1 Firing Tongs
- 1 Automatic Temperature Checking Set ATK2 (for P500, P510, P700, PCS, PCS2, EP 3000 and EP 5000)



Programat Firing Tray Kit 2

655 946

The Firing Tray Kit 2 was specially developed for the Programat furnaces. The support pins can be individually mounted to the firing tray if required.

Delivery Form:

- 2 Firing Trays
- 6 Support Pins A
- 6 Support Pins B
- 2 Fan Pins
- 2 Ceramic Sleeves
- 1 Hobbyceram Pen



Ceramic Sleeves for Metal Pins (Refill)

666 068

The ceramic sleeves can be placed on the metal pins to avoid that metal-free restorations come in contact with the metal pins.

Delivery Form:

- 10 Ceramic Sleeves



Metal Pins for Firing Tray (Refill)

626 698

The metal pins can be used for the Programat firing tray and can be mounted to the honeycomb firing tray if required.

Delivery Form:

- 4 Metal Pins A
- 4 Metal Pins B
- 4 Metal Pins C



Firing Tongs

530 665

These firing tongs are used to place objects in the firing chamber. The tongs are characterized by their excellent material quality.

Delivery Form:

- 1 Firing Tongs



Programat WLAN Kit

655 949

The WLAN Kit is used for the wireless data transfer (for P510).



Programat USB Flash Drive

598 601

Individual programs and data can be saved on the USB flash drive.



Programat®

A STORY OF SUCCESS



Furnace accessories ideally complement your furnace

Temperature Checking Set 2

544 903

The Temperature Checking Set is used to check and adjust the temperature in the Programat furnace chamber (e.g. for Programat P100, P200, P300 and P310).

Delivery Form:

2 Sample Holders
10 Silver Test Wires



Silver Test Wires (Refill)

545 634

Delivery Form:

5 Silver Test Wires



Programat Tong Holder xx10

655 945

The Programat tong holder can be mounted to the Programat housing. It can be used to conveniently store the firing tongs.



Automatic Temperature Checking Set ATK2

598 833

The Automatic Temperature Checking Set ATK2 is used to check and adjust the firing temperatures in furnaces with automatic calibration programs for the ATK2 system (for Programat P500, P510, P700, PCS, PCS2, EP 3000 and EP 5000). The automatic double-range calibration is a special feature of the ATK2.



USB Download Cable

597 183

The USB download cable can be used to update the furnace software via a PC/laptop. Additionally, the furnace can be operated or a protocol of the firing data can be generated with the computer using the PrograBase X10 software and the USB download cable.



Programat Display Cleaning Cloth

655 947

This special cloth can be used to clean the display and the glass of the IR sensor.

