

# Programat® CS3

## Software Revision Log

Version	Changes from previous version	Status
<b>V3.2</b> 11.03.2019	<b>Software optimized</b> <ul style="list-style-type: none"> <li>• Optimized parameters for IPS Empress CAD with IPS e.max CAD Crystall.</li> <li>• DSA analysis optimized</li> <li>• DSA protocol function implemented</li> <li>• Optimized vacuum test program</li> <li>• Optimized dehumidification-function</li> <li>• More special characters</li> <li>• Acoustic feedback of touchscreen</li> <li>• WLAN performance has been enhanced</li> <li>• Various languages have been improved</li> <li>• Various software optimizations</li> </ul>	<b>Currently delivered version</b>
<b>V3.1</b> 01.02.2018	<b>Software optimized</b> <ul style="list-style-type: none"> <li>• WLAN performance has been enhanced</li> <li>• Push function can be tested</li> <li>• Multimedia files can be directly deleted</li> <li>• Various languages have been improved</li> <li>• Various software optimizations</li> </ul>	<b>No longer available!</b>
<b>V3.0</b> 28.04.2017	<b>Software optimized</b> <ul style="list-style-type: none"> <li>• New program structure with new program parameters</li> <li>• Various software optimizations</li> </ul>	<b>No longer available!</b>
<b>V2.1</b> 25.09.2016	<b>Software optimized</b> <ul style="list-style-type: none"> <li>• Optimized text boxes for Japanese language</li> <li>• Various software optimizations</li> </ul>	<b>No longer available!</b>

# Programat® CS3

<b>V2.0</b> 08.04.2016	<b>Software optimized</b> <ul style="list-style-type: none"><li>• DSA - Digital Shade Assistant: Orientation of the reference teeth either “incisal edge to incisal edge” or “incisal edge to tooth neck” possible</li><li>• In-process dehumidification: Function that protects the vacuum system from contamination and its consequential damages</li><li>• Various software optimizations</li></ul>	<b>No longer available!</b>
<b>V1.0</b> 15.10.2015	<b>Basic version</b>	<b>No longer available!</b>

Version 11.03.2019