

SERVICE+ | OVERVIEW IMPLANT CONNECTION

Manufacturer*	Implant system*	Implant platform (ø mm)					Abutment level	Single-unit abutments	Implant-supported superstructures
Biomet 3i	Osseotite®	3.4		4.1	5.0	6.0		✓	✓
	Osseotite® Certain®	3.4		4.1	5.0	6.0		✓	✓
CAMLOG®	Camlog®	3.3	3.8	4.3	5.0	6.0		✓	✓
	Conelog®	3.3	3.8	4.3	5.0			✓	✓
	iSy®		3.8	4.3	5.0			✓	✓
Dentsply Implants	Astra Tech Osseo Speed™	3.5/4.0		4.5/5.0				✓	✓
	Frialit® + Xive®	3.4	3.8	4.5/5.5				✓	✓
	Ankylos®						Balance base C / small		✓
medentis medical	ICX	3.45	3.75	4.1	4.8			✓	✓
Nobel Biocare	Brånemark System®	NP 3.5		RP 4.1	WP 5.1			✓	✓
	NobelActive™	NP 3.5		RP 4.3/5.0				✓	✓
	NobelReplace™	NP 3.5		RP 4.3	WP 5.0	6.0		✓	✓
							Multi-unit		✓
Straumann	Bone Level	NC 3.3		RC 4.1	RC 4.8			✓	✓
	synOcta®				RN 4.8	WN 6.5		✓	✓
Zimmer	Screw-Vent®	3.5		4.5	5.7			✓	✓

* The implant systems and products listed above are not registered trademarks of Ivoclar Vivadent AG.